Mécanique

Banc de rotation

Mechanical

Rotation bench

Réf :
322 210

Français – p 1

English – p 8

Version : 4209

Banc de rotation

Rotation Bench

1. Description

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09 69 32 02 10 *

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1. Description

The rotation bench is an educational device designed to study the concepts of rotational mechanics, such as uniform and accelerated circular motions. It includes a mass-pulley system with three levels, connected by a red string, which triggers the rotation of the bench when the mass is released. This generates a torque that sets the rotating arm in motion. This system allows the exploration of rotational dynamics, particularly accelerated movements, and introduces Newton's law of gravitation.

The rotating arm, graduated is equipped with a displacement guide, allows the measurement and adjustment of the mass positions for various dynamic experiments, including the study of moments of inertia. Thus, the rotation bench provides an in-depth exploration of fundamental mechanical principles such as the moment of inertia, the conservation of angular momentum, and the relationships between force, torque, and angular acceleration.

- This rotation bench is:
 1. Robust, ensuring long life even with intensive use.
 2. High repeatability of experiments, crucial for precise and reliable mechanical practical work.
- Modularity and Options:

The bench is designed to be highly adaptable, compatible with various optional modules that extend its capabilities: see chapter "educational application"

These options broaden the scope of study, allowing for:

- Uniform and non-uniform circular motions
- In-depth analysis of moments of inertia
- Detailed exploration of rotational dynamics laws

2. Components

The bench consists of the following elements:

1. **Base platform:** this robust component serves as the support for the rotation bench.

- 2. Main fixation axis:** the axis for rotational fixation, to screw in the center of the platform.

- 3. Multi-level pulley and main rotating arm:** This versatile element is modular with three different diameter levels. This feature allows testing a variety of moments of inertia by adjusting the mass distribution. The rotating arm is equipped with an ingenious displacement guide, allowing the masses to be fixed at different positions.

- 4. Return pulley on the support:** Connected to the mass, this pulley plays a key role in triggering the rotation when the mass is released.

- 5. Notched mass:**

- 6. Cylindrical masses:** These additional masses can be fixed on the main axis for various experiments, with wing screws.


7. **Cord:** used for operating the system, it ensures the transmission of force between the notched mass and the rotational mechanism, enabling efficient conversion of potential energy into rotational kinetic energy.

3. Assembly

The rotation bench is easy to assemble, by simply screwing onto the anodized aluminium plate, dispose a spirit level and adjustable feet to adjust horizontality.

4. Educational Applications

The rotation bench offers a wide range of practical activities, allowing for an in-depth exploration of rotational mechanics concepts:

Fundamental module

Study of Non-Uniform Rotational Motions:

- Objective: Analyze angular acceleration and its effects on motion.
- Method: Use of different masses to create variable torques.
- Key Concepts: Moment of force, angular acceleration, rotational motion equation.

Several optional modules (not provided) allowing for an advanced exploration of concepts

1. Study of Uniform Rotational Motions with an Encoded Motor (322211):
 - Objective: Explore constant angular velocity and its implications.
 - Method: Precise control of speed via the control software.
 - Key Concepts: Angular velocity, centripetal force, dynamic equilibrium.
2. Analysis of Rotational Motions with an Accelerometer (ref 488060):
 - Objective: Measure and interpret accelerations in circular motion.
 - Method: Use of the accelerometer for real-time measurements.
 - Key Concepts: Centripetal acceleration, angular acceleration, angular velocity.
3. Study of Motions with a chronocine Fork (322212):
 - Objective: Precisely determine rotation periods and frequencies.
 - Method: Chronometric measurements using the fork.
 - Key Concepts: Period, frequency, angular velocity, conservation of angular momentum.

These practical activities allow students to develop a deep and practical understanding of the fundamental principles of rotational mechanics, essential in many fields of physics and engineering.

5. Technical Characteristics

-Dimensions : approximately 325*325*180mm

-Base platform: a base for attaching the rotation axis, the motor accessory (not included), and the chronocine fork accessory (not included).

-Assembly: rotation axis, stepped pulley (diameter: 30mm, 45mm and 60mm), graduated rotating arm (length 30cm, graduation in tenths), coloured dot for use with a high-speed camera.

-Return pulley: the return pulley is attached to the mass; releasing it triggers the rotational movement of the axis.

-Cord

-2 cylindrical masses of 50g

-2 cylindrical masses of 100g

-Notched mass of 200g

6. After-Sales Service

The warranty is for 2 years.

For any adjustments, please contact **Technical Support** at **09 69 32 02 10***.

The equipment must be returned to our workshops, and for any repairs or spare parts, please contact:

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