Ex.A.O

Capteur Foxy / Tooxy / AirNeXT

Data logging

Sensor Foxy / Tooxy / AirNeXT

Réf :
482 091

Français – p 1

English – p 8

Version : 5203

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Drop Counter sensor

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09 69 32 02 10*


** appel non surtaxé*

You have just purchased a differential pressure sensor and we thank you for your confidence. This document aims at guiding you in using your sensor.

1. Description

Combined with a pH sensor or conductivity meter, this dropper sensor makes your ph-metric or conductivity dosing and titration a success. Repeatability and traceability of measurements are greatly improved compared with visual determination on a graduated burette.

TP is fully or partially automated, guaranteeing fast results.

2. Use

This sensor can be used with any console in the Foxy range.

- Present the sensor facing the console connector.
- Slide the connector into place, following the console guide.
- Push until it clicks into place.

3. Measure

3.1 Using the sensor in Standalone Mode

3.1.1 Size selection

Select the size using the green size button (10 mL, 100 mL).

The selected size is indicated by a triangle in the green band opposite the size. The volume value is displayed.

3.1.2 Zero adjustment

First press and hold the tare button (purple). A triangle flashes in front of the word "Cal".

Pour a known volume of solvent, drop by drop, through the sensor.


Then turn the purple knob to indicate the volume of solvent poured.

Confirm the setting by clicking on the tare button (purple) again.


When the sensor is correctly calibrated, the arrow above "Cal." goes out.
If the setting does not conform to the factory setting (too great a discrepancy between the volume entered and the number of drops detected), an arrow appears above the Caution symbol.

3.2 Using the sensor with the Scientific Workshop Software


- Run the Standard Scientific Workshop software


- Drag and drop the differential pressure sensor in ordinate:


- The tabs specific to the sensor appear:

- In the caliber tab, select 100mL or 10mL:

- Then click on the Measure tab:

If the value displayed is very different from 0, click on the Settings tab.

3.2.1 Zero adjustment on computer

To set the zero, simply follow the instructions in the Atelier Scientifique software:

Volume

Caliber Quantity Measure Color Regulation

Click on start button to launch calibration process.
All the previous setting will be cancelled .
Dont calibrate your probe during data measurements.

NB : Calibration is needed only to enhance measure accuracy.

Start

Once you have read the instructions, click on Start to bring up the following window:

Volume

Caliber Quantity Measure Color Regulation

Enter the value of calibration,
then click 'Adjust'

mL

Adjust

Pass a known volume through the sensor, then enter the volume value. Click on "Adjust".

The zero setting confirmation window then appears.

Volume

Caliber Quantity Measure Color Regulation

Your probe is calibrated!

4. Technical characteristics

Measurement channels: volume and number of drops

Calibration (volume channel): 10 mL / 100 mL

Display accuracy: 0.01 mL / 0.1 mL

Max. speed: 20 gts /s

Max. temperature probe diameter: 6 mm

Max. electrode diameter: 16 mm

5. After-Sale Services

The warranty is of 2 years.

For all settings, contact our **Technical Support** on **+33232294023**.

The material must be returned to our workshops and for all repairs or spare parts, kindly contact:

JEULIN – S.A.V.

468 Jacques Monod Street
CS 21900
27019 EVREUX CEDEX France

+33 2 32 29 40 23