# Didactique

Kit de chimie

***Didactic***

*Chemistry Kit*

Réfs :  
102 190 – 102 164

Français – p 1

English – p 7

Version : 3111

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# 1 Description

## 1.1 General

Thanks to its innovative and original packaging, the lavender aroma synthesis kit in phials means the teacher can substantially reduce preparation time PT, by minimizing the quantity of glassware required. In addition, the phials are pre-dosed and ready for use. So the reagents can be distributed individually to the students without any need to measure their volume.

The quantities of reagents in the phials have been reduced to a minimum to limit risks linked to handling of the products. The sulfuric acid traditionally used as a catalyst has been replaced by para-toluene sulfonic acid which is less dangerous to handle.

Made of plastic, the break-open phials can be opened by students without any danger or risk of cuts.

## 1.2 Composition and description

The lavender aroma synthesis kit in phials comprises a batch of 10 phials containing 5 ml pure ethanoic acid, 10 phials containing 5 ml of linalool and 10 phials containing paratoluene sulfonic acid which acts as a catalyst.

This kit allows 10 linalyl ethanoate synthesis operations.


Photo 1

- 1 – Breakable tip of phials
- 2 – Label complying with the regulations
- 3 – Phial containing ethanoic acid

- 4 – Phial containing linalool
- 5 – Phial containing paratoluene sulfonic acid

Table 2

## 2 Additional equipment

Synthesizing ethyl propanoate with the phial lavender aroma synthesis kit requires the following accessories:

Designation	Quantity
Reaction chamber with overhead air refrigerant	1
Bath of hot water or bain marie	1
Glass agitator or rod	1
Pill box or glass pot, 60 ml	2
Plastic pasteur pipette	1
Strip of filter paper	1
Pair of protective vinyl or latex gloves	1
Pair of safety goggles	1

Table 3

### For reference 102190

Saturated sodium chloride solution	1
Saturated sodium hydrogeno-carbonate solution	1

Table 3bis

## 3 Precautions for use and safety

### 3.1 Warning

JEULIN cannot be held liable in the case of accidents occurring during use of the lavender aroma synthesis kit under conditions other than those set out in this notice.

Also, JEULIN cannot be held liable in the case of accidents which occur on account of non-compliance with the safety instructions described herein.

### 3.2 Labeling and Safety Data Sheet

The phials containing the reagents are individually labeled in compliance with the regulations in force.

Before any operation, have the students read the label so they become aware of the risks and the protective measures to be adopted.

For additional information, the material material safety Data Sheets for the two reagents can be downloaded from the site [www.jeulin.fr](http://www.jeulin.fr).

Product	FDS Ref.
Lavender aroma synthesis kit in phials – Ethanoic acid	102190-A
Lavender aroma synthesis kit in phials – Linalool	102190-B
Lavender aroma synthesis kit in phials – Paratoluene sulfonic acid	102190-C

Table 4

### 3.3 Individual protection

The phial lavender aroma synthesis kit was formulated to minimize risks linked to the use of chemicals. However, handling pure ethanoic acid requires the wearing of goggles. Wearing gloves is also recommended.

### 3.4 Collective protection

Given the small quantities of reagents used during the operation, it is not necessary to work under an extractor hood.

### 3.5 Waste

The lavender aroma synthesis kit does not contain any compounds which present a particular danger to the environment. However, for a responsible approach to environmental protection the products must be considered as laboratory chemical waste and treated as such. The quantities of reagents used in the kit have been reduced to minimize the volume of waste and the costs linked to its destruction.

## 4 Operating method

### 4.1 Principle of synthesis

The lavender aroma synthesis kit allows synthesizing linalyl ethanoate by esterification. The reaction on which synthesis relies is esterification by acidic catalysis of ethanoic acid on linalool (Figure 5). This balanced reaction results in the formation of linalyl ethanoate and water. The catalyst used is para-toluene sulfonic acid presented in a solid form which is less dangerous to handle than the sulfuric acid normally used.


Figure 5

### 4.2 Preparation of reaction mixture

At each work station assemble the recommended accessories (*paragraph 2*) and an phial of each reagent.

Photo 7

Photo 8


Open carefully and pour the phials containing ethanoic acid and linalool into a reaction chamber (receptacle or test tube) (*photos 7 and 8*).

*Photo 9*

Then pour in the contents of the phial containing the catalyst and agitate using the glass agitator. (*Photo 9*).

### 4.3 Reaction

Next place the air refrigerant on the reaction chamber and immerse in a bath of hot water (water which has brought to boiling point) for 20 to 30 minutes (*Figure 10*).

*Figure 10*

Remove the reaction chamber from the bath of hot water and remove the air refrigerant.

### 4.3.1 Washing the ester

#### Case 1: With water-soluble bags (102 164)

Water-soluble bags used to fill the experimental downtime due to reflux heating and discover original polymer.

##### **Saturated solution of sodium chloride:**

fill 2/3 glass container or a small 100 mL beaker. For a saturated solution of sodium chloride, NaCl dissolve soluble bag (photo A) marked using a rod glass (photo B).

Photo A

Photo B

#### Case 2 : Prior manufacturing solutions

Half fill a glass receptacle or small 50 to 100 mL beaker with a saturated solution of sodium chloride.

##### **In both cases:**

Pour the contents of the reaction chamber into the receptacle containing the saturated sodium chloride solution (*photo 11*).

Agitate using the glass agitator and leave to settle. Two phases will appear, an aqueous phase and above, an organic phase rich in synthesized ester.

Photo 11

Photo 12

### 4.3.2 Neutralization of excess acid

The organic phase retains traces of acids of which the sharp strong odor interferes with identification of the aroma by smelling. It is necessary to neutralize the excess acid to eliminate the odor.

#### Case 1: With water-soluble bags (102 164)

##### **Solution saturated with bicarbonate of sodium:**

Fill 2/3 glass container or a small 100mL beaker. For a saturated solution of sodium bicarbonate dissolve soluble bag marked  $\text{NaHCO}_3$  using a glass rod.

#### Case 2 : Prior manufacturing solutions

For this, one-third fill a glass receptacle or a small 50 to 100 mL beaker with a saturated solution of sodium hydrogen-carbonate.

#### In both cases:

Using a plastic pipette, draw off the organic phase containing the ester (repeating as necessary) and **drop by drop, add** the organic phase to the receptacle containing the sodium hydrogeno-carbonate solution (*photo 12*). There will be strong effervescence indicating the neutralization of the excess acid. Do not add the organic phase too rapidly, to avoid excessively violent effervescence and overspill from the beaker.

### 4.3.3 Characterization of the aroma

Once the excess acid has been neutralized the aroma can be characterized by smelling. For this purpose, take a strip of filter paper and place a few drops of the synthesized aroma on one end. Smell the aroma while waving the strip of filter paper under the nose as if in a perfume shop.

**Do not smell the contents of the bottle directly, since the quantity of ester is substantial and there could be risk of irritating the airways.**

## 5 After sales service

For any questions or complaints, please contact :

#### **JEULIN - TECHNICAL SUPPORT**

468, rue Jacques Monod

CS 21900

27019 EVREUX CEDEX FRANCE

+33 (0) 2 32 29 40 50

Didactic

**Kit for synthesizing lavender aroma in phials and soluble bags**

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# Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediately to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

**Free service\***

**+33 2 32 29 40 50\*\***

\* Call cost not included.  
\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux