Immunologie

Kits

Immunology

Kits

Réf :
106 235

Français – p 1

English – p 5

Version : 8110

Kit Ouchterlony de simulation
Ouchterlony simulation kit

1. Description

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0 825 563 563*

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1. Description

The purpose of the Ouchterlony gel immuno-diffusion test is to demonstrate the antigen-antibody bond and its specific nature.

The actual test uses serum from a small mammal (e.g. a rabbit), previously immunised by injection of a known antigen (generally, BSA beef serum albumin) and various antigen solutions including the one injected. The presence in the serum of antibodies directed against the previously injected antigen and their specific nature are tested using the so-called double diffusion or Ouchterlony method, which is based on two principles:

- diffusion in the gel of serum antibodies from a central well, where the deposit is made, and that of other antigens from peripheral wells, to allow the antigens-antibodies to meet.
- the recognition and fixing of the antibodies on the antigens specific to them. This stage is manifested by the formation of immune complexes (antigen-antibody complexes) which are too large for further diffusion and are visible to the naked eye in the form of a whitish arc.

This type of simulation reproduces the principle of double diffusion.

The results observed are similar to those of the actual test. The meeting in the gel of the NaOH solution (stimulating the serum containing the antibodies) and the solution of zinc sulphate (simulating the specific antigen) is manifested by the formation of a white precipitate simulating the formation of immune complexes.

Advantages of the simulation version:

The use of substitutes has several advantages:

- Since diffusion of the products in the gel is quicker, the results appear after approximately 40 minutes, compared with 24 hours for the actual test. This means the test and analysis of the results is possible during the same session.
- The students do not handle any product of animal origin.
- Preservation is at room temperature.

2. Composition

The kit allows handling up to 20 pairs with 2 sets per pair.

- 15 Petri dishes, diameter 55 mm,
- 1 tube labelled A containing 0.85 g of agar-agar powder,
- 1 tube labelled B containing 1 mL of soda solution,
- 1 tube labelled C containing 1 mL of zinc sulphate solution,
- 4 cryotubes labelled D to G,

A bag of small equipment:

- 1 set of 5 reusable patterns for locating wells,
- 5 punches with suction bulb, reusable.

3. Additional equipment required

To prepare for the operation:

- a microwave oven or bain-marie,
- a beaker, 250 mL,
- a glass agitator or spatula,
- a pipette 5 mL (or more) with pipeteur
- a heatproof glove,
- distilled water (90 mL).

For the operation:

- 6 polyethylene Pasteur pipettes, 1mL (or 3 mL) for depositing the drops,
- permanent felt tip markers (1 per pair),
- a test tube holder.

4. Experiment Application - procedure

4.1 Preparation for the operation

Allow 1h30 for this stage which can be completed a few days in advance.

4.1.1 Preparation of the agar-agar gel

- Pour the contents of tube A (0.85 g of agar-agar) into a 250 mL beaker.
- Add 90 mL of distilled water while shaking the beaker to dissolve the agar-agar.
- Heat the mixture almost to boiling point, in the microwave or bain-marie.

Heat gradually in the microwave, stopping regularly to shake the mixture. Stop heating when the mixture is totally clear, that is at the very start of boiling.

In the bain-marie heat the mixture while stirring until it becomes totally clear and stop at the very start of boiling (approximately 20 minutes).

4.1.2 Filling the dishes with the warm gel

- With the pipette take 5 mL of warm fluid agar gel and place 4 mL (no less) in each dish.
- Immediately shake the dish gently to level the gel and remove any bubbles. Place the cover on the dish.
- Leave the dishes to cool on a flat surface without touching them for the first five minutes necessary for the agar gel to set. Full cooling takes approximately thirty minutes.

Caution: The gel begins to solidify fairly quickly, so do not wait too long before filling the dishes. If the gel is not sufficiently fluid, simply reheat until it is again totally translucent.

Remark: The dishes obtained can be used from the time of solidifying of the gel and can be stored for four to five days in the refrigerator (4 °C) in a damp container to avoid dehydration of the gel.

4.1.3 Preparing the reagents

- No preparation is necessary for tubes B and C which are ready to use. However, they can be relabelled 'tested serum' for tube B and 'BSA' for tube C to hide their true composition from the students.
- Fill the four cryotubes with distilled water, to simulate the other antigen solution. These tubes can also be renamed ovalbumin, lactalbumin, goat, sheep or rabbit serum.
- Using a felt tip pen mark the six Pasteur pipettes with the letters B to G so each tube has its dedicated pipette.

4.2 Handling by pairs

Each pair has two dishes of gel, a sight for marking the wells, a punch and a permanent felt tip marker.

4.2.1 Making wells in the gel

- Place the dish (cover upwards) on the sight which is used as a model and open.
- While firmly gripping the balloon of the punch, pierce the gel to the base of the dish using the sight mark then release the balloon and remove the punch. The disk of gel formed remains in the punch and a well appears in the dish
- Follow the same method for the five other wells.

Remark: If the gel disk does not lift up when the punch is removed, it can easily be removed using a toothpick or a pin.

Caution: It is important to avoid any slits in the gel when making the wells and to comply with the model in terms of the arrangement. An increase in the distance between the central well and the peripheral wells will adversely affect the legibility of the results.

4.2.2 Deposit of reagents (“serums”)

- Before making any deposit, use the felt tip to mark under the dish and at each well, the letter corresponding to the reagent deposited to remember the arrangement of the products. The central well will always be marked B. The peripheral wells will be marked C to G.
- In the central well deposit a drop of solution B ‘tested serum’ using dedicated pipette B without making any bubbles.
- Using the same method, deposit in the peripheral wells a drop from tubes C to G (‘antigens’) using their dedicated pipettes.

Caution: It is important that the wells should be **totally** filled without any bubbles. If the solution runs over onto the gel, it must be removed immediately using a small piece of absorbent paper and if necessary, fill up the well.

4.2.3 Obtaining the results

Leave the dishes flat with their covers on for approximately forty minutes. The wells will gradually dry up, a sign diffusion is occurring, and the results will appear. The results will remain visible for a few days if drying out of the gel is avoided by storing the dishes in a cool damp container.

Suggestion: It is possible to obtain more rapid results by reducing the distance between the central well and the peripheral wells. A distance of one cm between the respective centres of the wells will provide very satisfactory results.

Reading and analysing the results

Reading: A white arc appears in the gel between the central well which contained the tube B solution (‘tested serum’) and the peripheral well which contained the tube C solution (‘BSA’). For improved observation of this arc look at the dish using a transparent or against a black background.

Analysis: The tested serum and the antibodies it may contain are diffused in an aureole around the central well. The same applies to the various antigen solutions around the peripheral wells. The white precipitation arc indicates that at the meeting point between the tested serum and the BSA, a new complex is formed (antigen-antibody complexes). Hence, in the tested serum there are antibodies capable of recognising and fixing themselves on the beef serum albumen (BSA). The absence of the arc of precipitation between the tested serum and the other albumens indicates the antibodies were not able to fix on it. The antibodies are thus specific to one antigen only, that is BSA.

The presence of anti-BSA antibodies in the serum indicates the animal in which the BSA had been injected had produced antibodies to defend itself and the latter are specific to the antigen attacking the animal.

Preservation, sorting and recovery

The kit can be stored at an ambient temperature protected from direct sunlight.
The gel in the dishes can be thrown in the dustbin after possible decontamination with bleach if colonies of bacteria or mould have appeared.
The punches and dishes can be reused after washing.

5. Teaching hints

The Ouchterlony test can be incorporated in various teaching initiatives concerning the study of immune reactions.

- The students know the serum contains antibodies directed against BSA and the test allows demonstrating 'the antigen-antibody bond' and the specific nature of the antibody.
- Since the students know the white arc corresponds to an antigen-antibody bond and that the rabbit had been injected with BSA, the 'synthesis of antibodies by the body following contact with an antigen can be demonstrated' as well as 'the specific nature of antibodies vis-à-vis the antigen received.
- For evaluation purposes, it allows detecting antibodies and/or their specific nature in the serum.

In addition, the two dishes per pair allow performing two different tests such as for example, a test on an immune serum (solution B) and a non-immune serum (mimicked by the tube of distilled water).

6. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux