# **Didactique**

Kit de chimie

***Didactic***

*Chemistry kit*

**Ref :  
107 105**

Français – p 1

English – p 7

Version : 1109

**Kit protection et corrosion du fer**  
***Iron protection and corrosion kit***

# 1 Description

## 1.1 Généralité

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# 1 Description

## 1.1 General

The experiments contained in the iron protection and corrosion kit make it easy to view how iron is attacked and protected. The appearance of colours around the nail can be observed to evaluate progress with the experiment. These experiments are carried out with safe materials.

Firstly, the pupil uses the kit to determine how a nail can be protected using metals. To achieve this, the pupil is provided with two metals that he turns around the nail before they are put into the reactional solution. The appearance of the colour shows which metal protects the nail and which accelerates its oxidation.

A second experiment demonstrates cathodic protection of the metal. The pupil dips two nails into the second reaction solution and makes the same observations as before. Finally, he connects one of the nails to the positive terminal and the other to the negative terminal of a 9 V battery and then observes an immediate colour change.

## 1.2 Composition and description

This kit can thus be used to make 2 complete experiments:


Photo 1

1 – 2 x Agar-Agar 2 g packaging	4 – 1 aluminium wire
2 – 2 x 50 mL reaction solution	5 – 1 copper wire
3 – 2 x Petri dishes	6 – 5 x iron nails

Table 1

## 2 Additional equipment

It is recommended that the following equipment should be obtained to complete these two experiments:

Item	Quantity
Small 100 mL beaker or glass pot	2
Glass stirrer or rod	1
Spatula	1
9V battery	1
Crocodile clip	2
Heating element	1
Pair of vinyl or latex protective gloves (1)	1
Safety goggles (1)	1

Table 2

(1) Personal protection equipment recommended but not compulsory.

## 3 Usage and safety precautions

### 3.1 Warning

The Jeulin Company cannot be held responsible for an accident that occurs after the iron protection and corrosion kit has been used under conditions other than those described in these instructions.

Similarly, the Jeulin Company cannot be held responsible for an accident that occurs because the safety instructions given in these instructions are not respected.

### 3.3 Individual protection

The iron protection and corrosion kit was formulated to minimise risks related to use of the product. Although not compulsory, we recommend that safety goggles shall be worn in addition to the lab coat during handling. It is also recommended that gloves should be worn.

### 3.4 Collective protection

No collective protection equipment is required during storage (safety cabinet, etc.) or during use (hood, etc.), because reagents do not contain any flammable volatile or toxic solvents.

### 3.5 Waste

The iron protection and corrosion kit does not contain any components that create any particular danger for the environment. However, these products must be considered as being chemical laboratory waste and treated as such, in the context of a citizen environmental protection policy. Reagent quantities used in this kit are kept low in order to minimise waste volumes and disposal costs.


## 4 Procedure

### 4.1 Principle of the experiments

The iron protection and corrosion kit demonstrates the attack or protection of iron by means of a solution conducive to corrosion, namely an alkaline medium. Oxidation/reduction couples can be positioned on the redox potentials scale to determine if the metal surrounding the iron will protect it or accelerate its corrosion. This kit can be used for experiments with cathodic protection of iron using a 9 V battery.

#### 4.1.1 Why galvanise metals?

The following couples are involved in these reactions:


It can be seen that copper is less reducing than iron but aluminium is more reducing than iron.

#### 4.1.1.1 Iron alone

A blue colour is observed around the iron nail.

This colour is due to the presence of potassium hexacyanoferrate (III) in the medium that forms a blue complex with  $\text{Fe}^{2+}$  ions.

A pinkish red colour appears at the end of the nail due to the oxidation of iron oxide ( $\text{FeO}$ ). The colour appears due to the presence of cresol red in the medium that changes to pinkish red when the pH exceeds 8.8.



#### 4.1.1.2 Iron surrounded by copper wire

In this case, the blue complex is concentrated around the nail, because iron is more reducing than copper and it will therefore be oxidised by the solution.

As time progresses, a pinkish red colour will be seen around the copper wire.


Ambient oxygen dissolves in the medium at this location and produces  $\text{HO}^-$  ions, which explains this colour change because the pH becomes basic at this location.

#### 4.1.1.3 Iron surrounded by aluminium wire

In this case, there is no blue complex because aluminium is more reducing than iron, therefore it is the aluminium that is oxidised by the solution. On the other hand, a pink colour is seen around the iron. Oxygen is consumed at this location.

As time progresses, a white deposit can be seen around the copper wire, because  $\text{Al}^{3+}$  ions are formed and oxidised by air to create the white  $\text{Al}_2\text{O}_3$  compound.

These experiments demonstrate that a metal can be protected by covering it with another metal more reducing than the first. This is the conclusion of this first experiment.

#### 4.1.2 Principle of cathodic protection

The second experiment demonstrates cathodic protection of iron.

Once the nails are in the reactional medium, they are connected to the terminals of a 9 V battery. Two phenomena can be seen after they have been connected.

1<sup>st</sup> phenomenon:

At the battery – terminal: The nail connected to this terminal is protected against corrosion (surrounded by a pinkish red colour). Water is reduced at this terminal.

2<sup>nd</sup> phenomenon:

At the battery + terminal: The nail is attacked (surrounded by a blue colour), electrons are extracted from the nail. Iron is oxidised at this terminal.

This demonstrates that vehicle bodywork can be protected by connecting it to the battery negative terminal.

## 4.2 Preparation for the manipulation

Place the complete kit and recommended accessories on each counter top.

### 4.2.1 Preparation of the reactional medium

Put the 50 mL solution into a beaker and then place it on a heating plate, taking all necessary precautions. During this time, turn an aluminium wire around one nail and a copper wire around another nail (*photo 2*). Position the nails in a Petri dish (*photo 3*). Take the temperature of the reactional medium. **The temperature of the solution must not exceed 75°C for the remainder of the experiment** (*photo 4*).

Photo 2

Photo 3

Photo 4

### 4.2.2 Preparation of reactions

#### 4.2.2.1 Preparation for demonstration of iron protection and corrosion by the use of metals

When the temperature of the medium is equal to 75°C, withdraw the beaker from the heating element and add 2 g of Agar-Agar (*photo 5*). Mix the solution (*photo 6*) and then pour the entire medium into the Petri dish, distributing it over the entire surface of the dish (*photo 7*). The reaction begins as soon as the mix is added to the nails.

Photo 5

Photo 6

Photo 7

#### 4.2.2.2 Preparation of the cathodic protection

Repeat the preparation described in 4.2.1. During the preparation, place two nails in the second Petri dish and connect the crocodile clips to the nail and to **a single terminal of the 9V battery so as not to discharge the battery** (*photo 8*).

When the temperature of the medium is 75°C, remove the heating element and add 2 g of Agar-Agar (*photo 5*). Mix the solution and then pour the entire medium into the Petri dish, distributing it over the entire surface of the dish (*photo 7*).

Once the mix is distributed, connect the second crocodile clip to the 9 V battery (*photo 9*).

Photo 8


Photo 9

#### 4.2.3 Cleaning

Rinse the beakers and other equipment used as soon as possible to prevent the need for the painstaking cleaning that would be necessary if the remainder of the reactional solution is allowed to stick to the walls of the glass and the tooling. Rinsing with water is all that is necessary.

## 5 After-Sales Service

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your  
disposal from Monday  
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries