Optique

Optics

Réf :
202 001

Français – p 1

English – p 4

Version : 8110

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Optical Bench

1. Composition

1.1 Dans le tube

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1. Composition

1.1 In the tube

One guide bench made of epoxy-coated sheet steel with a silkscreen-printed millimetric scale.

Bench dimensions: U-shaped member with outer dimensions 50 x 32 mm, thickness 3 mm, length 2 m.

Length of the silkscreen-printed scale: 1,800 mm.

1.2 In the accessory box:

- 1 crossbeam made of cadmium-plated steel with 2 rubber shoes and 1 screw for rail mounting.
- 1 crossbeam made of cadmium-plated steel with 2 locking screws and 1 screw for rail mounting.
- 1 lamp on spring shoe fitted with:
 - 220 V - 40 W lamp, B 22 socket, flat filament type.
 - Mains cord
 - Kodatrace rectangle located behind the letter F to avoid direct images of the lamp filament.
 - One lamp cover with a $\varnothing 38$ opening for use when the filament image is to be studied or when another object than the supplied F is used.
- One lens or mirror holder with spring shoe, slewable in any direction (which allows departing from Gaussian conditions, as required) and adjustable in height. The $\varnothing 40$ lenses or the mirror are secured by means of a metal clip also supplied.
- One lens holder similar to the previous one but neither slewable nor adjustable.
- One spring-loaded screen holder.
- One cardboard millimeter-scale screen.
- One set of diaphragms.

2. Additional accessories required

Ordering Code 203 232: Kit comprised of a set of 6 lenses $\varnothing 40$ and one plane mirror $\varnothing 40$, i.e.:

- convergent lens $\varnothing 40$: focal distance $f = + 200$
- convergent lens $\varnothing 40$: focal distance $f = + 100$
- convergent lens $\varnothing 40$: focal distance $f = + 50$
- divergent lens $\varnothing 40$: focal distance $f = - 100$
- divergent lens $\varnothing 40$: focal distance $f = - 200$
- divergent lens $\varnothing 40$: focal distance $f = - 500$
- plane mirror $\varnothing 40$.

3. Set-up

- On both bench ends, mount the 2 cadmium-plated steel crossbeams used to support the bench and for horizontal adjustment.
- Remove the rear cover from the lamp (wire side), remove any cushioning around the bulb and refit the cover.
- Place the lamp on the bench so that the index on the lamp points onto division « 0 » of the millimetric scale.
- According to the experiments to be performed, install the required number of lens or screen holders so that the indices move above the millimetric scale.
- Mount the lenses, diaphragms and the required screen.
- Plug into a 220 V power outlet.

4. Use

The Optical Bench allows the study of reflection and of the plane mirror, the qualitative study of spherical convergent and divergent lenses, and the study of the magnifying glass principle. It is particularly suited for quantitative experiments on lenses.

4.1 Verification of thin lens formulas

The procedure can be as follows:

- Move the lens on the bench in order to give the object-to-lens distance p ten values or so allowing the image to be recovered at an object distance at most equal to the bench length.
- Let p' the lens-to-image distance; read the respective distances p and the image-to-object distances ($p+p'$); on the screen, measure the dimension $A'B'$ of the image of the object F with dimensions $AB = 18 \text{ mm}$.
- Repeat the same steps with a virtual object obtained by means of an auxiliary lens.
- The table of results obtained can be presented as follows:

p	$p + p'$	p'	$\frac{1}{p} + \frac{1}{p'} = k = \frac{1}{f}$	$\frac{p'}{p}$	$\frac{A'B'}{AB}$
-----	----------	------	--	----------------	-------------------

The experiment can be performed with several lenses. Their focal distance can be evaluated by measuring the lens-to-image distance of the lens of an object located at a great distance. Then the following equation is verified:

$$K = \frac{1}{f}$$

4.2 Focimetry

❖ Approximate method:

See above.

❖ **Autocollimation method:**

On the bench fitted with its lamp and with the cardboard screen, install a plane mirror a few centimeters behind the lens whose focal distance is to be measured. Move the lens-mirror assembly until a very sharp image of the object is obtained on the screen, reversed with respect to the object and of the same size, then the object and its image are spaced apart by the distance f .

This allows the verification of the theorem of convergence: successively measure the focal distances focal F_1 and F_2 of two convergent lenses then the focal distance F of the set comprised of the two attached lenses.

❖ **Bessel's method:**

Place the screen towards the bench end at a distance D . Observe that for two positions symmetric with respect to the middle O of the object-to-screen distance and spaced apart by a distance d , a sharp image is obtained on the screen. These positions are verified to correspond to:

$$F = \frac{D^2 - d^2}{4 D}$$

Measurements can be taken for different values of D and recorded in the following table:

D	d'	$D^2 - d^2$	$F = \frac{D^2 - d^2}{4 D}$
-----	------	-------------	-----------------------------

❖ **Silbermann's method:**

This is a particular case of the previous method where $D = 4 f$: in this case, the two lens positions are merged; look for the screen and lens positions yielding a real image of the same size as the object.

5. Compliance with standards

This equipment complies with Standard EN61010-1, Installation Class II, pollution degree 2.

6. After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux