Optique

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Optics

Geometric optics

Réf :
202 012


Français – p 1

English – p 38

Version : 9101

Fibroptonic

Fibroptonic

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
This set of optoelectronic accessories is used to illustrate and study:

- Transmission of signals
- Their transmission by optical fibre
- Their reception.

It is also used to determine the frequency of various vibrations.

The suggested experiments are simple and often spectacular.

Composition :

- 1 transmitter module
- 1 optical fibre (about 5 meters long)
- 1 receiver module
- 1 photoresistor module
- 1 instructions manual including 25 experiments
- 1 storage kit

1 Principle - Description

During experiment sessions, it is possible to have 2 groups of students working simultaneously: for example, one can experiment with the transmitter and receiver whereas the other can perform measurements of periods with the photoresistor module.


Figure 1: the transmitter module


Figure 2: The receiver module


Figure 3: the Photoresistor module

2 Commissioning

❖ Power supply:

Both, transmitter and receiver modules must be powered by a constant DC generator providing between 6 and 12 volts. It is thus possible to power them both with a single voltage generator even if it is educationally better to separate the input circuit from the output circuit by using either 2 voltage generators or 2 9 V batteries.

❖ Transmitter module test:

Simply connect its red and black terminals to the power supply and operate the on/off switch located beside the optical transmitter. When the green LED comes on, an infrared signal having a wavelength of 0.85 micrometers is emitted by the optical transmitter.

❖ Receiver module test:

Simply provide a voltage of between 6 and 12 V between its red and black terminals (it is not necessary to connect the optical fibre). The red LED must be off. As soon as we illuminate the photoreceptor with a torch, the red LED comes on; the receiver has received an optical signal.

❖ Optical fibre:

VERSION No.1

To connect the 2 transmission and reception units with optical fibre:

- Loosen the 2 blue and black plastic screws
- Pass the fibre until its ends touch the emitting and receiving diodes located inside the blue and black ends of the transmitter and receiver units
- Tighten once again the 2 screws without forcing them (figure 4).



Figure 4: installing the optical fibre

VERSION No. 2

To connect the 2 transmission and reception units with optical fibre, simply insert the two ends of the fibre in each of the units.


2.1 How to cut an optical fibre:

- The 2 ends of the fibre have been carefully stripped and scraped. But if, for any reason, we wish to cut it, the ends must be recut carefully.
- Using a cutter cut off 5 to 10 mm of the transparent end of the black protective casing.

- Cut in such a way that only 3 mm of the fibre is outside the casing.
- To increase the transmission factor of the fibre, it is recommended to polish the end with very fine sandpaper. The polishing can also be done with a cloth soaked in mirror.

2.2 Spectral characteristics of the light transmitter and photoreceptor

3 After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

+33 (0)2 32 29 40 50

FIBROPTONIC	PROPAGATION OF AN INFRARED SIGNAL IN AIR	SHEET 01
-------------	--	-------------

1 Purpose

Experimentation on the propagation of infrared signals in air and in transparent mediums

The light transmitter provides an infrared signal of about hundred nanometers of spectral bandwidth, centred at the wavelength of 850 nm. The hatched area on figure 5 corresponds to the spectral emissivity of the light transmitter. The optical fibre of 1 mm diameter transmits this signal without much absorption. The photoreceptor has a wide spectral sensitivity ranging from near UV (300 nm) to infrared (1100 nm), with a maximum at 820 nm (figure 5).


Figure 5 : Spectral characteristics of the light transmitter and photoreceptor

2 Experiment No. 1

2.1 Equipment:

- Transmitter
- Receiver
- 6, 9 or 12 V DC power supply (the optical fibre is not used)

2.2 Set-up and experiment

- Power the transmitter and the receiver (connect the red terminals to the + of the generator and the black ones to -).
- Place the switch located beside the transmitter connector on O (on), the green LED must come on and if we look inside the blue connector we can see a red glow in version No. 1 which corresponds to the visible part of the emitted signal.
- The red LED of the receiver must remain off, as it does not receive any signal.
- Place the 2 connectors facing each other (figure 6). The red LED of the receiver must come on. The infrared signal has crossed the space between the light transmitter and the photoreceptor. Beyond a gap of 20 mm (when powered by 12 volts) and the receiver's LED goes off.


Figure 6: Propagation of an infrared signal in air

2.3 Additional experiments

- Inserting a piece of paper interrupts the propagation of the signal.
- Inserting a glass slide between the light transmitter and the photoreceptor does not stop the propagation (try with various materials).

2.4 Conclusion

A signal in the near infrared propagates through air and glass (caution, certain mediums that are transparent in the visible range are not so in infrared at certain wavelengths, glass for example).

FIBROPTONIC	OPTICAL AND GEOMETRICAL PROPERTIES OF AN INFRARED BEAM	SHEET 02
-------------	---	-------------

1 Purpose

To study the directivity of an infrared beam and its optical properties. Acoustical detection.

2 Experiment no. 2

Rectilinear propagation of an infrared beam

2.1 Equipment:

- Transmitter
- Receiver
- 6 V or 12 V DC power supply
- Signal generator (delivering a sinusoidal voltage with a frequency between 1000 and 2000 Hz).
- The optical fibre is not used.

2.2 Set-up:

- Place the transmitter and receiver face to face (figure 7)
- Power the transmitter and receiver with 6 V or 12 V DC (in this case 12 V is preferable)
- Place the switch located beside the transmitter connector to On (O), the green LED must come on.
- Place the second switch, which controls the microphone to EXT (see figure 7)
- Connect the signal generator between the yellow and black terminals (black to the ground of the signal generator). The amplitude of the sine-wave signal is of the order


Figure 7: Set-up indicating the rectilinear propagation of an infrared beam

2.3 Experiment:

After adjusting the potentiometer of the receiver unit to the maximum gain, try to align the light transmitter and photoreceptor. When they are perfectly aligned, an audible signal is heard even if the distance between the blue and black connectors is 10 to 15 cm. On the other hand, the red LED comes on only if the distance is less than 20 mm.

2.4 Conclusion:

The critical alignment of the transmitter and receiver indicates the significant directivity of the infrared beam emitted. This experiment also illustrates the transformation of an electrical signal into an infrared signal by the transmitter unit, the propagation of an infrared signal in air and the opposite operation, which transforms an infrared signal into an audible signal emitted by the loudspeaker on the receiver.

3 Experiment no. 3

Reflection of an infrared beam. Acoustical detection.

Repeat the previous experiment but without aligning the transmitter and receiver. The 2 units are separated by a few centimetres. We then insert a glass slide or a small plane mirror, which we adjust until an audible signal is heard (figure 8). We can thus illustrate Snell's law for the reflection of an electromagnetic wave, which does not belong to the visible region.


Figure 8: Reflection of an infrared beam

4 Experiment no. 4

Refraction of an infrared beam. Acoustical detection.

In relation to the previous experiment, simply replace the plane mirror with a prism or a half cylinder in Plexiglas or glass. The prism's apex angle must be 30 or 45° to prevent the occurrence of total reflection of the beam too easily. Rotate the prism with small movements around a vertical axis until an audible signal is heard (figure 9). Although it may not be possible to carry out accurate measurements of angles, we can however show that the refraction with angular incidence corresponds to a deviation of the beam during a change in medium.


Figure 9: Refraction of an infrared beam

FIBROPTONIC	PROPAGATION OF AN ELECTROMAGNETIC SIGNAL IN AN OPTICAL FIBER	SHEET 03
-------------	--	-------------

1 Purpose

To convert signals of different nature, carry information, display.

The transmitter emits an infrared signal of 0.85 micrometers wavelength, but the initial information that has created this signal can be of a different nature: electrical, acoustic, light, etc. At the receiver end, the infrared signal is converted into an audible acoustic signal or an electrical signal, which can be viewed on an oscilloscope. The spectral response of the fibre is quite wide. The fibre also transmits radiations located in the visible region.

2 Experiment no. 5

Propagation through optical fibre. Detection of a light signal.

2.1 Equipment

- Fibre
- Receiver
- 6, 9 or 12 V DC power supply
- Pocket torch (or helium-neon laser)

2.2 Set-up and experiment

- Power the receiver (red and black terminals)
- Connect the fibre to the receiver (the transmitter unit is not used)
- Illuminate the end of the fibre with a pocket torch (figure 10).


Figure 10: Detection of a light signal by the receiver

2.3 Observations and conclusion:

- The red LED of the receiver must come on.
- The optical fibre has transmitted the light information to the photoreceptor. The light signal whose maximum emission is at about 0.57 micrometers has been detected. Thus the fibre does not transmit only infrared signals.

3 Experiment no. 6

Optical transmission of logical information.

3.1 Equipment:

- Transmitter
- Receiver
- Optical fibre connected to the transmitter and receiver
- 6, 9 or 12 V DC power supply
- Oscilloscope

3.2 Set-up:

- Power the transmitter and the receiver (red terminal to +, black terminal to -)
- Connect the fibre properly to the 2 units (see figure 4)
- Place the microphone switch on the EXT position.

3.3 Experiment:

- Toggle the other switch alternately (On – Off positions).
- When the switch is on ON, the green LED comes on, an infrared signal is sent through the fibre. This signal is detected by the photoreceptor, the red LED comes on.
- If we disconnect the optical fibre, the infrared signal cannot be detected: the red LED goes off.
- These 2 binary logic states (state 1: LED is on / state 0: LED is off) can be viewed on the oscilloscope. For this, connect the yellow and black terminals of the receiver to the oscilloscope (yellow: channel A and black to the ground of the device). The selector switch of the oscilloscope must be on the DC position. When the LEDs are on, the amplitude of the signal on the oscilloscope reaches 0.4 to 0.7 V (depends on the supply voltage of the transmitter and receiver and on the length of the optical fibre).
- Figure 11 shows the appearance of a signal detected on the oscilloscope when the switch of the transmitter unit is quickly toggled (On – Off).


Figure 11: Oscillogram corresponding to the rapid switching to change from logic state 0 to state 1.

FIBROPTONIC	OPTICAL FIBER TELEPHONY	SHEET 04
-------------	-------------------------	-------------

1 Purpose

To discover the principle of optical fibre telephony (transmission through wire).

2 Experiment no. 7

2.1 Equipment:

- Transmitter
- Receiver
- Fibre
- 6, 9 or 12 V DC power supply

2.2 Experiment

- Unwind the optical fibre as much as possible: the transmission and reception units are separated by 5 m. Place the microphone switch on MIC and the transmitter switch on On. The 2 green and red LEDs must come on.
- A student speaks in the transmitter microphone (figure 12)
- The receiver is placed on the table amongst the students. The sound can be perfectly heard if the sound level (which we adjust using the potentiometer located on the receiver unit) is sufficiently high.


Figure 12: Optical fibre telephone

2.3 Conclusion:

The transmitted audible signal modulates the infrared signal, which propagates through the fibre. At the receiver end, the signal is “decoded” and the audible signal reconstructed.

2.4 Alternative:

The student whistles regularly, the signal received can be viewed by connecting an oscilloscope to the yellow and black terminals of the receiver.

This experiment therefore illustrates: the transmission, propagation, reception and modulation of an infrared signal.

FIBROPTONIC	OPTICAL TRANSMISSION OF A PERIODIC SIGNAL	SHEET 05
-------------	---	----------

1 Purpose

To illustrate the transmission of an electrical signal through an optical fibre.

2 Experiment n° 8

2.1 Equipment:

- Transmitter
- Receiver
- Fibre
- 6, 9 or 12 V power supply
- Sinusoidal or rectangular signal generator
- Oscilloscope

2.2 Set-up and experiment

- Adjust the signal generator so that it gives sinusoidal or rectangular signals with peak-to-peak amplitudes of between 50 and 300 mV, and frequencies of 1000 to 2000 Hz.
- Connect the signal generator to the yellow and black terminals of the receiver (black to ground). Place the microphone switch on the EXT position and the other switch on On. The green and red LEDs must come on (signal received).
- The sound level is adjusted with the receiver's potentiometer. It is possible to set it on zero: there is no sound in this case.
- To view the signal received, the oscilloscope must be connected to the yellow and black terminals of the receiver unit (black to the ground of the oscilloscope). Refer to figure 13.


figure 13 : Transmission of a periodic signal

2.3 Observations:

- The detected signal is also periodic. It may be slightly deformed. The deformation depends on the nature of the signal (rectangular or sinusoidal), on its amplitude, its frequency (between 10 Hz and 50 kHz) and on the value of the supply voltage of the transmitter and receiver.
- An electromagnetic signal (infrared) modulated by a periodic signal of a completely different nature passes through the optical fibre.
- By rapidly toggling the on/off switch located on the transmitter unit, we can create bursts of easily audible sine waves.

FIBROPTONIC	DETECTION AND TRANSMISSION THROUGH AN OPTICAL FIBER	SHEET 06
-------------	--	-------------

1 Purpose

To use the optical fibre to transmit a signal detected by a sensor.

2 Experiment no. 9

2.1 Equipment:

- Transmitter
- Receiver
- Photoresistor module
- 6 V (or 9 V) DC power supply
- Optical fibre
- Pocket torch
- Oscilloscope

2.2 Set-up:

- Connect the optical fibre to the transmitter and receiver (figure 4)
- Place the microphone switch located on the transmitter unit on EXT
- Connect the photoresistor (yellow and black terminals) in series with the 6 V DC power supply to the red and black terminals of the transmitter (figure 14).
- Power the receiver with 6 V.


Figure 14: Set-up to detect and transmit a variation in illumination.

2.3 Experiment

- When we toggle the receiver's switch to On, the green LED comes on only if we illuminate the photoresistor sufficiently. An infrared signal passes through the optical fibre and the red LED also comes on.
- When the photoresistor is lighted with a pocket torch, the signal received when viewed on the oscilloscope has an amplitude of 0.3 to 0.6 V for a supply voltage of 6 V.
- If we light the photoresistor with a pocket torch in an oscillating motion, the green and red LEDs go on and off alternately. The signal observed on the oscilloscope is almost periodic.

Various applications

3 Experiment no. 10

Remote alarm triggered by light

This set-up could serve as an alarm clock sensitive to daylight.

3.1 Equipment:

- Transmitter
- Receiver
- Optical fibre
- 6, 9 or 12 V Power supply
- Pocket torch
- Photoresistor module
- LF sine-wave signal generator (amplitude 0.5 V and frequency 1 to 2 kHz).

3.2 Set-up and experiment

- Repeat the set-up of figure 14 by connecting the signal generator to the black and yellow terminals of the transmitter (black to the ground of the LFG). One of the 2 devices (power supply or LFG) must have a floating ground because of the set-up on the transmitter side.
- Place the 2 switches on EXT and On (transmitter unit)
- Set the sound level to the maximum (potentiometer of the receiver)
- Illuminate the photoresistor with a pocket torch and the alarm is triggered (the sound level is greater and the triggering threshold lower if a 12 V DC power supply is used). The torch can be moved quickly from left to right and right to left alternately over the photoresistor.

4 Experiment no. 11

Examining a spectrum with an optical fibre.

4.1 Equipment:

- Slide projector
- Slit of 1 mm cut out into a 5x5 cm cardboard.
- A diffraction grating of 140 lines/mm or a glass prism
- Receiver
- Fibre
- 6, 9 or 12 V DC power supply (the transmitter is not used)
- Oscilloscope.

4.2 Set-up and experiment:

- Power the receiver
- Place the vertical slit in the slide projector
- Focus the image of the slit on the cardboard screen located at 1.5 m (approximately) from the projector.
- Fix the diffraction grating (vertical lines) against the projector's lens. Many coloured bands are visible; these are white light spectra in different orders of the diffraction grating (single spectrum with a prism)
- Move the screen in such a way that many spectra are half cut by the upper edge of the screen; examine the 1st order spectrum by moving the end of optical fibre tangentially to the screen.
- To view the received signal, connect an oscilloscope (switch on DC) to the yellow and black terminals of the receiver.


Figure 15: Examining a white light spectrum with an optical fibre

4.3 Observations

- The spectral response of the assembly: optical fibre + photoreceptor is weak in blue-violet; it is maximum in red and infrared.
- This experiment shows that there exists radiation outside the visible region, beyond the red: the infrared.
- The amplitude of the signal received is dependent on the optical quality and the state of the surface of the end of the fibre. It must be perfectly polished and transparent. The amplitude of the signal also depends on the lighting power of the white light source and on the distance between the projector and optical fibre.
- To increase the signal we can disconnect the fibre, which absorbs a small fraction of the electromagnetic radiation and then directly examine the spectrum with a single receiver module
- On the other hand, with a diffraction grating as a dispersive system, a large part of the light emitted is lost in the different orders of the fibre which are not examined. Furthermore, when we examine the 1st order infrared (wavelength greater than 800 nanometers), it superimposes itself on the 2nd order violet-blue radiation (wavelength greater than 400 nm). Thus, it is better to use a prism, which gives a single very bright spectrum and avoid the problems of overlapping of orders, but it has the disadvantage of further dispersing the blue instead of the red.

FIBROPTONIC	MEASUREMENT OF PERIODS	SHEET 07
-------------	------------------------	-------------

1 Purpose

To determine the period of a periodic system using three fibroptonic modules and an oscilloscope.

2 Experiment no. 12

2.1 Equipment:

- Transmitter
- Receiver
- Optical fibre connected to the previous units
- Photoresistor
- 6, 9 or 12 V DC power supply
- Oscilloscope
- Light source: pocket torch, laser, etc.
- Various vibrating systems: vibrator, etc.

2.2 Set-up:

- Place the receiver switches on On and EXT
- Power the transmitter and receiver with DC voltage; insert the photoresistor into the emitter's power supply circuit between the black terminals of the transmitter (figure 16).
- If we illuminate the photoresistor, the 2 LEDs must come on.
- Connect the oscilloscope between the yellow and black terminals of the receiver. The oscilloscope must be able to detect the variable AC and DC components (selector switch on DC).
- For a supply voltage of 6 V, the signal attains an amplitude of 0.5 V when the photoresistor is illuminated by a pocket torch.


Figure 16: Set-up to determine the period of low frequency vibratory mechanical phenomena.

2.3 Experiment:

- Illuminate the photoresistor with a torch or better still with a laser beam.
- Intercept the light beam with the vibrator whose period needs to be determined (saw blade, spring, etc.)
- Move the vibrator from its equilibrium position and observe the signal received on the oscilloscope (figure 17)
- The vibrator cuts the beam twice every period. Thus the period of the vibrator is double of the signal period viewed on the oscilloscope.
- If in the resting position, the blade of the vibrator does not cut the laser beam, the appearance of the signal can be modified (figure 18).


Figure 17: At every half period a vibrating blade intercepts the light beam.


Figure 18: The vibrating blade intercepts the light beam twice every period, but at rest, the vibrator does not cut the beam.

2.4 Note:

To determine the frequency delivered by a flash stroboscope or a rotating disk with slits, simply project the light emitted by the stroboscope on the photoresistor and view the signal on an oscilloscope.

2.5 Using optocouplers to electrically separate 2 circuits:

In electrical set-ups, we are often faced with problems of grounding (0 volts on certain points of the circuit). The use of optoelectronic components solves these problems by ensuring an optical communication between 2 circuits, which cannot have common electrical points.

In this set of accessories, the transmitter and receiver are not connected electrically if each one of them is powered by a 9-volt battery.

Before acquiring these components, disconnecting transformers were used.

FIBROPTONIC	DETECTION OF VARIATIONS IN LIGHTING	SHEET 08
-------------	-------------------------------------	----------

1 Purpose

To detect a variation in the light intensity using a photoresistor.

2 Precautions for use

The resistance of the photoresistor reduces when it is illuminated. In sunlight, the resistance is a dozen ohms; in darkness the resistance can exceed thousands of ohms. The photoresistor R is installed in series with the potentiometer RV of 2.2 kilo ohms and a resistance R0 of several hundred ohms. If U is the supply voltage on the terminals of the assembly, the electrical power for the photoresistor is:

$$P = \frac{R \cdot U^2}{(R + R_0 + R_V)^2}$$

This power is maximum for RV = 0 and R = R0 thus Pmax = U²/4 R0. As the maximum power indicated by the manufacturer is 0.3 watts, we deduce that U must be less than (1.2 R0)^{1/2} that is 12 V for R0 = 120 ohms.

3 Light sources whose brightness varies periodically. Measurement of periods

This part can be the subject of experiment sessions. The experiments suggested are easy to implement and not very common; the results are immediate.

This part can be the subject of experiment sessions. The experiments suggested are easy to implement and not very common; the results are immediate.

3.1 Equipment:

- Photoresistor module
- 6, 9 or 12 V DC power supply
- Oscilloscope
- 40 or 60 W light bulb operating on 220 V mains supply
- Pocket torch (battery powered)
- Neon indicator (stair opening type) operating on mains supply
- Fluorescent tube light (classroom lighting, for example)
- Mechanical stroboscope or a DC motor fitted with a cardboard disk pierced with holes or slits.

or cinecolor Part no. 201 019

3.2 Set-up:

- Power the photoresistor with 6 or 9 V DC between the 2 black and green outer terminals. The 3 elements of the module being symmetrical dipoles, the connecting

direction does not matter. We will therefore connect the black terminal to the negative and the green to the positive of the power supply.

- Connect the oscilloscope to the terminals of the photoresistor: black terminal to the ground of the oscilloscope and the yellow terminal to channel A (figure 19).


Figure 19: Connecting the photoresistor module

4 Experiment no.13

Fluctuations in the brightness of a light bulb operating on the mains supply.

- Illuminate the photoresistor with a light bulb held at about ten centimetres.
- Set the oscilloscope in the AC position (only AC component) and amplify.
- We observe a sine curve, whose amplitude can be adjusted with the potentiometer located on the unit (the amplitude depends on the distance between the bulb and the photoresistor).
- The period determined on the oscilloscope is 10 milliseconds (thus a frequency of 100 Hz, i.e. double of the mains frequency because the bulb comes on at every positive or negative half-wave of the AC voltage of the mains).
- If we repeat the same experiment with a pocket torch (operating on batteries), there are no fluctuations.

5 Experiment no. 14

Periodic variations in brightness of a neon indicator operating on mains

- Repeat the previous experiment and place a small neon lamp (orange light) powered by mains AC voltage, very close to the photoresistor.
- Adjust the amplitude of the signal observed on the oscilloscope with the potentiometer (see figure 20). Determine the period, that is 20 ms.


Figure 20: Signal observed on the oscilloscope when the photoresistor is lit with a neon indicator powered by mains current. The shape of the signal depends on the experimental conditions.

6 Experiment no. 15

Variations in the brightness of a fluorescent lamp powered by mains current.

- Repeat the previous experiment and illuminate the photoresistor with a single fluorescent lamp (this gives a white light, its spectrum contains certain mercury lines during transmission).
- Figure 21 shows the shape of the signal detected on the oscilloscope. Note the asymmetry of the signal between the lighting and quenching phase of the lamp.


Figure 21: Signal corresponding to periodic fluctuations in brightness of a fluorescent lamp powered by mains current. The shape of the signal depends on the experimental conditions.

7 Experiment no. 16

Determining the rotational speed of the mechanical stroboscope disk and frequency calibration.

- The oscilloscope is in the AC position.
- Illuminate the photoresistor with a pocket torch or a DC powered light source.
- Intercept the beam, which lights the photoresistor with the rotating disk of a mechanical stroboscope.
- We observe a non-sinusoidal alternating signal whose appearance depends on the rotational speed of the disk, on the number of slits or holes, on the width of these openings, etc.
- Even so the signal is periodic. If the disk has only one opening, the frequency measured corresponds to the rotational speed of the motor in revolutions per second. If the disk has 4 evenly spaced out identical openings, the frequency of the flashes by the stroboscope is equal to 4 times the rotational speed of the motor.
- We can thus use this experiment to determine the rotational speed of a motor. Simply fit it with a disk with one opening.

FIBROPTONIC	MECHANICAL VIBRATIONS AND OSCILLATIONS	SHEET 09
-------------	---	-------------

1 Purpose

To determine the period of a vibrator or mechanical oscillator using a photoresistor module.
To use a photoresistor to study periodic or pseudo-periodic mechanical movements.

THE EXPERIMENTS DESCRIBED IN THIS PART ARE EASY TO PERFORM.

1.1 Equipment:

- Photoresistor module
- 6, 9 or 12 V DC power supply
- Oscilloscope
- Helium-neon laser
- Small piece of mirror measuring a few square millimetres on a adhesive foil (mirror splinters are very convenient).
- Miscellaneous:
 - . Tuning fork
 - . Melde's vibrator
 - . Saw blade
 - . Spring
 - . Loudspeaker of about 10 cm diameter (recovery).
 - . Simple pendulum consisting of a ball attached to a thread
 - . LFG to drive the loudspeaker.

1.2 Set-up:

- Power the photoresistor with 6, 9 or 12 V DC (black terminals connected to the minus of the power supply and the green to the positive)
- Connect the oscilloscope between the yellow and black terminals of the photoresistor (black to the ground and yellow to channel A).

2 Experiment no. 17

Determining the period of vibration of a tuning fork

- Fix the small adhesive mirror on one of the branches of the tuning fork (not useful if the metal is fairly reflecting)
- Project the laser beam on this mirror and intercept the reflected beam with the photoresistor (figure 22). Strike the tuning fork.
- Adjust the amplitude of the signal detected by the oscilloscope with the potentiometer on the photoresistor module (oscilloscope selector switch on AC).
- Strike the tuning fork again and determine its period on the oscilloscope. Deduce the frequency (440 Hz for "A 3" tuning fork). Note that the signal transmitted by a tuning fork that is not sustained is pseudo-periodic (and not periodic).

3 Experiment no. 18

Determining the period of a mechanical vibrator

- Repeat the previous set-up and replace the tuning fork with Melde's vibrator or a saw blade; the small adhesive mirror will be fixed to the vibrator.
- We thus show that the frequency of a Melde's vibrator is 100 Hz when it is powered by the mains AC current (220 V – 50 Hz).
- We can also check that the frequency of the fundamental mode (the one that is generally driven) varies as the inverse square of the length of the vibrating blade. In this case the excitation is not sustained.

4 Experiment no. 19

Vibrations of a loudspeaker diaphragm

- Repeat the set-up of figure 22 and replace the tuning fork by a loudspeaker whose visible membrane is of a large diameter (for example. 10 cm). Fix the small adhesive mirror mid way between the centre and the edge of the membrane.
- Drive the loudspeaker with a low frequency generator, which gives a sinusoidal signal. Adjust the frequency so that the amplitude of vibrations of the adhesive mirror is significant.
- The signal detected by the photoresistor module and viewed on the oscilloscope is periodic. Measure its period and compare it with the indications of the low frequency generator.
- Use another loudspeaker fitted with a mirror, which also intercepts the laser beam. It is now possible to obtain Lissajous figures by carefully combining the 2 periodic functions. Another LF generator will be required.

5 Experiment no. 20

Oscillations of a spring. Measuring the period.

- Project the laser on the photoresistor and intercept the beam with the small mass hanging at the end of the spring (figure 23).


Figure 23: Set-up to determine the period of a spring.

- Adjust the amplitude of the signal detected by the oscilloscope using the potentiometer on the photoresistor module (oscilloscope in AC).
- Make the spring oscillate.

- Depending on the relative position of the laser and the mass at rest, we observe 1 or 2 interruptions of the beam per period of the spring (refer to figures 17 and 18).
- With this experiment we can study the influence of the mass attached to the spring on the period. To study the influence of the spring constant, the spring must be changed.
- We can also intercept the laser beam by the turns; in this case it is necessary to know the number of turns, which intercept the beam.

Alternative version: replace the spring with a vibrating wire which cuts the laser beam and study the action of the tension of the wire on its period of vibration.

6 Experiment no. 21

Determining the period of a pendulum

- The pendulum consists of a small mass attached to a thread with a length of a few centimetres.
- In the set-up of figure 23, we replace the spring with the pendulum; the mass when at rest intercepts the laser beam. Thus the signal detected by the oscilloscope has 2 minima per period.
- This experiment helps verify the relation between the period of a simple pendulum and its length (T varies with $L^{1/2}$).
- If the oscilloscope has a high-performance time base: 0.1 s/cm or even 0.5 s/cm or 1 s/cm, it is possible to determine much longer periods and study all sorts of pendulums.

FIBROPTONIC	SPECTROSCOPY	SHEET 10
-------------	--------------	-------------

1 Purpose

To use the properties of the photoresistor to examine a white light spectrum.

2 Experiment no. 22

2.1 Equipment:

- Slide projector
- 1mm wide slit cutout in 5x5 cm size cardboard.
- Diffraction grating of 140lines/mm
- Photoresistor module
- White screen (of cardboard, for example)
- Multimeter, which can be used as an ohmmeter (to measure the resistance)

2.2 Set-up and experiment:

- Place the vertical slit in front of the slide projector and adjust the focus so that a clear image of the slit is obtained on the screen located at about 1 meter from the projector.
- Fix the diffraction grating well against projector (lines from the diffraction grating must be parallel to the slit).
- Intercept the 1st order spectrum using the photoresistor (figure 24). Attach a small slit cut in a black adhesive tape (width: 0.5 to 1 mm depending on the lighting power of the source and the distance between the projector and photoresistor) on the photoresistive cell.
- Connect the ohmmeter to the black and yellow terminals of the photoresistor.


Figure 24: Examining a white light spectrum with a photoresistor.

2.3 Observations and results:

- We observe that the resistance of the photoresistive cell varies when it is moved in the 1st order spectrum. In blue-violet, the resistance is high because the photocell is less sensitive in this wavelength range and the light emitted by the projector's halogen lamp has very little blue-violet radiation. The resistance is minimal in yellow and orange (figure 25). Note the sensitivity of the photoresistor to near infrared radiation (760 to 800 nm).
- The curve no. 25 is a result of the combination of the spectral response of the photoresistor and the spectral emissivity of the projector's white light lamp. The inverse of the resistance is on the ordinate and the wavelength on the abscissa.
- To establish such a curve, the 1st order spectrum obtained with a diffraction grating must be calibrated using the dispersion relation of diffraction gratings working under normal incidence, that is:

$$\sin i = n.k.\lambda$$

With λ as the wavelength, k the order (in this case 1), n the number of lines 140mm⁻¹ and i the angle between order 0 and order 1 for the selected wavelength (see figure 24). We can also use the helium-neon laser to calibrate.
- Note: beyond 800 nm in the 1st order, there is superposition with the 2nd order radiation of wavelength greater than 400 nm.


Figure 25: Curve resulting from the combination of the spectral response of the photoresistor and the emissivity of the source of white light.

FIBROPTONIC	PHOTOMETRY	SHEET 11
-------------	------------	-------------

1 Purpose

Use a photoreceptor to measure the illumination.

2 Experiment no. 23

The purpose of the experiment described here is to specify the behaviour of the resistance of the photoresistive "cell" in function of the illumination. It gives simple and relatively accurate measurements.

2.1 Equipment:

- Photoresistor module
- Multimeter used as an ohmmeter
- Light source: 40 W to 100 W light bulb

2.2 Set-up and experiment:

- Connect the ohmmeter to the black and yellow terminals of the photoresistor.
- Illuminate the photoresistor with the bulb; D is the distance between the bulb and the photoresistor (figure 26).
- Vary D by every 10 or 20 cm between 0.20 m and 2.5 m and measure the value of the resistance R of the photoresistor.


Figure 26: Measuring the resistance R in function of the illumination.

2.3 Results:

Figure 27 shows the shape of the resistance R in function of the distance D for 2 bulbs of different wattages (40 and 75 W). We observe that R is a linear function of D. Since the illumination on the photoresistor varies with $1/D^2$, we deduce that R varies with the inverse square root of the illumination thus with $1/E^{0.5}$.


Figure 27: Curve representing the resistance R of the photoresistor in function of the distance D from the light source (A: 40 W light bulb, B: 75 W bulb).

2.4 Note:

- We obtain similar results if the light sources are fitted with colour filters.
- The wattage of the bulbs (for example 40 W) is an electric power and not a lighting power.
- It is also possible to work at a given distance D with bulbs of different electric power (25 - 40 - 60 - 75 - 100 - 150 W). In this case we measure the resistance R in function of the electric power of the bulb and the flux ϕ emitted.
- While taking a photo, the quantity of light that exposes the film is adjusted by the diaphragm; it varies with the area of the diaphragm thus with the square of the aperture. If the quantity of light is measured using a photoresistor, we deduce that its resistance is proportional to the value of the diaphragm (a reminder of these values: 16 - 11 - 8 - 5.6 - 4 - 2.8 - 1.4). This property is used in photography.

FIBROPTONIC	COLORIMETRY	SHEET 12
-------------	-------------	-------------

1 Purpose

To approximately determine the concentration of a coloured solution from a standard solution of known concentration.

2 Experiment no. 24

2.1 Equipment:

- Photoresistor module
- Ohmmeter
- Directive source of white light or laser
- Set of identical glass tanks
- Different solutions of potassium permanganate of concentrations between 10^{-4} and 10^{-1} mol/litre.

2.2 Set-up and experiment:

- Connect the ohmmeter to the black and yellow terminals of the photoresistor.
- Illuminate the photoresistor with the light beam
- Place a glass tank filled with a decimolar potassium permanganate solution (control solution) in the beam; measure the resistance R of the photoresistor (figure 28).
- Restart the operation with the other control solutions of up to 10^{-4} mol/l. Each time, check that the tank is perpendicular to the beam (for this, observe the laser's reflected spot) and that the distance from the source to the photoresistor is unchanged.
- Plot the curve of the resistance R in function of the concentration C .
- Pour the solution whose concentration needs to be determined into the glass tank, measure R , deduce C from the curve established earlier.


Figure 28: Determining the concentration of a coloured solution by absorption.

2.3 Note:

- The accuracy of the concentration is 5%.
- We also obtain good results with copper sulphate II solutions.
- If we rotate the tank around a vertical axis, the optical thickness as well as the absorption increases (Lambert's cosine law).
- If we use a laser, the distance between the laser and the photoresistor must be sufficient so that it illuminates a significant part of the photoresistor.
- Be careful of ambient light interference.

FIBROPTONIC	CHEMICAL KINETICS	SHEET 13
-------------	-------------------	-------------

1 Purpose

To study the kinetics of the chemical reaction between potassium iodide and sodium peroxodisulfate solutions using photometric measurements.

2 Experiment no. 25

2.1 Equipment:

- Photoresistor module
- Ohmmeter
- Directive white light source or laser
- Set of identical glass tanks
- Potassium iodide solution of 0.5 mol/litre
- Sodium peroxodisulfate solution of 0.05 mol/litre.
- Stopwatch

2.2 Set-up:

- Mix 250 ml of potassium iodide with 250 ml of sodium peroxodisulfate, the solution becomes yellow and then slowly turns brown with the formation of iodine (I_3 ions). Wait for about one hour (this depends on the ambient temperature). The solution obtained does not change anymore. We now prepare different control solutions of known concentration by dilutions: $a - 0.8 a - 0.6 a - 0.4 a - 0.2 a - 0.1 a - 0.05 a - 0.025 a$ etc. (a is the concentration of I_3 ions in the solution obtained earlier, that is 0.025 mol/litre. The volume of each solution is 50 ml.
- Connect the ohmmeter to the black and yellow terminals of the photoresistor
- Illuminate the photoresistor by the light beam (adjust the distance between the source and the photoresistor in such a way that the lighted area is not too small).
- Introduce the tank containing the control solution in the beam, for example one with a concentration of: $a = 0.025$ mol/l (refer to figure 28); measure the value of the resistance R of the photoresistor.
- Restart the operation with each of the control solutions by checking that there is no change in the relative arrangement of the beam, the tank or the photoresistor.
- Plot the curve representing R in function of the iodine concentration.

2.3 Experiment:

- Take a glass tank and pour 25 ml of peroxodisulfate solution in it.
- Prepare 25 ml of potassium iodide in a beaker.
- Place the tank filled with the solution under a light beam; reset the stopwatch.
- Quickly pour the contents of the beaker into the tank and start the stopwatch.
- As the solution becomes brown, the resistance read on the ohmmeter increases. Note down the times when the resistance passes from one of the values obtained in the absorption curve of the control solutions earlier.
- Thus, at these instances the concentration of iodine is well known. Now just plot the curve of the iodine concentration in function of time. As the iodine ions are initially in excess, the final iodine concentration tends towards a.

2.4 Note

- The equation of the reaction is:

- The temperature is a kinetic factor: the lower it is, the slower the change (try between 10° and 40°C), even if the critical value "a" is unchanged.
- We can also study the influence of the initial concentration of iodine or peroxodisulfate ions.
- This set-up also helps study kinetics of the reaction between hydrochloric acid and the thiosulphate ions.

or the action of oxalate ions on a potassium permanganate

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux