Optique

Optics

Réf :
202 171

Français – p 1

English – p 22

Version : 8111

Mini labo. d'optique

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1. Description et principe

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
0 825 563 563*

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1 Principle and description

The Jeulin Optics Mini Lab is a complete set of optical equipment dedicated to experimentation at the college level. It was planned and designed with an objective to enable easy experimentation by young students (light, small dimensions and efficient attachment of components by magnetization), without the need to darken the classroom.

It consists of:

- A black box with a metal bottom used to fit optical accessories installed on magnetised supports
- A white point source installed on a magnetic support
- 5 magnetised slide supports made of black plastic
- A white loop installed on a magnetised support
- A small white ball installed on a rod
- A small plastic file containing 19 small plates of 50 x 50 mm:
 - 2 white opaque screens
 - 2 black opaque screens
 - 2 pink and yellow fluorescent screens
 - 2 red and green filters
 - 4 translucent slides made of polypropylene
 - 3 single slits centred on black and opaque slides
 - 2 black opaque screens with a square and round opening at the centre
 - 2 transparent slides with black opaque square and round at the centre.

2. Technical characteristics

Black box:

- Dimensions: 320 x 170 x 90 mm
- An observation window fixed at one of the ends
- One of the sides may be opened or closed as required, enabling the observation of optics experiments over the entire length of the box from the side
- A removable cover providing the student with easy access from the top to place optical accessories.

Light source:

- Box dimensions: 60 x 30 x 50 mm
- LED position: 30 mm for perfect optical alignment with the accessories
- Magnetic soleplate
- On/Off slide switch
- Power supply: 1 6F22 9 V battery (supplied)
- Source available separately under the part number 202 170.

Accessories:

- Large white ball: Ø 25 mm approximately on magnetic support
- Small white ball: Ø 10 mm approximately installed on a rod
- Small plates: 50 x 50 mm.

3. Some examples of observation

To prepare a non-exhaustive list of observations that can be made using the optics mini lab, let us imagine some actual situations, all taken as examples from an amateur astronomer's concerns, and suggest experiments that can be performed to experimentally check the possible answers.

Why this bias for astronomy? The main thing most probably lies in the interest that the students have in this type of questions, and to the direct observation that can be suggested to them to support these discoveries affirmed in class.

3.1 Question 1 : Under what conditions can the Moon be seen ?

Research: Why do we see an object?


3.1.1 Experiment 1 : without light, vision impossible

Objective: Demonstrate two types of light sources: primary sources and diffusing objects.

Equipment: single source in the box closed by its cover.

Description: first, the observer looks inside the box, with the light turned off, then the source is turned on and placed in different positions.

Observation:


3.1.2 Experiment 2 : the eye must receive the light emitted by the object

Objective: Understand that to see an object, the eye must receive light from it.

Equipment: source, black opaque screen pierced with an opening.

Description: the source illuminates the pierced screen; the student places his eye on the other side, near the opening of the lab.

Observation:


3.1.3 Experiment 3 : an obstacle can cover the source

Objective: Understand that to see an object, the eye must receive light from it. Demonstrate the fact that a light source defines two zones: an illuminated zone from where the observer sees the source and a shadow zone from where the observer does not see the source.

Equipment: source, transparent slide with a black round (or a square) at the centre or a ball on magnetised support, translucent screen (is used to dim the source if the luminosity is too bright for the eye).

Description: Previous experiment “in negative”: when the observer reports that he cannot see the source, the monitoring student sees the shadow of the obstacle stand out against the open eye of his friend. The observer cannot see the source when the eye, the obstacle and the source are aligned.

Observation: The round (or square) shadow of the obstacle stands out against the observer's face, and particularly against his eye when he says that he cannot see the source. A monitoring student looking at the observer, as in the previous experiment, will confirm this. Replacing the slide by an opaque screen of any colour will confirm that an object can be seen if the light diffused by it reaches the eye.

3.1.4 Conclusion:

Certain objects are visible because they emit light (the Sun, for example) and others, most in fact, can be seen only because they diffuse the light that illuminates them (like the Moon acting as a mirror).

We see the Moon, above the horizon, when our eyes can receive light from it.

3.2 Question 2 : Why do we find our way more easily on a full Moon Night ? Is there Earth light on the Moon ?

Research: Can an object illuminate without producing light?

3.2.1 Experiment 1 / : A white opaque object is diffusing

Objective: Define the experimental conditions so that a diffusing screen illuminates an object: illuminating the screen by a source, relative positions of the screen, the object and the source, etc.

Equipment: opaque object (ball), source, white screen.

Description: the source is placed next to the object (without illuminating it); the white opaque screen is placed facing the ball and the source.

Observation:

	<p>To illuminate an object not lighted directly, the white screen must be placed within the light beam so that it returns the light towards the object.</p>
---	---


3.2.2 Experiment 2 : the object must be opaque and diffusing


Objective: Compare the illumination obtained by a diffusing object according to its optical properties (colour, transparency, etc).

Equipment: source, white opaque object (ball), opaque black screen, coloured opaque screens, translucent slide.

Description: configuration of the previous experiment wherein the white screen is replaced by a black opaque screen with an opening, or a translucent slide, or a coloured screen.

Observation: The screens studied here are either non-diffusing because the black screen absorbs light, or partially diffusing because the coloured screens do not allow the ball to appear white. With these experiments, it is therefore possible to specify the concept of "secondary source". The use of a translucent screen, which diffuses only a part of the light that it receives, helps refine the definition.

	<p>With the black screen, the ball is not illuminated.</p>
	<p>With the green screen, the white ball appears green.</p>

3.2.3 Experiment 3 : The Moon and “Earth shine”


Objective : Experimentally reproduce the phenomena of diffusion responsible for earthshine.

Equipment : source, large and small white opaque balls (the Earth and the Moon)

Description : The Moon is placed between the source placed at the far end of the box (the Sun) and the Earth, almost aligned, in such a way that the Moon is seen as a crescent from the Earth, and thus from the opening of the mini lab.

This experiment may be repeated for the study of “lunation”.

Observation: Notice that the actual shadow of the Moon – ball, which appears black next to a clear crescent in the absence of the Earth, appears grey when it is illuminated by the light of the large white ball.

3.2.4 Conclusion

The Moon is an object diffusing the light of the Sun, and it sends back less light when it is quarter than when it is full.

As for "Earth light" on the Moon, it can often be observed three days after the new moon: it is referred to as "earth shine". The Earth also has the power to diffuse in all directions the light that it receives from the Sun, and the Moon "as backlight", almost in the same direction as the Sun, receives a part of it.
The planet along with the Moon in this photograph is Jupiter.

3.3 Question 3 Why does the planet Mars appear red ?

Research: How to receive red light in the eye?


3.3.1 Experiment 1 : The light is red

Objective: Receive in the eye red diffused light from a red light source.

Equipment: source, red filter, white ball.

Description: we place the red filter in front of the source and illuminate the white ball located at the far end of the box.

Observation: the ball can send back only the light that it receives and therefore appears red.

The white ball is illuminated by red light (white LED + red filter).

3.3.2 Experiment 2: the diffusing object is red

Objective: Receive in the eye red diffused light from a white light source and a red diffusing object.

Equipment: source, red filter, red object (if this is not available, a red diffusing screen)

Description: It is possible, firstly, to replace in the previous experiment the white ball by a red object. Then with the filter removed, the red object is observed in white light.

Observation: Analysis of light from the source will show that red radiation is only one of its components, it will be easy to show that the red object diffuses only a part of the light that it receives.


An eye in place of the ball will perceive the red light diffused by the coloured screen.

3.3.3 Conclusion

The planet Mars, because of the nature of its red-orangy soil containing numerous iron oxides, sends reddish light to the Earth. At the pole (on the top left of the planet), the ice cap sends white light from the Sun.

On this photograph, the stars are tracked in their apparent motion around the pole star for 10 minutes, with a lens of 28 mm focal length. In Taurus, to the right of Pleiades, the planet Mars appears red-orangy.

3.4 Question 4 : Why can we measure the distance between the Earth and the Moon with laser beam ?

Research: How does light travel ?

3.4.1 Experiment 1: marking the light path

Objective: Demonstrate the rectilinear propagation of light using sightings.

Equipment: source, 3 black slides with slit, half sheet of A4 paper.

Description: we place a half sheet of A4 paper at the far end of the black box; we then approximately align with the source 2 slits about twenty centimetres apart; we place the

eye in alignment with the system and place a third slit between the two previous ones so that the source is still seen. We then note the position of the 3 slits and the source.

Observation: the observer can see the source only when the three slits are aligned.


3.4.2 Experiment viewing a light beam

Objective: Demonstrate the rectilinear propagation of light by visualising beams.

Equipment: source, black slide with round (then square) opening, translucent slides.

Description: place the slide with the opening a few centimetres away in front of the source to obtain a beam of small cross-section, then place the translucent slides parallel to the first one with the hole.

Observation: We will observe on the translucent screens “circular light spots” whose diameter increases with the distance from the source. We will note that these translucent slides are partly transparent as the light reaches the next one, and partly diffusing as they each illuminate the previous screen.

3.4.3 Conclusion

Laser light propagates in a straight line. Thanks to this property, we can measure the Earth-Moon distance. The laser beam marks the line passing through the telescope that has emitted it; the mirror placed on the Moon reflects the laser beam and sends it back to the firing point on the Earth.

The laser is fired from the observatory at Côte d'azur, on the Calern plateau near Grasse, and it comes under the control of the Nice observatory. Additional information can be found on the websites of these organisations: "wwwrc.obs-azur.fr/cerga/laser/" or "www.obs-nice.fr"

For further study on sundials, information can be viewed on the website "http://fr.wikipedia.org/wiki/cadran_solaire"

3.5 Question 5 : Why is the shadow of the Earth on the Moon round during a Lunar eclipse ?

Research: On what do the position, shape, size and colour of shadows depend?

3.5.1 Experiment 1: how to produce a shadow?

Objective: Predict the position and shape of shadows in the case of a point source.

Equipment: source, transparent slide with a black square at the centre, white opaque screen.

Description: we observe on the screen variations in the shape of the shadow according to the object (black square) that we move or pivot along its vertical axis.

Observation: using the black square without any thickness helps concentrate on the properties of the shadow cast without the risk of confusion with the shadow on the obstacle. We can first study the variations in the size of the shadow according to the distance of the slide to the source. Using the round obstacle instead of the square shape will help define

the shadow in the form of an ellipse, which could then simplify the description of the terminator² on the Moon at the time of the study of phases.


3.5.2 Experiment 2: The color of the shadow

Objective: Notice that the shadow cast remains black even in the case of a coloured source

Equipment: source, colour filter (green or red), transparent slide with a black square at the centre, white or coloured opaque screen.

Description: To show that the shadow cast remains black whatever the colour of the light or the screen. It is not useful during these experiments to modify the distances from the obstacle to the source or to the screen or to vary the shape of the shadow.

Observation:

Note: It is possible to suggest as an exercise to place at an angle a white or yellow opaque screen close to the source, to note the appearance of a second shadow cast and ask why these shadows are no longer actually black but rather grey or dark yellow. A black shadow is a zone where there is a total absence of light.

3.5.3 Experiment 3 : Natural shadow and cast shadow

Objective: Interpret natural and cast shadows.

Equipment: source, white ball on rod, white opaque screen and colour filters.

Description: The natural shadow on the ball and the cast shadow on the screen will be all the more clearer when the source is point-shaped and the ball is away from the source.

Observation: We could observe the shadows on the obstacle and on the screen. The shadows remain black in coloured light (filter placed after the source).

² The terminator is the line separating the shadow zone and the lighted zone of the Moon or more generally a planet.


3.5.4 Experiment 4 Observing the shadow cone

Objective: Demonstrate the existence of the shadow cone.

Equipment: source, transparent slide with a black round at the centre, translucent slides.

Description: The four translucent slides are placed parallel to the one with the obstacle and they divide the shadow cone into fine black sections.

Observation:

	<p>We observe on the translucent slides round shadows whose diameter increases with the distance to the source.</p> <p>This experiment can be compared to the one suggested to mark the light beam, and thus reaffirm the rectilinear propagation of light.</p>
---	---


3.5.5 Experiment 5 : Analysis of the shadow cone

Objective: Interpret the shadow cone from a plot.

Equipment: source, white ball on magnetic support, coloured opaque screen, two black slides with slits crossing each other so that they only create a small hole, a half A4 sheet to make plots.

Description: The two cross-slits are moved in front of the screen, perpendicular to the beam, or else, after having removed the screen, an eye is placed behind these two assembled slides (behind the hole).

Observation:

Depending on whether the hole is in the shadow zone (right-hand photo) or in the illuminated part (left-hand photo), a small white luminous point (a small square "pixel") may or may not appear on the coloured opaque screen, or else, if the screen is not there, we could see the source from the small hole if this is not situated in the shadow cone. We can note on a sheet of paper placed below, the position of the small hole after which the source disappears or becomes visible again. By doing this several times at different distances from the ball, we can plot the projection of the shadow cone on the paper.


Note: the two cross-slits may be replaced by the black slide with a square central opening: the analysis of different zones will be less accurate but we would be able to study better the shadow/light limit.

3.5.6 Experiment 6: Is the Earth flat or Round


Objective: Interpret shadows.

Equipment: source, white large ball (Earth) and small ball (Moon), transparent slide with a black round, white opaque screen.

Description: The experiment can start by obtaining two identical round shadows on the flat screen, the obstacles being hidden by strips of paper.

According to you, what are the obstacles that create this same shadow on the screen?
Remove the masks.

On the left, the shadow is created by the transparent slide with a black round at the centre, whereas on the right it is the large opaque ball that gives the same result on the screen.


If we now pivot slightly each obstacle along a vertical axis, the shadow of the ball does not change whereas that of the slide becomes an ellipse.

We then replace the white screen by the small Moon-ball.

Observation: The Earth-ball must be placed quite away from the source so that its shadow cast is of a reasonable size. Only a round Earth can cast a shadow with a circular outline, whatever the position of the source or the obstacle. In ancient times, this was one of the arguments put forward to reject the model of a “flat Earth”.

The phenomenon of a Lunar eclipse is described in greater detail in the last paragraph.

	<p>Shadow of the Earth on the Moon.</p>
	<p>Shadow of the black round, positioned perpendicular to the axis of the beam, on the small white ball on the rod (Moon).</p>

Shadow of the black round, positioned at an angle with respect to this axis on the small white ball on the rod (Moon).

3.5.7 Conclusion

Lunar eclipse.

We will note that on Earth, the half night side corresponds to a natural shadow and on the Moon emerges as a shadow of large diameter: shadow of the Earth (cast shadow).

3.6 Question 6 : Why does the Moon change its appearance every day

Research: Do the relative positions of the Sun, the Earth and the Moon help explain the phases?

3.6.1 Experiment 1: description of lunation with two balls exposed to light

Objectives: Simulate the phases of the Moon over a complete lunation.

Equipment: source, large and small balls

Description: source at the far end of the box, large ball at the centre and small ball held by hand.

Observation: By rotating the small ball (Moon) around the larger one (Earth), in the light of the source (Sun), counter clockwise, we can observe the appearance of the small sphere in its different positions by viewing it just above the fixed ball. We could even imagine the observer placed on the Earth (marked by a lightly fixed sticker) by trying to assess the time corresponding to this observation point (knowing that it is noon at the centre of the lighted face of this large Earth ball).


We could use the sheet on page 40 with two concentric circles in this more thorough research:

- The inner circle graduated in hours (to place the observation time),
- The outer circle, trajectory of the Moon, graduated in days.



The pattern on page 40 may be cut out or reproduced, and placed at the far end of the box (noon towards the Sun). The diameter of the trajectory of the Moon must measure about 12 centimetres on the reproduction.

Note that the sidereal period of the Moon (complete revolution on its trajectory) is a little less than 28 days while the duration of a lunation is 29.5 days (because of the movement of the Earth around the Sun during this period of almost one month, which delays the new alignment between the Sun, the Moon and the Earth).

The black background of the drawing will help eliminate stray light for observing the phases as shown in some of the examples given below:


We will note that the terminators (lines separating day and night) on the Earth and on the Moon will have to be identical with the Sun very far away (at 150 million kilometres). In fact, on these photos they are not exactly identical because the source is too close to the two balls. We could, for greater clarity, propose an observation outside the box.

	<p>Gibbous moon on the eleventh day, reaching culmination at the beginning of the night.</p>
	<p>Last quarter, 22 days after the new moon, reaching culmination at Sunrise.</p>

3.6.2 Conclusion

The phases of the Moon depend only on the position of the Moon with respect to the Earth and the Sun.

3.7 Question 7 : Why in France are Solar eclipse much rarer than Lunar eclipse ?

Research: Can we mark on a model the visibility zones of each of the eclipses?

3.7.1 Experiment 1 : Simulation of a Solar eclipse

Objective: Establish in which configuration of the Sun-Earth-Moon system can a Solar eclipse occur.


Equipment: source, small and large balls.

Description: source (Sun) at the far end of the box, ball on support (Earth) close to the opening, ball on rod (Moon) held in front of the Earth but quite close to it so that the shadow cast is quite small.

Observation: We can move the small ball slightly with respect to the Sun – Earth line to show that every year, one or two eclipses may occur, but only small regions on the Earth, different every time, will be able to view them.

We could rotate the Earth/ball to introduce the “visibility band” concept.

We could finally position the Earth on the diagram graduated in hours and in lunation to estimate the visibility time of the eclipse.

3.7.2 Experiment : Simulation of a Lunar eclipse

Objective: Establish in which configuration of the Sun-Earth-Moon system can a Lunar eclipse occur.

Equipment: source, Earth-ball and Moon-ball held in the shadow cone.

Description: source (the Sun), white ball on magnetised support placed quite far away from the source (the Earth), small white ball on rod maintained behind the first and quite close (the Moon).

Observation: We will interpret (with estimations on the diameters, possibly) the shadow zones. We could show that the eclipse is visible at night for inhabitants of more than half the Earth (obscured by its natural shadow). We could also specify that the Sun, even though very far away, has a very large diameter and that the shadow cone shrinks moving away from the Earth.

3.7.3 Conclusion

Solar eclipses are visible around noon, and Lunar eclipses occur in the middle of the night.

4. Photo credits

Astronomy photographs appearing in sections 3.2.4 page 29, 3.3.3 page 31, 3.5.7 page 38 were created by Jean-Luc SINGER, astronomy animator and holder of a BEATEP in astronomy and computer science – <http://www.astropixel.fr>.

The other photographs and diagrams were created by Françoise CLIPET and Jean-Luc FOUQUET, professors of physical sciences at the Saint-Martin de Ré college. M. FOUQUET is also an active member of the Comité de Liaison Enseignants et Astronomes (CLEA) (Teachers and Astronomers Liaison Committee).

5. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux