# **Optique FR**

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# ***Optics GB***

*Geometrical Optics GB*

Réf :  
204 084

Français – p 1

English – p 5

Version : 8105

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
***Reflection-refraction set GB***

## 1. Description

Image 1

### 1.1 Composition

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
### 1.1 Composition

*Product composition*

- 10. Triangular metal support, screen printed
- 11. Rotating arm
- 12. Graduated rotating disc
- 13. Transparent tank
- 14. Semi-cylindrical lens made of polished Plexiglas
- 15. Cover for transparent tank
- 16. Fixing screw for tank and lens
- 17. Plane mirror on support
- 18. Rod for vertically placing the triangular metal support (see image 1)

The device is intended for use in horizontal or vertical position.  
 Set delivered in a storage box.

### **Recommended Additional Equipment (not provided)**

	- Multi-beam laser	<b>Ref. 201 039</b>
or	- Light source 12 V / 30 W	<b>Ref. 211 066</b>
or	- Laser diode	<b>Ref. 201 032</b>
	+ Kit laser plane	<b>Ref. 202 076</b>

## **2. Implementation**

Place the device horizontally on a table, the rubber buffers ensure its secure fit during handlings.

Place the plane mirror, the semi-cylindrical lens or the transparent tank on the graduated disk depending on the experiment to be carried out (the plastic film present on the Plexiglas lens should be removed before use).

Place the graduated disk in null position: graduation 0 of the disk facing the red arrow on the triangular support.

Position the light source or the laser on the most suitable graduation (see below) and then align the light beam with the device central axis.


Use the rotating arm to find the reflected or refracted beam when the latter is not clearly visible on the graduated disk. The arm makes it easier to read the angles of the reflected or refracted beams for several disk positions (angle of reflection of 15° minimum).

Using the rod supplied, place the whole on an optical stand base accepting 10 mm diameter rods.

The experiments listed below are achievable in vertical position for classroom demonstrations. Think of using the magnetic pad (ref. 291196) to set your light sources vertically.


### **Configuration 1 : use the multi-beam laser**

Place the laser in such a way that the beam exit is on index n°6.


**Configuration 2 : use with the light source 12 V / 30 W**

Place the source as close as possible to the graduated disc, in such a way that the beam exit is on index n°1.

**Configuration 3 : Use with the laser diode and the plane laser kit**

Place the laser diode at the end of the triangular support and tilt the laser diode forwards, and then place the laser plane kit on index n°5.

### 3. Technical Specifications

<b>Graduated disc :</b>	Ø 230 mm
<b>Plane mirror surface :</b>	about 50 x 10 mm
<b>Plexiglas lens :</b>	Ø 200 x 10 mm
<b>Transparent tank :</b>	Ø 200 x 20 mm
<b>Total dimensions :</b>	550 x 230 x 40 mm
<b>Rod :</b>	Ø 10mm
<b>Weight :</b>	about 1 kg
<b>Warranty :</b>	2 years

## 4. Handlings Examples

### 4.1 Plane Mirror: Verification of Descartes Law

Place the mirror at the center of the graduated disc in such a way that the reflective side corresponds to the diameter  $- 90^\circ / + 90^\circ$  (configuration 2).

The light beam being properly centered, change the angle of incidence on the mirror by turning the graduated disc starting from its original position (graduation 0 facing the red arrow).

For each position, measure the angle of incidence  $i$  and the reflected angle.

### 4.2 Refraction: Determination of the Index of a Medium

#### - With the Plexiglas semi-cylindrical lens ( $n = 1.5$ )

The beam being centered, place the lens at the center of the graduated disc in such a way that the beam arrives normally to the plane side, and attach the lens to the disc using 2 fixing screws. The spots of the reflected and refracted rays are then in the extension of one another ( $i = 0$ ).

Turn the tray to change the angle of incidence  $i$  from  $2^\circ$  to  $2^\circ$  until reaching  $15^\circ$ , and then from  $5^\circ$  to  $5^\circ$  until  $80^\circ$ . Note the angle of incidence  $i$  and the refracted angle. Deduce  $\sin i$  and  $\sin r$ .

Verify the relationship  $n_1 \sin i_1 = n_2 \sin i_2$ , with  $n_1$  (air) = 1 and  $n_2 = 1.5$ .

Verify the reversibility principle of light rays.

Verify the total reflection phenomenon and determine the limit experimental angle.

#### - With the transparent tank

Fill the tank with water until it is almost half-way to the top. Place the tank with the same care as the semi-cylindrical lens. Position the semi-cylindrical cover on the front of the tank to remove the refracted beam in the plastic part of the tank (see figure 2).

**Proceed as described above with the semi-cylindrical lens to note the angles. Verify the relationship  $n_1 \sin i_1 = n_2 \sin i_2$  for the water with  $n_1$  (air) = 1 and  $n_2 = 1.33$  (water).**

## 5. After-Sale Services

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

### Free service\*

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux