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Geometrical optics

Optical instrument

Réf :
204 350

Français – p 1

English – p 7

Version : 4208

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09 69 32 02 10 *

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1. Introduction

This model represents a simplified camera designed for educational purposes to illustrate the basic functioning of a camera. It is suitable for an optical bench. The goal is to study the formation of an image in a camera.

Educational Applications:

- **Illustration of Optical Principles:** This setup is ideal for demonstrating the principles of light reflection in a reflex camera.
- **Simulation of Reflex Camera Operation:** Real-time visualization of the light path inside a camera, from entering through the lens to exiting through the viewfinder.
- **Understanding Image Formation:** The setup also shows how the image is formed through the lens and perceived through the viewfinder by the user.

2. Composition

This product consists of the following elements:

1. Pentaprism: A component used to correctly reflect the image towards the viewfinder. It plays a crucial role in reflex cameras to orient the image perceived by the user.


2. Eyepiece Lens: Located behind the pentaprism, the user can observe the formed image through it.

3. 45° Removable Mirror: Located inside the housing, inclined at 45°, designed to reflect light from the lens. This mirror gives the system its 'reflex' name as it reflects light upwards and then lifts during the photo capture to allow light to reach the sensor.

4. Ground Glass Screen: Simulates the sensor during photo capture. When the mirror is raised, the image is formed on this screen.

5. Rod: Used to mount the model on an optical bench.


6. Transparent Housing: A rectangular housing made of transparent plastic allows observation of the internal components of the model.

3. Reflex Photography:

When a scene is observed with a reflex camera, the light from the scene passes through the lens. Inside the camera, a 45° inclined mirror is positioned directly behind the lens. This mirror reflects the light towards the pentaprism, which corrects the image (initially inverted by the mirror). The pentaprism redirects this corrected image towards the viewfinder. The photographer looks through the viewfinder and sees the image as it will be captured. This direct visualization is possible thanks to the mirror and pentaprism system, which is at the heart of reflex camera operation.

When the photographer is ready to take the picture and presses the shutter button, a complex mechanism is triggered inside the camera. The 45° mirror quickly flips up out of the light path, allowing light to reach directly the sensor or film located behind the mirror. With the mirror raised, the shutter opens to expose the sensor or film to the light, capturing the image. The exposure time depends on the settings chosen by the photographer (shutter speed). After exposure, the shutter closes, and the mirror quickly returns to its original position, restoring the view in the viewfinder.


Exploded view of a reflex camera

4. Model Operation Principles:

In our setup, the scene is represented by the lantern and the F object. The beam passes through the lens, which acts as the objective, then reaches the 45° mirror. The light reflected by the mirror then reaches the pentaprism, which redirects the image towards the viewfinder. Finally, the observer views the F object through the eyepiece lens. To simulate photo capture, simply tilt the mirror to the horizontal position so that the image is formed on the ground glass screen.


4.1 Scheme

4.2 Set-up

The set-up consists of the following elements:

- **Lantern + object F:** The device located at the far left of the setup, this source produces a beam of light directed towards the camera, simulating the role of an object.
- **Lens (Objective):** Positioned in the center of the setup, this lens acts as the objective in a reflex camera.
- **Camera Model:** Located at the far right, this model is the one previously described. It illustrates the principle of how a reflex camera works.

5. Technical characteristics and complements

5.1 Technical characteristics

Eyepiece Lens (reference 202219) ; $f=+150$ mm
Pentaprisim : optical surface 25*28mm
Mirror: 50*50mm
Screen
Rod

5.2 Complements

Lantern
Object F
Lentille : $f=+300$ mm
Optical bench
Optical mounting

6. After-sales service

The warranty is for 2 years.

For any adjustments, please contact **Technical Support** at **09 69 32 02 10***.

The equipment must be returned to our workshops, and for any repairs or spare parts, please contact:

JEULIN – After-Sales Service
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX France

09 69 32 02 10 *

*Local call rate, not surcharged.