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Didactic Physics

Optical Measuring Instruments

Réf :
207 502

Français – p 1

English – p 6

Version : 0201

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1. Presentation

1.1 Introduction

This model, supplemented by a parallel beams source (of type mirrors source ref.201038 or 211066) allows to conclusively replicate Eratosthenes' experiment on a student's table.

It will allow the gradual implementation of the link between physics and mathematics and provide the students with the essentials tools to build a problem solving strategy through the experiment.

With the measurements, the student will then be able to calculate and deduce a comparable length of the meridian and deduce the theoretical value of the Earth's radius.

1.2 A bit of History

Eratosthenes (275-194 BC) had read that on the summer solstice day in Syene (Egypt), sunlight reaches the bottom of a vertical well at noon, while on the same day and time, in Alexandria, a vertical obelisk always has a shadow. Eratosthenes also knew that in order to go from Syene to Alexandria, a camel caravan needed 50 days to cover a distance of 100 stadions per day (the stadion was the unit of length at that time and equaled around 160m).


2. Description

The model consists of a tray on which is represented the Earth in cross-section by longitude 30° .

Syene and Alexandria are indicated as well as the distance separating them (5000 stadions, which is around 800km according to the information known on Eratosthenes).

An arm graduated in mm rotates about the center of the Earth and projects perpendicularly to the Earth's surface. Illuminated by a parallel light beam, it projects a shadow that we can directly measure on a small graduated tray, fixed perpendicularly to the arm.

This model is especially designed to work with the mirrors light source that can produce a correctly parallel light beam. The model and the source are placed on the table. The source is mounted on the tray edge and reproduces the sunlight during the summer solstice.

To carry out the measurements easily, you may slide the small tray along the arm while remaining strictly perpendicular to the arm. The shadow extends, and the measurements to calculate the arm rotation angle are even more precise.

3. Experimental Protocol

Use the 2 parallel lines represented on the model to set the parallelism of the light source beam.



Set the arm on Syene, place the light source facing Syene and measure the length of the shadow on the small tray: it must be zero.

Set the arm to Alexandria (for example), move the light source along the model and measure the shadow on the small tray.

These 2 measures allow calculating the angle α between the sun's rays and the vertical to Alexandria at the summer solstice, and the Earth's radius, knowing the distance Syene - Alexandria.

The manipulation is also possible between 2 different latitudes provided that the angle remains slight (less than 8°).

With the theorem of internal/external angles, 2 angles α seem equal.


4. Results and Digital Applications

4.1 Experimental Values


We are looking for the value of angle α , corresponding to angle \widehat{CDE} .

Measure of the shadow in Syene:

The observation shows the absence of shadow on the moving ruler.
 Angle $\widehat{CDE} = 0^\circ$ -- The points C and E coincide.

Measure of the shadow in Alexandria:

We note a shadow of length = 6 mm
 The angle between the rotating arm and the ruler are constant and still: 90°
 Since angle $\widehat{DCE} = 90^\circ$, it is then a right triangle and the trigonometric relationships apply:

$CE = 0.6 \text{ cm}$ (side opposite to angle \widehat{CDE})

$CD = 4.7 \text{ cm}$ (side adjacent to angle \widehat{CDE})

Based on the relationship $\tan \widehat{CDE} = \frac{[\text{opposite}]}{[\text{adjacent}]}$ we then have $\tan \widehat{CDE} =$

$$\frac{[CE]}{[CD]} = \frac{0.6}{4.7}$$


$$\tan^{-1} \frac{0.6}{4.7} = 7.27^\circ$$

We have $\widehat{CDE} = \widehat{ABC} = 7.27^\circ$

4.2 Results & Calculations

This angle corresponds to $\frac{360^\circ}{7.27^\circ} = 49.5^{\text{th}}$ of the length of a meridian M.

We know that this angle also corresponds to 5000 stadions, being 800km:

Therefore, we deduce that travelling 800 km on the Earth's surface means travelling an angle of 7.27° .

7.27° corresponds to $1/49.5^{\text{th}}$ of a full rotation of the meridian.

Therefore, 800 km also corresponds to 49.5^{th} of a full meridian length.

We then deduce that the length of the meridian M:

$$[M] = 800 \times 49.5 = 39600$$

Length of a meridian = 39 600 km

The relationship linking a meridian to the radius of the Earth is such that:

$$M = 2\pi r$$

We deduce:

$$r = \frac{M}{2\pi}$$

Hence:

$$r = \frac{M}{2\pi} = \frac{39\,600}{2\pi} = 6302 \text{ km}$$

The actual theoretical data corresponds to:

Length of the meridian = 40 075 km

Radius of the Earth = 6371 km

4.3 Error Calculation

$$Error = \frac{Theoretical\ value - Experimental\ value}{Theoretical\ value} \times 100$$

$$\text{Error on the measurement of the meridian: } Error = \frac{40075 - 39600}{40075} \times 100 = 1.2 \%$$

$$\text{Error on the measurement of the radius: } Error = \frac{6371 - 6302}{6371} \times 100 = 1.1 \%$$

Eratosthenes' method made it possible to calculate the length of a meridian and the Earth's radius with only 1.2%!

5. Technical Specifications

Tray dimension: 300 x 420 x 35mm
Rotating arm dimension: 260 x 15 x 10 mm
Moving ruler dimension:
Mass: 1 kg

6. Recommended Accessories

Light source with mirrors – halogen	ref. 211066
Light source with mirrors – LED	ref. 211038

7. After-Sale Service

The warranty is of 2 years.

For all settings, contact out **Technical Support** on **0 825 563 563**.

For all repairs, spare parts or return to our workshops, kindly contact:

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux