Ultrasons

Ultrasounds

Réf :
221 001

Français – p 1

English – p 3

Version : 8111

Paire de diapasons
Pair of tuning forks

1. Objectifs

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0 825 563 563*

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1. Purpose

The tuning fork pair is used to demonstrate the properties and characteristics of sound waves.

2. Composition and assembly

- 2 sounding boards made of varnished wood,
- 2 steel tuning forks (A³), 440 Hz,
- 2 sliding pendulum weights,
- 1 hammer.

Place each tuning fork on its respective sounding board.

One of the tuning forks may receive the pendulum weights allowing the tuning fork frequency to be varied.

3. Experimental examples

3.1 Sound pitch

This experiment uses the tuning fork fitted with the pendulum weights.

- Install and secure the two pendulum weights at the upper end of the tuning fork.
- Strike and listen.
- Then place the weights at the lower end of the tuning fork.
- Strike, listen, and comment.

In order to visualize the "physical difference" between the two sounds, you can use a microphone connected to a storage oscilloscope (or datalogging).

Then perform frequency measurements in the two cases described above.

3.2 Beats

Place the two boards with the tuning forks side by side.

For different weight positions, strike each tuning fork (using the same force) in quick succession.

Sound intensity maxima and minima can be heard. This is the so-called "beat" phenomenon, which is all the stronger as the difference in tuning fork frequency is small (weights approx. at mid-height on tuning fork).

3.3 Resonance

- Place the two boards with the tuning forks facing each other.
- Strike a tuning fork then quickly hold it with your fingers.

If the tuning fork frequencies are close, or better still identical, the second tuning fork will start to vibrate.

Repeat the experiment with different weight positions and also without any weight.

4. After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux