# **Mécanique**

Etude du pendule

# ***Mechanics***

*Study of the pendulum*

**Ref :  
222 008**

Français – p 1

English – p 11

Version : 7005

**Kit pendule simple et pesant**  
***Simple and weight pendulum kit***

## SOMMAIRE

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**NOTES**

## CONTENTS

1	Description.....	12
2	Composition and technical characteristics.....	13
3	Accessories (not supplied).....	13
4	Experiment: Study of the simple pendulum by video recording	13
4.1	Purposes of the experiment	13
4.2	Equipment used	14
4.3	Set-up of the assembly	14
4.4	Video recording	14
4.5	Processing the video	14
4.6	Analysing the measurements obtained	15
5	Experiment: damping, aperiodic and critical regimens .....	16
5.1	Purposes of the experiment	16
5.2	Equipment used	17
5.3	Set up and adjustment	17
5.4	Video recording	17
5.5	Processing the video	17
5.6	Analysis of measurements obtained	17
6	Experiment: Energy studies.....	18
7	After-Sales Service .....	19

## 1 Description

This assembly was designed specifically for basic study of the mechanical oscillation of a simple or weight pendulum.

**It allows:**

- Studying the specific period of a simple pendulum
- Studying the damping of a weight pendulum
- Verification of the law of isochronism of the minor oscillations of a simple pendulum
- Learning about the concepts of amplitude, period, pseudo-period, aperiodic and critical regimens

The composition of the simple and weight pendulum kit allows varying one single parameter at a time: weight, length of the pendulum or damping.


The assembly can be modified in just a few seconds thanks to an ingenious system for attaching the various components:

- System for hooking-on the pendulum using a bracket clipped onto a stepped pulley
- Weight with screw which can be attached to the simple weight pendulum to adjust its height
- A damping disk which can be clipped onto the weights and orientated to vary the damping
- A removable blade to study aperiodic and critical regimens with the weight pendulum

The colour marker on the simple pendulum is intended to facilitate computer processing during video operation.

## 2 Composition and technical characteristics

- ① A simple pendulum support: an attaching bracket on a stepped pulley, wire length 1 metre, 1 accessory for attaching the adjustable height weights, 1 colour marker.
- ② Weight pendulum support: an attaching bracket on a stepped pulley, graduated rigid rod  $L = 600$  mm
- ③ 2 weights, 100 and 150 g with attaching system for clip and graduated rod.
- ④ Damping disk (50 g) with system for hooking onto the 100 g weight.
- ⑤ Blade to study aperiodic and critical regimens, attached to the end of the graduated rod.
- ⑥ A stepped pulley with high quality ball bearing.

## 3 Accessories (not supplied)

- A support base in cast iron, vice grip, tank.
- A simple chronometer for an initial approach to quantifying pendulum laws.

For more detailed study, a webcam and video processing software such as the science workshop, for computer input and processing of the data.

It is also possible to computerise the data using your Ex.A.O. console with an electronic stepped pulley and an angle module which supplies voltage proportional to the angle of the pendulum. This voltage can be used directly on any Ex.A.O type interface.

## 4 Experiment: Study of the simple pendulum by video recording

### 4.1 Purposes of the experiment

To study the specific period of the simple pendulum: influence of the weight and length of the pendulum.

To demonstrate the influence of the angle on amplitude.

Verification of the law of isochronism of minor oscillations of the simple pendulum.

## 4.2 Equipment used


- The simple pendulum in the “simple and weight pendulum kit”,
- A webcam,
- Data input and processing software.

## 4.3 Set-up of the assembly

- Attach the stepped pulley on the cast iron support base using the vice grip.
- Clip the pendulum on to the pulley at the 2 notches taking care to add the colour marker.
- Connect a webcam to your PC and position it opposite, perpendicular to the pendulum.

## 4.4 Video recording

In the video workshop, accessible using the icon , go to the “fast recording and video tab”, parameter 10 or 15 seconds duration for the sequence.

Start the pendulum then start recording by clicking on the webcam icon.

## 4.5 Processing the video

Click on the automatic processing tab.

### Stage 1: Calibration


- Fix the origin of the marker
- Identify a scale

### Stage 2: Work frame

- Select the green frame in one corner using the mouse and re-size it so it covers solely the zone corresponding to the use zone.

### Stage 3: Parametering

- Select the object to be followed.


- Start processing by clicking on the green light.

## 4.6 Analysing the measurements obtained

Analysing the data acquired is performed directly with the science workshop and the modelling workshop:

### ❖ Influence of weight on the period:

#### Experiment Conditions


Different weights

Identical angles

Identical lengths

The curves are exactly superimposed. Hence it can be shown weight has no effect on the period.

### ❖ Influence of length of the pendulum on the period:

#### Experiment Conditions

Identical weight

Identical angles

Different lengths

It can be observed the period varies with the length of the pendulum.

It is possible to obtain the value of the experimental period from the graph and after modelling compare it with the theoretic intrinsic period of the pendulum. This gives:


#### Solid line curve

$$T_{\text{exp}} (\text{modelling}) = 1.14 \text{ s} \quad L = 32\text{cm} \quad T_{\text{theo}} = 2\pi\sqrt{\frac{L}{g}} = 1.13$$

#### Dotted line curve

$$T_{\text{exp}} (\text{modelling}) = 1.61 \text{ s} \quad L = 64\text{cm} \quad T_{\text{theo}} = 2\pi\sqrt{\frac{L}{g}} = 1.60 \text{ s}$$

#### ❖ Verification of the law of isochronism of the minor oscillations of the simple pendulum:


#### Experiment Conditions

Identical weights

Identical lengths

Different angles

This study indicates the influence of the angle on the period.

A close-up of the 10<sup>th</sup> period of oscillation shows that for small angles, the period is identical (principle of isochronism) and then for larger angles, the angle has an influence on the period.

## 5 Experiment: damping, aperiodic and critical regimens

### 5.1 Purposes of the experiment

Study of the shock absorption of the weight pendulum.  
Concept of pseudo period, a periodic and critical regimens.

## 5.2 Equipment used

- The weight pendulum, the damping disc and the blade from the "Simple and weight pendulum kit",
- A webcam,
- Data processing and input software.

## 5.3 Set up and adjustment

- Attach the stepped pulley to the cast iron support base using the vice grip.
- Clip the pendulum on to the pulley at the 2 notches taking care to position the colour marker.
- Connect a webcam to your PC and position it opposite, perpendicular to the pendulum.

## 5.4 Video recording


In the video workshop, accessible using the icon , go to the "fast recording and video tab", parameter 10 or 15 seconds duration for the sequence.

## 5.5 Processing the video

See section 4.5

## 5.6 Analysis of measurements obtained

❖ Concept of pseudo period and damping:

Modelling the curves allows comparing the pseudo periods and indicating they are identical (here  $T = 1.42$  s).


This also gives the value of the damping coefficient.

### Experiment Conditions


Identical weights

Identical lengths

Variable damping (the damping disc is inclined more or less to vary the damping)


❖ Concept of aperiodic and critical regimens:


The blade is attached to the end of the weight pendulum.  
The tank is filled with a more or less viscous liquid (example: oil, water, glycerine)  
The movement is recorded as previously.

Example of recording for damping in water.  
Recording of critical regiment

## 6 Experiment: Energy studies

To go further in the analysis of the data, the Generis 5 software allows performing all the calculations necessary to study the energy aspects.

Example of use in Ex.A.O with an ESAO 4 interface, an electronic stepped pulley P/N: 453109 (optional accessory not supplied) and the angle module P/N: 453110 (optional accessory not supplied).

## **7 After-Sales Service**

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27 019 EVREUX CEDEX FRANCE**  
**+33 (0)2 32 29 40 50**

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

**Free service \* :**  
**+ 33 (0)2 32 29 40 50**

*\* Call cost not included*