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Periodic phenomena

Réf :
222 033


Français – p 1

English – p 4

Version : 8111

Pendule à paramètres variables ***Pendulum with variable parameters***

1. Description

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
1. Description

- | | |
|--------------------|--------------------------------|
| 1. Winder | 6. Recess |
| 2. Suspension cone | 7. Pendulum retaining screw |
| 3. Threads | 8. Face of the pendulum |
| 4. Mass support | 9. Lugs for storing the masses |
| 5. Suspended mass | |

This device respects the conditions of a simple pendulum model: non-extensible thread of sufficient length and negligible mass, control of the oscillation amplitude.

The pendulum is set up on a standard laboratory support equipped with a 6 to 16 mm diameter rod or a 15x15 mm square-section rod.

The device helps:

- 1) Control the tilt angle of the pendulum from 0 to 40°
- 2) Modify and measure the length of the thread easily: On every rotation of the winder (1), the thread's (3) length varies by 10 cm
- 3) Have 3 objects (5) of the same volume and known mass set up on the support (4) and arranged on the 3 lugs (9):
 - PVC, 8 g
 - Aluminium, 14 g
 - Brass, 40 g

We can easily study the influence of these parameters on the value of the period and show that the only parameter that influences it at a given location is the length of the thread.

By varying the amplitude of the oscillations, we can demonstrate the limits of validity of a simple pendulum model.

2. Initial adjustment of the length of the thread

The device is set up and adjusted at the factory. However, in case the thread accidentally breaks, it can be readjusted to the correct length in the following manner:

1. Wedge the thread winder (1) on its notch
2. Position the mass support (4) (without the mass) in the recess (6)
3. Pass the thread through the mass support (4), stretch it and fix it to the support with a screw.

The length of the thread measured from the suspension cone (2), up to the centre of inertia of the object is then exactly 10 cm.

3. After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

+33 (0)2 32 29 40 50

Experiment sheet : Simple pendulum

1. Objectives

To build a simple pendulum that oscillates every second.

2. To reflect

- What is a simple pendulum? How to make one easily?
- What does the expression “mark every second” mean?
- What is an oscillation of a pendulum called? The period of an oscillation? The amplitude of an oscillation?
- How to measure the amplitude of an oscillation?
- What are the parameters that can be varied to modify the period of a simple pendulum?
- When does the period of a simple pendulum remain constant for a sufficient duration?
- Why is the study of a simple pendulum of historic importance?

3. Carry out

1) Carry out the experiments that will help determine the parameters (physical quantities) on which the period of a pendulum depends:

- ❖ Select a parameter and determine its influence on the period
 - Description of the experiment
 - Results
 - Conclusion
- ❖ Select another parameter and determine its influence on the period
 - Description of the experiment

- Results

- Conclusion

- ❖ Select a third parameter and determine its influence on the period
 - Description of the experiment

 - Results

 - Conclusion

2) Write down the overall conclusion here

4 Process

- 1) Draw a curve linking the period T of a pendulum and its length L . For example, draw $T=f(L)$ or $T = f(\sqrt{L})$ or $T^2=f(L)$.
- 2) What should be the length of the thread so that the pendulum marks every second?
- 3) Make a pendulum of this length and check with a stopwatch that it oscillates every second.
- 4) What must be done so that this pendulum oscillates without stopping?

5 Conclude

Imagine that we are in the 17th century and you need to sell pendulum clocks, the latest novelty. What are the arguments you will use to convince the buyers?

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux