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## *Periodic phenomena*

Réf :  
222 037

Français – p 1


English – p 4

Version : 8111

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### *Melde's membrane vibrator*

## 1. Description

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
**0 825 563 563\***

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## 1. Description

The vibrator consists of a loudspeaker membrane connected to a rod that vibrates vertically.

The end of the rod is fitted with a tightening screw to fix a cord or a spring.

Power is supplied by a low frequency generator providing 20 V DC (not supplied) which helps adjust the frequency of vibration of the vertical rod.

**The following accessories are required for the experiments** (not provided with the device):

- LFG 2 MHz, part no. 293 047
- Modumontage support leg, part no. 701 293
- 1 square-section rod 15 mm, length. 70 cm, part no. 703 459
- 2 tightening nuts, part no. 703 099
- A thin cord of visible colour (25 metres), part no. 323 011
- A spring for standing waves  $\varnothing$  10 mm, length. 800 mm, part no. 222 002
- A graduated metre rule, part no. 313 263
- A lab-jack, part no. 701 238
- A guide pulley on the rod, part no. 323 017
- A box of masses with hooks, part no. 703 016

## 2. Technical characteristics

Melde's membrane vibrator complies with the EN61010-1 standard, installation category II, pollution degree 2.

- **Maximum voltage:** 12 V
- **Connection:** safety sockets ø 4 mm
- **Protection:** auto-reset fuse
- **Dimensions:** 135 x 92 x 80 mm
- **Mass:** 640 g
- **Guarantee:** 2 years

## 3. Principle

This is the basic device for the study of transverse and longitudinal standing waves in vibrating motion.

Melde's experiments consists in observing transverse standing waves along a string, the vibrator generates vibrations sustained electrically.

Longitudinal standing waves can easily be obtained using a spring.

All these experiments help verify the phenomenon of high-pitched resonance. Melde's experiment also helps verify the formula of vibrating

strings:

$$f = \frac{k}{2l} \sqrt{\frac{F}{\mu}}$$

## 4. Experiments


### 4.1 Transverse vibration of a string

Place the Modumontage® support at approximately 1 metre from Melde's vibrator. Fix the pulley on the rod at the same level as the vibrator using the tightening nut.

Attach one end of the string to the rod of the vibrator (using the screw), and the other to the mass with a hook. Pass the string on the pulley.

The maximum amplitude is obtained between 40 and 50 Hz for a given value of the tension of the string: for example, 2 antinodes for an additional mass of about 20 g.


## 4.2 Longitudinal vibration of a spring

Place the vibrator on the table in line with the square-section rod.

Attach one end of the spring to the vibrator and the other to the horizontal rod, as shown in the diagram below.

Adjust the frequency of vibration to obtain antinodes with maximum amplitude.

## 5. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.  
For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
468 rue Jacques Monod  
CS 21900  
27019 EVREUX CEDEX FRANCE

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

### Free service\*

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux