**Méthode EXAO®**

***EXAO® method***

**Réf :  
222 040**

Français – p 1

English – p 5

Version : 6006

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**0 825 563 563\***

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
## 1. Description

This assembly helps study free oscillations: influence of the mass, the spring constant and the damping in air.

It was specially designed for a video study.

The acquisition and video processing software Cineris (or Generis 5 Plus) helps film and easily process all the experiments that can be carried with this system using its colour reference mark.

Thus, you have access to the period or pseudoperiod measurements and to the amplitude of the movement. All this data can then be processed and modelled.

## 2. Technical characteristics

The composition of the vertical elastic pendulum helps vary a single parameter at a time:  $m$ ,  $k$  or the damping, for a simplified study of the oscillations.

The spring is fixed with lock-grip pliers.

The shaft of the mass + shaft system is pierced at its end to attach the device to the spring.

The damper disc is added (between the washers) or removed using its slot.

The colour reference mark is designed to facilitate the data processing of the video.

The threaded rod, screw and washer system ensure the retaining of the elements, their adjustment and prevent their loss.

### 3. Experimentation : Oscillations of a mass-spring system

#### 3.1 Objectives of the experiment

To study the characteristics of the static opposing force.

To identify and measure the amplitude and period of oscillation.

To display the curves: displacement  $x(t)$ , velocity  $v(t)$  and acceleration  $a(t)$  of the motion.

Comparison of displacement curves in different cases:

- Same stiffness and different masses
- Same mass and different stiffness of springs
- Same mass and stiffness with or without damping.

Check the consistency of some experimentally obtained values, for example  $(2\pi/T)$  in the equation  $a(t) = -(2\pi/T)^2 \cdot x(t)$ .

#### 3.2 Equipment used

Vertical elastic pendulum  
Cinéris software  
ToUCam Philips webcam

To establish proportionality between the restoring force of a spring and its displacement using mass  $m_1$  of 150 g and the two available springs.

Once the influence of the spring is determined, study:

- The influence of the mass: use a spring of stiffness  $k_2$  and carry out the experiment successively with  $m_1$  then  $m_2$
- The influence of damping: use a spring of stiffness  $k_2$  and carry out the experiment successively with mass  $m_2$  then with mass  $m_1$  + the disc.

#### 3.3 Setting up and adjustments

Set-up the experiment with elements that interest you.

Connect the webcam to the computer

Open Cinéris

Video record the motion (refer to the Cinéris instructions)

#### 3.4 Processing of measurements obtained


The processing of the video is done directly by Cinéris.

The motion curves can be obtained and modelled.

##### ❖ Free oscillations without damping

Example of acquisition obtained with  $k_2$  and  $m_1$

The measurement of the period is carried out from oscillation to oscillation (to check its consistency), then by measuring the duration of the maximum number of visible oscillations divided by this number.


If slight damping occurs, select a lower acquisition duration on two or three recorded periods so that the phenomenon is barely visible.

From this recording move on to study the equation of the curve, its modelling, the lines of the curves relating to the velocity and acceleration of the motion.

Compare the curves obtained with  $m = \text{constant}$  and  $k$  variable then with  $k = \text{constant}$  and  $m$  variable.

### ❖ Free oscillations with damping

Example of acquisition obtained with  $k_2 m_1$  and the disc

The measurement of the pseudo-period is carried out as in the previous case.

## 4. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.  
For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
468 rue Jacques Monod  
CS 21900  
27019 EVREUX CEDEX FRANCE

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

### Free service\*

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux