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Periodic phenomena

**Ref :
222 041**

Français – p 1

English – p 5

Version : 9001

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
1 Description

This system completes the Vertical Elastic Pendulum system for the study of mechanical resonance.

It is used for:

- Comparing the period of the exciter and the resonator.
- Effect of damping on the amplitude of oscillations.

Composition

- 1 – Eccentric motor (powered by security sockets \varnothing 4 mm, 12 V max)
- 2 - Wire + rod to connect to the mass spring assembly.
- 3 – Grooved pulley + wire guide
- 4 – Tube with base
- 5 - A 150 g mass, \varnothing = 25 mm (to study the effect of damping in water)

Supplements (not supplied)

- Polynux® Nut
- Square-section rods 15 mm
- Modumontage® support legs
- R15F5/2A power supply

2 Technical characteristics

This assembly, just like the vertical elastic pendulum, was designed for a video study.

The pulley is equipped with a groove and a guide wire to maintain the wire at resonance.

The mass provided in this assembly has the same mass as the mass m_1 provided with the “vertical elastic pendulum” but a different diameter to enable the study of the effect of damping on the amplitude of oscillations.

The tube is used to:

- Study the fluid damping
- Guide the mass at resonance

The eccentric motor drive has a colour reference mark to easily measure the excitation period by video and compare it with that of the resonator.

3 Experimentation: Study of mechanical resonance

3.1 The educational objectives of the experiment

The phenomenon of resonance must be tackled by an experimental presentation. One must know the conditions when the phenomenon of resonance occurs and in particular also know that resonance occurs when the period of the exciter is close to that of the resonator's natural frequency.

3.2 Performing the experiment

❖ Equipment used

- Vertical elastic pendulum
- The excitation system for vertical elastic pendulum
- The Cinéris acquisition and video processing software (or Génériss 5 plus)
- ToUCam Philips webcam
- 4 Modumontage[®] legs
- 2 lock-grip pliers
- 2 square-section rods 15 mm $L = 1000$ mm
- 1 square-section rod 15 mm $L = 400$ mm
- 1 polynux[®] tightening nut
- R15F5/2A power supply

❖ **Setting up and adjustments**

Connect the webcam to the computer

Open Cinéris

Prepare the set-up (see photo on the 1st page)

Immerse the mass in a tube filled with water.

Acquisition of the movement is done with a webcam and the Cinéris software

Note: Water acts as a *damper* for the system and limits the amplitude of the movement during resonance. The level must be adjusted well so that the mass is always in water.


❖ **The experiment**

Vary the frequency of rotation of the eccentric drive using the power supply, therefore the excitation frequency, and observe the movement of the resonator.

The phenomenon of resonance appears clearly and can be studied on the recorded curves.

3.3 The records obtained and their processing

The recording is done over a duration of 30 s (600 images). On the recording, the experimental points are sufficiently close to be able to connect them without obtaining a curve with a *broken line*.

The student can:

Measure the period of the movement and its amplitude, using the *Cinéris* pointer

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Model to obtain the equations of curves, etc.

Furthermore, thanks to the colour reference mark on the motor, it is also possible to obtain the excitation curve and compare the exciter's period with that of the resonator.

4 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries