Ultrasons

Ultrasounds

**Ref :
222 047**

Français – p 1

English – p 5

Version : 8009

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
1 Description

- | | |
|------------------------|----------------------|
| 1. Wooden sounding-box | 4. Bridges |
| 2. Opening 5 x 5 cm | 5. Pulley on support |
| 3. Tension bolt | 6. Metal strings |

The sound level meter consists of a wooden sounding-box (1) on which are fixed 2 metal strings of different diameters (6) and a nylon string. The tension in the 2 metal strings is created by the tension bolt (3). The sounding-box has a screen-printed marking from 0 to 1 m on its top face.

Two bridges (3) are placed at each end of the marking. The third bridge is positioned in this area to define 2 portions of the string.

There is an opening (5 x 5 cm) at one end to place the microphone inside the sounding-box in order to study the sound emitted during the vibration.

Delivered with a nylon string fitted with a hook at one end.

2 Objectives

- To measure the position of nodes and antinodes in function of the frequency of vibration of a string
- To generate sound
- To study of the vibration of a metal string through which an AC current passes close to a magnet
- To determine the velocity of propagation of a wave


3 Technical characteristics

- Graduation: 0 to 1 m with a precision of 1 mm
- Effective length: 1 m
- Overall dimensions: 120 x 12 x 15 cm
- Mass: 3,8 kg
- Guarantee: 2 years

4 Experiments

4.1 Study of the vibration of a plucked string

Place two bridges on the device at the 0 and 100 markings of the two rulers. These bridges limit the string length to 1 m, see figure below.

Tighten a metal string placed on 2 bridges with a tension bolt until the string generates a sound when plucked.

Place the third bridge at the 66.6 cm marking; it defines two portions of the string with same tension but with a length of one double than the other.

Place paper riders (not supplied) on the 16.6 cm, 33.3 cm and 50 cm markings on the longer portion of the string.

Make the shorter portion of the string vibrate by plucking with the nail or by bowing (bow not supplied): The longer portion begins to vibrate through resonance. We observe that the paper riders on the 16.6 cm and 50 cm markings are thrown off whereas the 33.3 cm marking remains in place.

We thus show that the longer part of the string vibrates in two zones, at its first harmonic and at the same frequency as the first portion of the string.

We can repeat by placing the third bridge at the 75 cm marking to show that the longer part of the string vibrates this time in three zones, at its second harmonic, etc ...

Theoretical values

The modes of vibration are defined by the relation $2.L = n. \lambda$. The plucked portion of the string l_c determines the wavelength of vibration that is $\lambda = 2 \times l_c$.

For $l_c = 33.3$ cm, $\lambda = 2 \times 33.3 = 66.6$ cm, which gives the number of nodes on the total length of the string: $n = 2 \times 66.6 / 66.6$ that is $n = 2$ nodes.

For $l_c = 25$ cm, we thus have $\lambda = 2 \times 25 = 50$ cm, which gives the number of nodes on the total length of the string: $n = 2 \times 75 / 50$ that is $n = 3$ nodes, etc.


4.2 Velocity of propagation a wave generated by LFG

Additional equipment required (not supplied):

- 2MHz amplified generator
- Rheostat $33\ \Omega$
- U-shaped magnet
- Insulated alligator clips

The sound level meter is used over a maximum length of 1 m with a bridge on each end of the metal string (connected to the circuit by two alligator clips).

The U-shaped magnet is placed at the 25 mm marking. An AC voltage of 140 Hz frequency is applied between the 2 bridges: The string vibrates in 2 zones with a node at the 50 mm marking, see the diagram below.

The modes of vibration are defined by the relation $2L = n\lambda$.

The mode of vibration observed corresponds to 2 (n = number of zones), where $\lambda = L = 1\text{ m}$.

We also have the relation $\rightarrow v = \lambda.f$, which gives : $\rightarrow \mathbf{v = 140\text{ m/s}}$.

4.3 Measuring the frequency of the sound generated by a string

Additional equipment required (not supplied):

- 2 MHz amplified generator
- 600Ohms Microphone
- Loudspeaker
- 50MHz Oscilloscope

Method 1:

Connect a loudspeaker to a LFG with a digital display of the frequency.

By listening, find unison between the sound generated by the loudspeaker and the sound generated by the string being studied.

Method 2:

Place a microphone in the sounding-box of the sound-level meter.

Connect the loudspeaker to a storage oscilloscope to measure the periodicity of the sound captured. Deduce the frequency from this.

4.4 Calculating the tension of a string in relation to the frequency of vibration

A string attached at both ends always presents a vibration node at each end and a certain number of intermediary nodes. This system of standing waves is represented by an integer number of zones spread over the length of the string.

If we see k zones, the length of each one will be $\lambda/2$, the total length of the string L is given by the expression: $L = k \cdot \lambda / 2$

In addition, we have: $\lambda = v / \gamma$ and $v = \sqrt{F/\mu}$

Which leads to: $L = k \cdot \sqrt{F/\mu} / 2 \gamma$

With

- γ : frequency in Hz
- v: velocity of transverse waves in m/s
- F: tension of the string, in N
- L: length of the string, in m
- k: integer number of zones
- μ : linear density of the string, in g/m

Example:

A 1 m string having a total mass of 5 g vibrating only in one zone generates a sound having a frequency of 130.5 Hz. Its tension F is expressed by the following relation:

$$F = 4 \cdot \mu \cdot L^2 \cdot \gamma^2 / k^2$$

That is $F = 4 \cdot \mu \cdot L^2 \cdot \gamma^2$ with $k = 1$

That is $F = 345 \text{ N}$

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
0825 563 563

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries