# Ultrasons

## *Ultrasounds*

Réf :  
**222 082**

Français – p 1

English – p 12

Version : 8110

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## 1 Description

These instructions describe the experiments that can be performed with the Moduson ultrasound modular system. The following is the optimum composition:

- 1 transmitter-telemeter, part no. 222 082, powered by 12V power supply part. 281 512 or a 15 V power supply part. 281 486. The frequency of the transmitter can be adjusted to about 40 kHz.  
The transmission of ultrasonic waves is triggered by an On/Off switch. A reversing switch helps obtain either a continuous transmission or a transmission in bursts. A second reversing switch helps select the duration of the bursts:  
Long bursts: duration of 6 ms, with a gap of 4 ms between bursts  
Short bursts: duration of 1 ms, with a gap of 4 ms between bursts  
This function is used to measure distances from 0 to 30 cm using the sonar and scan techniques. Actually, this technique requires that the echo reflected by the object does not get mixed up with the burst emitted, which requires the use of short bursts.
- 2 identical ultrasound receivers (Part no.222 078) equipped with an output designed to be connected to the vertical input of an oscilloscope.
- 1 dual track (Part no. 223 005) fitted with a slot for the transmitter and two guide rails designed for the receivers.
- 3 shielded cables for connection (with BNC socket and banana plugs) (Part no. 283 328).

A BNC output provides the signal emitted by the transmitter.

The following is the minimum composition:

- 1 transmitter-telemeter (Part no. 222 077)
- 1 receiver (Part no. 222 028 or 222 078)
- 1 dual track (Part no. 223 005)
- 1 shielded cable for connection (Part no. 283 328)

## 2 Operating instructions

1. The use of a shielded cable to connect the receiver part no. 222 028 or 222 078 to an oscilloscope is essential to obtain a stable and interference free measurement signal.
2. The transmitter-telemeter is fitted with an automatic rearming fuse; in case of a connection leading to a short circuit, the device goes into standby mode and does not emit any signals. In this case, switch off the power supply to the device, rectify the short circuit then wait for a few seconds and switch it back on.

## 3 Objective

1. To demonstrate the vibratory nature of an ultrasound wave
2. To measure the velocity, period and frequency of an ultrasound wave
3. To observe the existence of points vibrating in phase
4. To verify the relation  $\lambda = V/N$
5. To measure a distance using ultrasound waves

## **4 After-sales service**

The device is under a 2-year guarantee, it must be sent back to our workshops.  
For any repairs, adjustments or spare parts please contact:

**JEULIN – TECHNICAL SUPPORT**  
468 rue Jacques Monod  
CS 21900  
27019 EVREUX CEDEX FRANCE

**+33 (0)2 32 29 40 50**

## Phase velocity of sound with a two receivers

### 1 Demonstrating the vibratory nature of an ultrasound wave

This experiment requires a transmitter and a receiver. Connect the transmitter to the power supply and the receiver, which is placed on the axis of the transmitter to the vertical input of an oscilloscope (adjustment of the level: 100 mV/div and a scan time of: 5  $\mu$ s/div). The On/Off switch is placed in the On position and the Continuous – Bursts reversing switch in Continuous. Adjust the transmitter's adjustment button to obtain the maximum amplitude of the signal detected by the receiver.


We obtain an almost sinusoidal signal. The measurement of the length corresponding to one period on the screen gives the duration of the period from which the frequency of the ultrasound wave is deduced. We will note that the signal which powers the transmitter is a square wave and of a greater amplitude (connect the signal output to the second vertical output of the oscilloscope, level: 5 V/div).

### 2 Measuring the velocity of an ultrasound wave

#### 2.1 Equipment

The experiment requires all components of the assembly. The transmitter and the two receivers are placed on the rails: The transmitter in its slot on the left and the receivers A and B in the guide rails are positioned in such a way that their pointers are at the zero markings. Connect the transmitter to the power supply and the receivers to the vertical inputs of a dual-trace oscilloscope (level on each channel: 100 mV/div and a scan time of: 1 ms/div).

#### 2.2 The experiment

The On/Off button is placed in the On position and the Continuous – Bursts reversing switch in Bursts. The Long bursts – Short bursts reversing switch is placed in the Long bursts position. The two receivers are at the start at the same distance from the receiver.


Adjust the frequency adjustment button of the transmitter to obtain maximum amplitude of bursts detected by the two receivers on the oscilloscope screen.

The oscilloscope screen looks like the diagram opposite.


The bursts are received at regular intervals. They can almost be superimposed by a vertical translation of the curves.

The entire device is shown below:

## 2.3 The measurements

As soon as one of the receivers is moved (for example the receiver A to the right), the corresponding signal also moves to the right of the time axis (provided that the synchronisation channel is properly selected). The receiver A which was moved, receives the emitted signal with a delay  $t$  which depends on the distance  $d$  separating the receiver B, as shown in the diagram opposite.

The measurement of time  $t$  on the oscilloscope screen and distance  $d$  on the graduated ruler on the rails gives a direct measurement of the velocity of sound in air:  $V = d/t$ .

Measurement of the duration  $t$  is better if we select an angle of the burst signal waveform (such as the point M of the figure above) and if we increase the gain of the signal corresponding to the receiver A that moves away from the transmitter to compensate for the *attenuation of the signal because of the distance*. We can also increase the horizontal sweep speed and select a lower adjustment duration (0.5 ms/div).

## 2.4 The results


For a single measurement, if a time  $t = 1 \text{ ms}$  is read (that is one division on the screen with a selected measuring range), the distance is  $d = 34 \text{ cm}$ , which leads to a velocity of  $V = 340 \text{ m.s}^{-1}$ .

It is also possible to measure some values of the time  $t$  (by using the markings on the time axis of the oscilloscope) and obtain values of corresponding distances  $d$ :

<b>t: time (in ms)</b>	0.2	0.4	0.6	0.8	1
<b>d: distance (in cm)</b>	<b>7</b>	<b>14</b>	<b>20.5</b>	<b>27.5</b>	<b>34</b>
<b>Velocity (in m/s)</b>	<b>350</b>	<b>350</b>	<b>341</b>	<b>343</b>	<b>340</b>

This helps conclude that the velocity of ultrasound waves is constant.

### 3 Existence of points in phase and out of phase

This experiment is designed to show the existence of a same vibratory state at various points in space, in the transmission zone of an ultrasound wave.

The experiment requires the entire device. The transmitter and the two receivers are placed on rails: The transmitter in its slot on the left and the receivers A and B in the guide rails are positioned using their pointers placed at the zero markings. Connect the transmitter to the power supply and the receivers to the vertical inputs of a dual-trace oscilloscope (level on each channel: 100 mV/div and time base: 5  $\mu$ s/div). The On/Off button is placed in the On position and the Continuous – Bursts reversing switch in Continuous.

Adjust the frequency-setting button of the transmitter to obtain maximum amplitude of signals detected by the receiver on the oscilloscope screen.

After having rectified the zero position of one of the receivers if necessary, the two signals detected have the shape on the oscilloscope screen as shown below: they can be superimposed.

The points located on the sensors of the receivers are in the same vibratory state: *they are in phase*.

We move the receiver A slowly to the right; the corresponding signal also moves to the right of the screen and it is possible for multiple positions of A to obtain it in phase with the one detected by B.

### 4 Verifying the relation $\lambda = V/N$

The previous experiment is followed by this one. It introduces the definition of wavelength: distance separating two consecutive points vibrating in phase.

Move the receiver A while keeping B fixed. Measure the distance  $d$  separating A from B when the signals on the screen are  $n$  times in a row in coincidence. We deduce  $\lambda = d/n$ . We obtain a value close to **8.5 mm**.

The known values of **V** and **N** result in the theoretical value of  $\lambda = V/N$  which is very close to the previous measurement.


## Phase velocity of sound with a single receiver

### 1 Phase velocity of sound with a single receiver

#### 1.1 Experiment

- Set the transmitter on the «bursts» and «long bursts» positions as far away as possible from the source.
- View the receiver signal on channel A of the oscilloscope, and the "source" signal on channel B. (time base: 2 ms / Div).

Synchronise the oscilloscope on the "source" signal.

#### 1.2 Note

- We can observe the increase in response on the receiver signal. The offset with the signal sent to the transmitter is not significant. On the other hand, the significant reading is that of the displacement of the start of this signal on the screen (measurable lag or lead) in function of the receiver's displacement.
- Move the receiver by taking it closer to the source and note the displacement of the signal on the oscilloscope.
- For a receiver displacement of 170 mm, we find a signal displacement of 0.5 ms on the oscilloscope.
- For a receiver displacement of 340 mm, we find a signal displacement of 1 ms. This result corresponds to a phase velocity of 340 m / s.

#### 1.3 Lighter experiment

By heating the air between the transmitter and the receiver with a lighter or a lighted bulb, we observe the signal move: the phase velocity of sound increases with temperature.

$$v = \left( \sqrt{\frac{\gamma \cdot R \cdot T}{M}} \right)$$


## 2 Measuring the wave length of a single receiver

### 2.1 Experiment

- Set the transmitter to the "continuous" position
- View the transmitter signal on channel A of the oscilloscope, and the receiver signal on channel B. (time base: 10  $\mu$ s/div)

### 2.2 Synchronise the oscilloscope on the transmitter signal

- Locate the initial position of the receiver.
- Move the receiver away (or closer to) from the source and observe the scrolling of 10 periods on the oscilloscope.

### 2.3 Note

The relation  $\lambda = C.T$  helps calculate the phase velocity  $C$ . The period  $T$  is measured with the oscilloscope.

with  $T = 2.5 \times 10^{-5}$  s.


and  $\lambda = 8.29$  mm at 20° C.

we get  $C = 331$  ms $^{-1}$ .

## 3 Measuring distance using ultrasound waves

### 3.1 Objectives

To know the position of a small character enclosed in an opaque box. For this we will use ultrasound waves employed in the sonar and scan techniques.

The On/Off switch is placed in the On position and the Continuous – Bursts reversing switch on Bursts. The Long bursts – Short bursts reversing switch is placed in the Short bursts position.

### 3.2 Detecting ultrasound waves

The following experiment is performed on the desk. Bursts of ultrasound waves are emitted by the transmitter E and received by the receiver R. To determine the distance  $d$  that separates E from R, we study it on the oscilloscope:


The oscillogram 1 corresponds to E.

The oscillogram 2 corresponds to R.

The horizontal adjustment is 2 ms per division (one horizontal division is equal to a duration of 2 ms.).

Draw the diagram of the experiment.


1. Find the velocity of propagation of sound or ultrasound waves in air.
2. What is the duration that separates the start of the transmission from the start of the reception?
3. Calculate the distance  $d$ , then verify the result by measuring  $d$  with a graduated ruler.

**Partial answer:**

1.

Channel 1

2.  $v = 340 \text{ m/s}$

3. 0.77 ms;  $d = 26 \text{ cm}$ ; 25 cm with a graduated ruler (the distance of 1 cm between the sensitive tip and the E or R casing must be taken into consideration).

### 3.3 What happens when ultrasound waves encounter an obstacle?

Write an experimental protocol that can be used to know what happens when ultrasound waves encounter:

- A cardboard sheet
- A facial tissue paper
- A plastic character
- A hand.


The experiments are then performed on the desk by the teacher or the student. Note down the results in the notebook.

**Partial answer:**


Place the screen in front of the transmitter; the receiver is placed opposite the screen. Observe the oscillogram corresponding to the transmitter.

Results: A cardboard sheet, a plastic character or a hand reflect ultrasound waves; a facial tissue paper (open the tissue and keep only one layer) reflects and transmits ultrasound waves.

### 3.4 At what distance from point A is the small character standing?

1. Draw a diagram of an experiment using ultrasound waves to measure  $d$ .
2. Here are the results obtained:
  - The oscillogram 1 corresponds to E
  - The oscillogram 2 corresponds to R
  - The horizontal adjustment is 2 ms per division.


Determine the time separating a transmission from reception, then calculate the distance  $d$  and check the result by measuring  $d$  with a graduated ruler.

3. Find a practical application of this experiment.

**Partial answer:**

1. top view


2. 1.57 ms;  $d = (340 \times 1.57) / 2 = 267 \text{ mm} = 27 \text{ cm}$ ; 26 cm with a ruler.
3. Sonar.

### 3.5 What is the depth of the box?


The front face of a cardboard box was removed and replaced by a facial tissue paper. We propose to determine the depth  $l$  of the box using ultrasound waves.

Top view diagram of the experiment:

Here are the oscillograms obtained:

- Oscillogram 1 corresponds to the transmitter
- Oscillogram 2 corresponds to the receiver
- The horizontal adjustment is 2 ms per division.

1. Observe oscillogram 2 and explain what the first and second receptions correspond to.
2. Determine from this oscillogram the distances  $d_1$  and  $d_2$ , then find the depth  $l$  of the box.

Check the result using a graduated ruler.

**Partial answer:**


1. to  $t_1 = 0.61$  ms corresponds the reception of ultrasound waves reflected by the tissue paper, to  $t_2 = 1.52$  ms corresponds the reception of ultrasound waves transmitted through the tissue paper, then reflected by the bottom of the box.
2.  $d_1 = (340 \times 0.61) / 2 = 104$  mm;  $d_2 = (340 \times 1.52) / 2 = 258$  mm;  $l = 154$  mm ; 160mm with a ruler.

### 3.6 Where is the small character hidden in the box?

Place the character in the box.

1. After processing the oscillograms given below, draw the box by respecting its dimensions and place the character correctly inside.

- Oscillogram 1 corresponds to the transmitter,
- Oscillogram 2 corresponds to the receiver,
- The horizontal adjustment is 2 ms per division.

2. Find a practical application of this experiment.

**Partial answer:**

1. First echo corresponds to the front face of the box :  
 $t_1 = 0.59 \text{ ms}$ ;  $d_1 = 100 \text{ mm}$   
Second echo corresponds to the character :  
 $t = 1.01 \text{ ms}$ ;  $d = 172 \text{ mm}$   
Third echo corresponds to the rear face of the box :  
 $t_2 = 1.51 \text{ ms}$ ;  $d_2 = 256 \text{ mm}$   
The character is located 72 mm from the front face of the box.
2. Ultrasound scan.

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

### Free service\*

**+33 2 32 29 40 50\*\***

\* Call cost not included.  
\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux