Ultrasons

Ultrasounds

Réf :
223 008

Français – p 1

English – p 32

Version : 5103

Plateau Moduson

Moduson table

Table des matières

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Contents

1	Description	32
1.1	Moduson table.....	32
2	Installation.....	32
2.1	Required additional equipment	32
2.2	Recommended additional equipment.....	33
2.3	Setup	33
3	Operation	34
3.1	Ultrasonic reflection and absorption	34
3.2	Ultrasonic interference phenomena.....	34
4	Theoretical study of interference phenomena.....	35
5	After-Sales Service.....	38
EXPERIMENTS		40
1	Directivity of ultrasonic waves.....	41
1.1	Step 1: Setup.....	41
1.2	Step 2: Protocol.....	42
1.3	Step 3: Exploitation of the results.....	43
2	Reflection of ultrasonic waves.....	43
2.1	Step 1: Setup.....	43
2.2	Step 2: Protocol.....	44
2.3	Step 3: Exploitation of the results.....	45
3	Interference of ultrasonic waves	46
3.1	Step 1: Setup.....	46
3.1	Step 2: Protocol.....	47
3.2	Step 3: Exploitation of the results.....	48

LAB WORK SHEETS	52
1 Directivity of ultrasonic waves	53
1.1 Step 1: Setup.....	53
1.2 Step 2: Protocol.....	54
1.3 Step 3: Exploitation of the results.....	55
2 Reflection of ultrasonic waves.....	56
2.1 Step 1: Setup.....	56
2.2 Step 2: Protocol.....	57
2.3 Step 3: Exploitation of the results.....	58
3 Interference of ultrasonic waves	58
3.1 Step 1: Setup.....	59
3.1 Step 2: Protocol.....	60
3.2 Step 3: Exploitation of the results.....	61

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1 Description


1.1 Moduson table

This is a moulded PVC tray with two parts that fit into each other:

- The first part functions as a support for the ultrasonic transmitters: a recess allows the accurate positioning of the transmitters at various locations. These positions are identified by silkscreen-printed marks and they correspond to distances of 2, 4 and 8 cm between the two transmitters, respectively, centred with respect to the midline of the instrument.

In order to maintain these positions stationary, the bottom of the recess and the bottom of the transmitters are covered with an adhesive magnetic sole.

- The second part has two arc-shaped recesses with centre-zero graduations in degrees in which the ultrasonic receiver can be moved. The arc radii are 15cm and 30cm, respectively.

View of the MODUSON system (interference experiment)

2 Installation

2.1 Required additional equipment

- | | |
|---------------------------------------|----------------|
| - MODUSON ultrasonic transmitter | P/N 222 027 |
| - Or MODUSON 2 ultrasonic transmitter | P/N 222 076 |
| - MODUSON ultrasonic transmitter 2 | P/N 223 009 |
| - MODUSON ultrasonic receiver | P/N 222 028 |
| - Or MODUSON 2 ultrasonic receiver | P/N 222 078 |
| | |
| - Power supply model Evo 2 R15F5 | P/N 281 488 |
| - Oscilloscope | P/N 291 125 |
| - | Or P/N 291 123 |
| - | Or P/N 291 105 |

2.2 Recommended additional equipment

Banana-to-BNC adapter cable


P/N 283 328

2.3 Setup

❖ A few changes must be made prior to installation:

Adhere the adhesive magnetic sole on the bottom of the MODUSON transmitter (P/N 222 027 or P/N 222 076) and place the self-adhesive label (12 mm x 63 mm) representing a white line on the front and in the upper part of the transmitter.

On one of the MODUSON receivers (P/N 222 028 or P/N 222 078), adhere the self-adhesive label (8 mm x 32 mm) representing a white line on the front and in the lower part of the receiver.


In both cases, these labels are indexes for locating the position of the transmitter or the position of the receiver inserted in their recesses on the moulded PVC base.

❖ MODUSON requires the use of a 15V stabilised power supply.

Note: A built-in protection circuit prevents any instrument damage in case of wrong polarity connection of the power supply.

❖ MODUSON Complement can be used to perform two experiments:

In fifth form for the study of reflection. In this experiment, the mirror with a magnetic base is used together with one transmitter and one receiver. The transmitter (P/N 222 027 or 222 076) can be used directly but its base is too large and does not fit into the recess in the moulded PVC table. It is preferable to use the complementary transmitter 2, which has a suitable size, at the expense of a slightly bulkier setup.

In scientific upper sixth forms, the study of interference phenomena requires the simultaneous use of the two (main and secondary) transmitters. The transmitter 2 is connected by means of a shielded cable to the BNC output of the main transmitter that delivers a common and synchronous signal for the two transmitters. See the figure on page 1.

3 Operation

3.1 Ultrasonic reflection and absorption

These experiments are proposed as a complement to the various experiments that can be performed with the MODUSON system.

The main transmitter is used in the *Continuous* position. It must be set beforehand at resonance via the *Frequency* knob so that a signal with the highest possible amplitude is transmitted. The BNC connector of the transmitter is connected to the MODUSON transmitter 2 placed inside the circular recess on the moulded table. This transmitter sends an ultrasonic wave towards a metallic mirror located at the middle of the recess.

The MODUSON receiver is also placed inside the circular recess. It senses the reflected signal. The maximum amplitude of the received signal is found by observation on the screen of the oscilloscope to which the receiver is connected.

The setup demonstrates or verifies the equality of the angle of incidence and of the angle of reflection.


If the mirror is turned over, the reflective surface is replaced by plastic foam. Then the incident signal is entirely absorbed.

Note: The reflection of ultrasonic waves on the mirror goes along with a diffraction of the reflected signal. Owing to the small width of the mirror (4 cm), this phenomenon is not inconvenient: the centre peak is sharp enough and the secondary peaks are sufficiently weak and away from the main peak so that measurements can be easily performed in fifth forms.

3.2 Ultrasonic interference phenomena

The main transmitter (P/N 222 027 or 222 076) is assumed to be set at its maximum transmission level. It is connected to the MODUSON transmitter 2.

- The two transmitters connected to the same signal generator are synchronous sources. However, the following adjustments are required to give them an identical amplitude:
 1. Install the two transmitters in the magnetic recess on the moulded PVC table. Place them at a distance of e.g. 2 cm from each other (smallest distance). Place a receiver in one of the arc-shaped recesses opposite the zero-degree angle value (then the receiver is located on the median of the segment joining the two receivers and it is equidistant between them) and connect it to the vertical input of an oscilloscope.
 2. Disconnect the transmitter 2 and observe the signal transmitted by the main transmitter and sensed by the receiver. Record the amplitude value. Connect the transmitter 2 and switch the main transmitter *OFF* (in this position, the transmitter member is off but the generator is not and continues delivering the same signal on the BNC connector, hence to the transmitter 2), then only the transmitter 2 is active. Using the transmitter potentiometer, trim the amplitude of the received signal so that the amplitude value is identical to the previous reading.
 3. Switch the main transmitter *ON*: then the two transmitters operate simultaneously and they have the same amplitude. This produces a peak that is sensed by the receiver (a slight position adjustment by 1–2 degrees may be necessary to obtain this peak).

Note: Depending on the main transmitter and transmitter 2 used, the signal received at point 0 by the receiver when the two transmitters have identical amplitude settings may happen to yield a zero minimum instead of a peak. Then the two transmitters vibrate in opposition of phase instead of vibrating in phase as would be expected. To put them in phase, activate the phase selector switch.

- Once the receiver has been installed, its displacement inside the arc-shaped recess will produce a succession of peaks and minima having the amplitude of the signal sensed by the receiver and observed on the oscilloscope screen.

This is the series of the antinodes and nodes corresponding to the constructive and destructive interferences that appear in the vicinity of the two synchronous sources used.


- Select the distance a between the 2 transmitters (3 values are marked: $a=2$ cm, $a=4$ cm and $a=8$ cm, but any other value from 2 to 4 cm can be measured by means of a 20-cm rule) and the radius D of the receiver displacement arc (2 radii can be used: $D=15$ cm and $D=30$ cm). The measurement of the angular X-axis values θ of the consecutive minima (or peaks) yields the value i of the interfringe distance and allows the verification of its expression versus a and D with a fairly good accuracy.

$$i = \frac{\lambda \cdot D}{a}$$

4 Theoretical study of interference phenomena

❖ Ultrasonic interference with MODUSON

- Theoretical study
 - The theoretical top-view diagram is as follows. The receiver M moves on the arc of circle with radius $OM_0=D$. The two synchronous sources S_1 and S_2 are spaced by a distance $S_1S_2=a$

For an arbitrary receiver position M , the difference in the courses of the ultrasonic waves from of S_1 and S_2 is $\delta = d_2 - d_1$

d_2 and d_1 are the hypotenuses of two rectangular triangles. The following can be expressed:

$$d_2^2 = hM^2 + hS_2^2 = Ov^2 + (Oh + a/2)^2 = D^2 \cos^2(\theta) + (D \sin(\theta) + a/2)^2 = D^2 + a^2/4 + aD \sin(\theta)$$

$$d_1^2 = hM^2 + hS_1^2 = Ov^2 + (Oh - a/2)^2 = D^2 \cos^2(\theta) + (D \sin(\theta) - a/2)^2 = D^2 + a^2/4 - aD \sin(\theta)$$

$$d_2^2 - d_1^2 = (d_2 - d_1)(d_2 + d_1) = \delta(d_2 + d_1) = 2 aD \sin(\theta)$$

$$\text{or } \delta = 2 aD \sin(\theta) / (d_2 + d_1)$$


or else:

$$\delta = \frac{2aD \sin(\theta)}{\sqrt{D^2 + \frac{a^2}{4} + aD \sin(\theta)} + \sqrt{D^2 + \frac{a^2}{4} - aD \sin(\theta)}}$$

- Approximations
 - The computerised study of the denominator $\Phi(\theta)$ of δ by means of ACTILAB or of a spreadsheet, e.g. EXCEL, yields the denominator variations in the experimental conditions.
 - varies at most between -60 deg. and $+6$ deg. ($-\pi/3$ and $+\pi/3$)


Note: The graphical study can also be performed by using a graphic calculator.

In worst-case conditions, the values of a and D are $a=8$ cm and $D=15$ cm. Our working domain is comprised between these limits and in this case. The representative curve of $\Phi(\theta)$ versus θ yields the following graph:

- For the value $2D=30$ cm, the relative deviation $(\Phi(\theta)-2D)/2D$ is $1.05/30$ or 3.5% maximum at the angular abscissa $\theta=0$. The relative deviation drops to approx. $0.8/30$ or 2.7% for $\theta=+\pi/6$ or $-\pi/6$ and to $0.3/30$ or 1% for $\theta=+\pi/3$ or $-\pi/3$.

In worst-case conditions, the values of a and D are, $a=2$ cm and $D=30$ cm. Then the representative curve of $\Phi(\theta)$ versus θ yields the following graph:


- For the value $2D=60$ cm, the relative deviation is $0.033/60$ or 0.055% maximum at the angular abscissa $\theta=0$.
- The conclusion is that $2D$ can be reasonably substituted for $\Phi(\theta)$ with an error that remains acceptably low.

$$\delta \approx a \sin(\theta)$$

- Expression of the interfringe distance
 - Constructive interferences are obtained when the two waves are in phase, i.e. when $\delta = k \cdot \lambda$ (λ is the wavelength of the ultrasonic signal, here $\lambda = 0.84$ cm).
 - Therefore we have $k \cdot \lambda \approx a \sin(\theta)$ or $\sin(\theta) \approx \frac{k \cdot \lambda}{a}$. Since $\frac{\lambda}{a}$ is constant for a given value of a , $\sin(\theta)$ is an integral multiple of $\frac{\lambda}{a}$.
 - For $a = 4$ cm, we obtain $\frac{\lambda}{a} = \frac{0.84}{4} = 0.21$ and on a $\sin(\theta) = 0.21 \cdot k$.
 - Since $-\pi/3 \leq \theta \leq +\pi/3$, we have $-0.866 \leq \sin(\theta) \leq +0.866$, or
 - $-0.866 \leq 0.21 \cdot k \leq +0.866$ or else $-4.12 \leq k \leq +4.12$,
 - which gives the possible values for k and the corresponding amplitude peaks:
 $k = -4, k = -3, k = -1, k = 0, k = 1, k = 2, k = 3$ and $k = 4$.

On the diagram below, the peaks are obtained by a simple trace when Oh meets the following condition:

Since $\sin(\theta) = Mv/OM = Oh/D$, then we must have: $Oh = \sin(\theta) \cdot D \approx \frac{k \cdot \lambda}{a} \cdot D$.

If we assume $D = 15$ cm, $a = 4$ cm and $\lambda = 0.84$ cm, the possible values for k are: $Oh = 3.15 \cdot k$, or
 $Oh_{-4} = -12.6$ cm, $Oh_{-3} = -9.45$ cm, $Oh_{-2} = -6.3$ cm, $Oh_{-1} = -3.15$ cm, $Oh_0 = 0$.
 $Oh_1 = 3.15$ cm, $Oh_2 = 6.3$ cm, $Oh_3 = 9.45$ cm, $Oh_4 = 12.6$ cm.

These values allow the plotting of an array of straight lines parallel to OM_0 and equidistant by $\frac{\lambda D}{a}$. Their intersection with the arc of circle gives the locations of the peaks.

- Case of small angles

When $\sin(\theta)$ is small, i.e. a is large enough, we observe that the arcs M_0M_1 , M_1M_2 , M_0M_{-1} , $M_{-1}M_{-2}$ are equal if the deviation with respect to the position M_0 is not too great. In practice a value $a = 8$ cm can be used.

Let θ_i , which is the measure of the angle between two successive peaks. We have:

$$\theta_i \approx \sin(\theta_i) = \lambda/a = 0.105 \text{ or } \theta_i = 6^\circ.$$

This is well verified by the experiment.

Similar reasoning can be applied to the amplitude minima.

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
468 Rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

+33 (0)2 32 29 40 50

EXPERIMENTS

- ① The first part of this compendium of experiments is for teacher use only. It gathers a few simple experiments performed by using the MODUSON Complement. The numerical results contained therein are derived from actual measurements. However these measurements were taken in arbitrary conditions and they may differ from the experimenter's personal measurements, therefore the following measurements are not examples.
The first experiments are meant for fifth form students.
The last experiment on interferences is meant for sixth form science students.
- ② The second part of the compendium carries on with the experiments of the first part in the form of lab work sheets that can be used as such by the teacher in his/her class or easily modified. Questions are often asked in the student's text intended: these questions often correspond to the notes relevant to the observations, the experimental results and the conclusions of the first part.
- ③ Precautions of use: The preliminary setting of the transmitters consists in obtaining the signal with the highest possible amplitude received by the receiver. This must be done prior to any experiment.


1 Directivity of ultrasonic waves

Introduction: The purpose of the experiment is to demonstrate the directional effect of ultrasonic waves and to define the notion of « sound radiation ».

1.1 Step 1: Setup

① The system is set up as follows:

The transmitter is connected to a 15V power source. It is set at the transmission maximum and powered up in *Continuous* mode. It is positioned at the centre O of the arc-shaped recess and it transmits ultrasonic waves towards the graduation zero.


The receiver is positioned in the circular recess closest to the transmitter (radius = 15 cm) in which it can be moved. It is connected by means of a shielded cable to the vertical input of the oscilloscope.

② Material :

- MODUSON main ultrasonic transmitter
- One ultrasonic receiver
- MODUSON moulded PVC table (two elements)
- 0–15V stabilised variable power supply
- Oscilloscope
- Banana-to-BNC adapter cable

③ Setup schematic


1.2 Step 2: Protocol

Directivity of ultrasonic waves

① Determination of the reception peak:

The reception peak is obtained by finding the receiver position on the arc of circle that yields the highest amplitude for the received signal.

The oscilloscope settings are 0.1 V/div for vertical amplification and 0.1 msec/div for horizontal sweep. The periodical signal is then very narrow and its amplitude is better readable. The amplitude is 0.3 V.

The amplitude peak of the signal sensed by the receiver corresponds to the zero graduation.

② Sound transmission cone:


The ultrasonic waves are emitted by the transmitter in a zone to be determined.

Record the value a of the received signal amplitude for a variable value θ of the angle defining the position of the ultrasonic receiver (θ is negative on the left of 0 and positive on the right).

The following results are obtained:

a in V	0,09	0,12	0,15	0,19	0,23	0,27	0,28	0,3	0,3	0,28	0,26	0,22	0,20	0,17	0,13	0,10	0,08
θ in °	-40	-35	-30	-25	-20	-15	-10	-5	0	5	10	15	20	25	30	35	40

③ The corresponding graph is shown below:

④ Approximation of the measurement :

Determine the minimum variation of the angle for which the amplitude observed on the oscilloscope screen is the highest.

Observe that this value is about 5 degrees.

Therefore the determination of the amplitude peak is rather inaccurate.

❖ Conclusion :

We find that ultrasonic waves are rather directive (however much less so than light waves!)

The ultrasonic transmission zone (assuming that it is symmetrical about the transmitter axis), for an amplitude value decreased by half its maximum value, is a cone having an angle close to 60 degrees (angle α_{Ob} on the diagram).

Note: The experiment can be repeated by using the second recess for the receiver displacement.

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The so-called « sound radiation » is the direction of transmission that corresponds to the amplitude peak of the signal received by the receiver.

❖ Representation of sound radiation

This radiation is represented by an oriented line with an arrow pointing outwards from the transmitter. Its direction blends with the geometrical axis of the transmitter.

ER represents the sound radiation


2 Reflection of ultrasonic waves

Introduction: The purpose of the experiment is to verify the law of reflection. The influence of the absorption phenomenon on the ultrasonic waves can also be observed.

2.1 Step 1: Setup

❖ The system is set up as follows:

The reflective plane (mirror) fitted with a magnetic sole is positioned at the centre of the arc-shaped recesses. An ultrasonic transmitter and a receiver are positioned as shown on the diagram. The MODUSON main transmitter can be used directly but its base is too large and does not fit into the recess in the moulded PVC table. *It is preferable to use the complementary transmitter 2*, which has a suitable size, at the expense of a slightly bulkier setup. The transmitter 2 is positioned e.g. on the 30 deg. graduation on the left.


The transmitter 2 is connected to the main transmitter. The main transmitter is powered via a 15 V power supply and its switches are set on Continuous and OFF, respectively (then its transmitting member is cut off and only the transmitter 2 is powered).

The receiver is positioned in the circular recess closest to the transmitter (radius = 15 cm) in which it can be moved. It is located in the right portion. It is connected by means of a shielded cable to the vertical input of the oscilloscope.

❖ **Material :**

- MODUSON main ultrasonic transmitter and transmitter 2
- One ultrasonic receiver
- MODUSON moulded PVC table (two elements)
- 0–15V stabilised variable power supply
- Oscilloscope
- Banana-to-BNC adapter cable

❖ **Setup schematic**

2.2 Step 2: Protocol

❖ **Definitions**

On the above schematic, the symbols have the following meanings:

- On is the perpendicular to the mirror
- i is the angle of incidence between the incident sound radiation EO and the perpendicular to the mirror On
- r is the angle of reflection between the reflected sound radiation OR and the perpendicular to the mirror On

❖ **Measurements :**

Different values are assigned to i by moving the transmitter in its guide and the position to be given to the transmitter in order to obtain the direction of the radiation reflected on the mirror is determined.

This direction is obtained by determining the receiver position on the arc of circle that yields the highest amplitude for the received signal.

The oscilloscope settings are 0.1 V/div for vertical amplification and on 0.1 msec/div for horizontal sweep.

To perform correct measurements, proceed as follows:

- position the transmitter on the first value (10 degrees)
- slowly move the receiver from the lowest value (towards 0 degree)
- stop as soon as the first amplitude peak of the received signal is observed
- read the corresponding value
- fine-tune as required.


Note: The presence of other peaks (with lower amplitudes) is observed. This is due to complex phenomena called « diffraction phenomena ».

The results are summarised in the table below:

i in deg.	10	15	20	25	30	40	50	60
r in deg.	11	14	19	27	31	40	51	59

❖ **Plotting of the curve :**

Graphic representation of the function $i=f(r)$

❖ **Observation of the curve :**

The obtained curve is a straight line passing through the origin (within the limits of the domain under study).

❖ **Ultrasonic absorption :**

Complementary experiment.

Turn the mirror over so that its foam-coated side faces the ultrasonic beam from the transmitter.

For a few values of the angle i , try to obtain of the corresponding values of r .

Observe that the signal is no longer reflected: no observable amplitude can be obtained for the signal sensed by the receiver.

2.3 Step 3: Exploitation of the results

❖ **Verification of the law of reflection**

The measurements are rather tricky and they are limited to an approximate verification of the relationship

$$i=r$$

❖ **Ultrasonic absorption**

The ultrasonic waves are absorbed by the foam coating of the support used.


Note: The experiment can be repeated by using the second recess for the displacement of the transmitter or (and) of the receiver.

3 Interference of ultrasonic waves

Introduction: The purpose of the experiment is to conduct an interference experiment with ultrasonic waves and to compare the experimental measurements obtained with the theoretical results derived from the wave model of the phenomenon.

3.1 Step 1: Setup

❖ The system is set up as follows:

The transmitter 2 is positioned on a magnetic support inclined at 45 deg. This special position is necessary to achieve short distances between the two transmitters. Thus they contact each other in their smallest spacing (2 cm).


The transmitter 2 is connected to the main transmitter. The main transmitter is powered via a 15 V power supply and its switches are set on *Continuous* and *OFF*, respectively (its transmitting element is then cut off and only the transmitter 2 is powered).

The receiver is positioned in the circular recess closest to the transmitter (radius = 15 cm) in which it can be moved. It is positioned on the zero graduation. It is connected by means of a shielded cable to the vertical input of the oscilloscope.

❖ **Material :**

- MODUSON main ultrasonic transmitter and transmitter 2
- One ultrasonic receiver
- - MODUSON moulded PVC table (two elements)
- 0–15V stabilised variable power supply
- Oscilloscope
- Banana-to-BNC adapter cable

❖ **Setup schematic**

3.1 Step 2: Protocol

❖ Preliminary setting of the transmitters

The main transmitter is assumed to be set at its maximum transmission level. It is connected to the MODUSON transmitter 2.

The two transmitters connected to the same signal generator are synchronous sources. However, the following adjustments are required to give them an identical amplitude:

1. Install the two transmitters in the magnetic recess on the moulded PVC table. Place them at a distance of e.g. 2 cm from each other (smallest distance). Place a receiver in one of the arc-shaped recesses opposite the zero-degree angle value (then the receiver is located on the median of the segment joining the two receivers and it is equidistant between them) and connect it to the vertical input of an oscilloscope.
2. Disconnect the transmitter 2 and observe the signal transmitted by the main transmitter and sensed by the receiver. Record the amplitude value. Connect the transmitter 2 and switch the main transmitter *OFF* (in this position, the transmitter member is off but the generator is not and continues delivering the same signal on the BNC connector, hence to the transmitter 2), then only the transmitter 2 is active. Using the transmitter potentiometer, trim the amplitude of the received signal so that the amplitude value is identical to the previous reading.
3. Switch the main transmitter *ON*: then the two transmitters operate simultaneously and they have the same amplitude. This produces a peak that is sensed by the receiver (a slight position adjustment by 1–2 degrees may be necessary to obtain this peak).

❖ Measurements:

The oscilloscope settings are 0.1 V/div for vertical amplification and on 0.1 msec/div for horizontal sweep.

- Select the distance a between the 2 transmitters (3 values are marked: $a=2$ cm, $a=4$ cm and $a=8$ cm, but any other value from 2 to 4 cm can be measured by means of a 20-cm rule) and the radius D of the receiver displacement arc (2 radii can be used: $D=15$ cm and $D=30$ cm).
- Let for example the values $a=4$ cm and $D=15$ cm.
- Once the receiver has been installed, its displacement inside the arc-shaped recess will produce a succession of peaks and minima having the amplitude of the signal sensed by the receiver and observed on the oscilloscope screen.
- This is the series of the antinodes and nodes corresponding to the constructive and destructive interferences that appear in the vicinity of the two synchronous sources used.
- To obtain good measurements, proceed as follows:
 - position the transmitter opposite the zero-degree graduation
 - this corresponds to an amplitude peak (trim as required)
 - slowly move the receiver towards the increasing values of the angle θ (e.g. to the right)
 - stop when the first amplitude peak of the received signal is obtained, trim, read and record the corresponding value
 - continue until the extreme receiver position is reached
 - repeat these steps for the values of the angle θ to the left.
- Repeat the above steps for the angles corresponding to the amplitude minima.

Note : The minima are more easily identified than the peaks. However, the minima are not zero and they are less and less marked. The amplitudes of the interfering signals are not equal and the difference between them increases as the distances d_1 and d_2 become different.

❖ The results are summarised in the table below:

The angles are given in algebraic notation relative to the graduation zero.

θ minimum	Not measurable	-32	-18	-6		6	18	32	47
minimum difference			14	12	12	12	14	15	
θ maximum	Not measurable	-38	-24	-12	0	11.5	24	38	54
minimum difference		14	12	12		11.5	12.5	14	16

We observe that for angles in the range of -25 deg. to +25 deg., the variation of the angle θ from one peak to the adjacent peak (and from a minimum to the adjacent minimum) is constant and equal to:

$$\Delta\theta \approx 12^\circ$$

As for the number of peaks and of minima, we have:

- 8 peaks (the 9th one could not be read)
- 7 minima (the 8th one could not be read)

The differences between the peaks and the minima is identical and there are regularly alternated.

Note: The results are similar with other values. For example with $a=8$ cm and $D=30$ cm, the mean difference is $\Delta\theta \approx 6^\circ$

3.2 Step 3: Exploitation of the results

❖ Theoretical study of the interference phenomenon

① The superposition model of the waves from the two synchronous sources is used.

The wavelength is λ .

- At an arbitrary point M, the waves from the sources S_1 and S_2 exhibit a course difference $\delta=d_2-d_1$.
- If the two waves arrive *in phase* at point M, they have constructive a effect and their amplitudes are added: this yields a maximum amplitude point called *internode*.

This condition is rendered by $\delta=k \cdot \lambda$.

- If the two waves arrive *in opposition of phase* at point M, they have a destructive effect and their amplitudes are subtracted: this yields a point with zero amplitude (if the amplitudes are equal) or a low amplitude (equal to the amplitude difference), called *node*.

This condition is rendered by $\delta=(2k+1) \cdot \lambda/2$.

② Calculation of the course difference

In the above setup schematic, the course difference is expressed as $\delta=d_2-d_1$.

d_2 and d_1 are the hypotenuses of two rectangular triangles (MhS_1 and MhS_2).

The following relationships can be expressed:

$$d_2^2 = hM^2 + hS_2^2 = Ov^2 + (Oh + a/2)^2 = D^2 \cos^2(\theta) + (D \sin(\theta) + a/2)^2 = D^2 + a^2/4 + aD \sin(\theta)$$

$$d_1^2 = hM^2 + hS_1^2 = Ov^2 + (Oh - a/2)^2 = D^2 \cos^2(\theta) + (D \sin(\theta) - a/2)^2 = D^2 + a^2/4 - aD \sin(\theta)$$

The difference is written as follows:

$$d_2^2 - d_1^2 = (d_2 - d_1)(d_2 + d_1) = \delta(d_2 + d_1) = 2 aD \sin(\theta)$$

$$\text{or } \delta = 2 aD \sin(\theta) / (d_2 + d_1)$$

or else:


$$\delta = \frac{2aD \sin(\theta)}{\sqrt{D^2 + \frac{a^2}{4} + aD \sin(\theta)} + \sqrt{D^2 + \frac{a^2}{4} - aD \sin(\theta)}}$$

③ Approximations

The computerised study of the denominator $\Phi(\theta)$ of δ by means of ACTILAB yields its variations in the experimental conditions.

θ varies at most in the range from -60° to $+60^\circ$ (from $-\pi/3$ to $+\pi/3$)

For the selected values of a and D , the representative curve of $\Phi(\theta)$ versus θ gives the following graph:

For the value $2D=30$ cm, the relative difference $(\Phi(\theta)-2D)/2D$ is 0.27/30 or 0.9% maximum at the angular abscissa $\theta=0^\circ$. The relative difference drops to about 0.21/30 or 0.7% for $\theta=+\pi/6$ or $-\pi/6$, and it reached 0.6/30 or 0.2% for $\theta=+\pi/3$ or $-\pi/3$.

Therefore, $2D$ can be substituted for $\Phi(\theta)$.

Then the expression of the course difference is:

$$\delta \approx a \sin(\theta)$$

④ Number of peaks

For the peaks, we express $\delta = k \cdot \lambda \approx a \sin(\theta)$ or $\sin(\theta) \approx \frac{k \cdot \lambda}{a}$ ($\lambda=0.84$ cm).

Since $\frac{\lambda}{a}$ is constant for a given value of a , $\sin(\theta)$ is an integral multiple of $\frac{\lambda}{a}$.

For $a=4$ cm, we have $\frac{\lambda}{a} = \frac{0.84}{4} = 0.21$ and we have $\sin(\theta) = 0.21 \cdot k$.

Since $-\pi/3 \leq \theta \leq +\pi/3$, we have $-0.866 \leq \sin(\theta) \leq +0.866$, or

$$-0.866 \leq 0.21 \cdot k \leq +0.866 \text{ or else } -4.12 \leq k \leq +4.12$$

The possible values of k for the amplitude peaks are:

$k=-4, k=-3, k=-2, k=-1, k=0, k=1, k=2, k=3$ and $k=4$.

⑤ Number of minima

For the minima, we have $\delta = (2k + 1) \cdot \frac{\lambda}{2} \approx a \sin(\theta)$ or $\sin(\theta) \approx \frac{(2k + 1) \cdot \lambda}{2a}$

Since $\frac{\lambda}{2a}$ is constant for a given value of a, $\sin(\theta)$ is a integral multiple impair of $\frac{\lambda}{2a}$

For a=4 cm, we have $\frac{\lambda}{2a} = \frac{0,84}{8} = 0,105$ and on a $\sin(\theta) = 0.105 \cdot (2k+1)$.

Since $-\pi/3 \leq \theta \leq +\pi/3$, we have $-0.866 \leq \sin(\theta) \leq +0.866$, or

$-0.866 \leq 0,105 \cdot (2k + 1) \leq +0.866$ or encore $-4.63 \leq k \leq +3.63$

The possible values of k for the amplitude minima are:

$k=-4, k=-3, k=-2, k=-1, k=0, k=1, k=2$, and $k=3$

⑥ Case of small angles

When $\sin(\theta)$ is small (for θ = about 25 deg. about on both sides of position M_0),

- we have for the peaks: $\theta \approx \sin(\theta) = k \cdot \lambda / a = 0.21 \cdot k$, with $k=-2, k=-1, k=0, k=1$ and $k=2$

$\theta_{-2} = -24.2^\circ; \theta_{-1} = -12.1^\circ; \theta_0 = 0^\circ; \theta_1 = 12.1^\circ; \theta_2 = 24.2^\circ$

- we have for the minima: $\theta \approx \sin(\theta) = (2k+1) \cdot \lambda / 2a = 0.105 \cdot (2k+1)$, with $k=-2, k=-1, k=0, k=1$ and $k=2$

$\theta_{-2} = -18.4^\circ; \theta_{-1} = -6^\circ; \theta_0 = 6^\circ; \theta_1 = 18.4^\circ; \theta_2 = 31.7^\circ$

❖ Comparison of experimental and theoretical results

The experimental measurements obtained agree fairly well with the theoretical results.

- As for the number of peaks and the number of minima, the extreme negative values could not be obtained but the other 8 ones could be observed.
- With small angles, the difference between the experimental positions and the theoretical values is slight:

Experimental peak (in degrees)	-18	-6	6	18	32
Theoretical peak (in degrees)	-18.4	-6	6	18.4	31.7
Experimental minimum (in degrees)	-24	-12	0	11.5	24
Theoretical minimum (in degrees)	-24.2	-12.1	0	12.1	24.2

Moduson complément

LAB WORK SHEETS

Student's Documents


1 Directivity of ultrasonic waves

Introduction: The purpose of the experiment is to demonstrate the directional effect of ultrasonic waves and to define the notion of « sound radiation ».

1.1 Step 1: Setup

① The system is set up as follows:

The transmitter is connected to a 15V power supply. It is set at its maximum transmission level and operated in *Continuous* mode. It is positioned at the centre O of the arc-shaped recess and it transmits ultrasonic waves towards the graduation zero.

The receiver is positioned in the circular recess closest to the transmitter (radius = 15 cm) in which it can be moved. It is connected by means of a shielded cable to the vertical input of the oscilloscope.

② Material:

List the equipment used

③ Setup schematic: draw the schematic of the system

1.2 Step 2: Protocol

❖ Directivity of ultrasonic waves


① Determination of the reception peak :

It is obtained in determining the receiver position on the arc of circle that yield the highest amplitude for the received signal.

Specify the settings used for the oscilloscope:

- vertical amplification :
- horizontal sweep :

In the frame below, draw a diagram of what you see on the oscilloscope screen.

What is the value of the angle that produces the maximum amplitude of the received signal?

② Sound transmission cones:


The ultrasonic waves are emitted by the transmitter in a zone to be determined.

In the table below, record the amplitude values a of the received signal for the values θ of the angle that identifies the ultrasonic receiver position.

(θ is negative on the left of the graduation zero and positive on the right).

a in V																	
θ in deg.	-40	-35	-30	-25	-20	-15	-10	-5	0	5	10	15	20	25	30	35	40

③ Plot the corresponding graph below:

④ **Approximation of the measurement:**

Determine the minimum angular variation that produces the highest amplitude on the oscilloscope screen.

$\Delta\theta \approx \dots\dots\dots$

❖ **Conclusion :**

Are ultrasonic waves directional?

What is the value of the angle of a transmitter transmission cone?

Draw the transmission cone on the diagram below.

Note: The experiment can be repeated by using the second recess for the displacement of the receiver.

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We call « sound radiation » the direction of transmission from the transmitter to the receiver corresponding to the peak amplitude of the signal received by the receiver. An arrow symbolises its direction.

❖ **Representation of sound radiation**

Represent this radiation on the diagram below.

E ☐ Transmitter

R ☐ Receiver


2 Reflection of ultrasonic waves

Introduction: The purpose of the experiment is to verify the law of reflection. The influence of the absorption phenomenon of waves ultrasonic will also be observed.

2.1 Step 1: Setup

❖ The system is set up as follows:

The reflective plane (mirror) with magnetic sole is positioned at the centre of the arc-shaped recesses. An ultrasonic transmitter and a receiver are positioned as shown on the diagram. The MODUSON main transmitter can be used directly but its base is too large and does not fit into the recess in the moulded PVC table. *It is preferable to use the complementary transmitter 2*, which has a suitable size, at the expense of a slightly bulkier setup. The transmitter 2 is positioned e.g. on the 30 deg. graduation on the left.


The transmitter 2 is connected to the main transmitter. The main transmitter is powered via a 15 V power supply and its switches are set on *Continuous* and *OFF*, respectively (its transmitting element is then cut off and only the transmitter 2 is powered).

The receiver is positioned in the circular recess closest to the transmitter (radius = 15 cm).

❖ **Material:**

List the equipment used.

Setup schematic: complete the following schematic (use the definitions given on the next page)

2.2 Step 2: Protocol

❖ Definitions

On the above schematic, the symbols have the following meanings:

- On is the perpendicular to the mirror
- i is the angle of incidence between the incident sound radiation EO and the perpendicular to the mirror On
- r is the angle of reflection between the reflected sound radiation OR and the perpendicular to the mirror On

❖ Measurements :

Different values are assigned to i by moving the transmitter in its guide and the receiver position to be used to obtain the direction of the radiation reflected on the mirror is determined.

How do you identify this direction?

What are the oscilloscope settings?

- for vertical amplification :
- or horizontal sweep :

Describe how a measurement is performed.


Note : The presence of other peaks (with lower amplitudes) is observed. This is due to complex phenomena called diffraction.

Record the results in the table below:

i in deg.	10	15	20	25	30	40	50	60
r in deg.								

❖ Plot the corresponding curve on the graph below:

Graphic representation of the function $i=f(r)$

Observations of the curve:

What is the nature of the curve?

❖ **Ultrasonic absorption:**

Complementary experiment.

The mirror is turned over so that the foam-coated side faces the ultrasonic beam from the transmitter.

For a few values of the angle i , try to obtain the corresponding values of r .

What do you observe concerning the signal received by the receiver?

2.3 Step 3: Exploitation of the results

❖ **Law of reflection**

Write the law of reflection that can be derived from the measurements performed.

❖ **Ultrasonic absorption**

What conclusion can you make?


Note: The experiment can be repeated by using the second recess for the displacement of the transmitter or (and) of the receiver.

3 Interference of ultrasonic waves

Introduction: The purpose of the experiment is to perform an experiment on ultrasonic interference and to compare the experimental measurements obtained with the theoretical results derived from the wave model of the phenomenon.

3.1 Step 1: Setup

❖ The system is set up as follows:

The transmitter 2 is positioned on a magnetic support inclined at 45 deg. This special position is necessary to achieve short distances between the two transmitters. Thus they contact each other in their smallest spacing (2 cm).

The transmitter 2 is connected to the main transmitter. The main transmitter is powered via a 15 V power supply and its switches are set on *Continuous* and *OFF*, respectively (its transmitting element is then cut off and only the transmitter 2 is powered).


The receiver is positioned in the circular recess closest to the transmitter (radius = 15 cm) in which it can be moved. It is positioned on the zero graduation. It is connected by means of a shielded cable to the vertical input of the oscilloscope.

❖ **Material:**

List the equipment used.

❖ **Complete the setup schematic**

(plot the 2 sound radiations that reach an arbitrary point M)

3.1 Step 2: Protocol

❖ Preliminary setting of the transmitters

- The main transmitter is assumed to be set at its maximum transmission level. It is connected to the MODUSON transmitter 2.
- The two transmitters connected to the same signal generator are synchronous sources. However, the following adjustments are required to give them an identical amplitude:
 1. Install the two transmitters in the magnetic recess on the moulded PVC table. Place them at a distance of e.g. 2 cm from each other (smallest distance). Place a receiver in one of the arc-shaped recesses opposite the zero-degree angle value (then the receiver is located on the median of the segment joining the two receivers and it is equidistant between them) and connect it to the vertical input of an oscilloscope.
 2. Disconnect the transmitter 2 and observe the signal transmitted by the main transmitter and sensed by the receiver. Record the amplitude value. Connect the transmitter 2 and switch the main transmitter OFF (in this position, the transmitter member is off but the generator is not and continues delivering the same signal on the BNC connector, hence to the transmitter 2), then only the transmitter 2 is active. Using the transmitter potentiometer, trim the amplitude of the received signal so that the amplitude value is identical to the previous reading.
 3. Switch the main transmitter ON: then the two transmitters operate simultaneously and they have the same amplitude. This produces a peak that is sensed by the receiver (a slight position adjustment by 1–2 degrees may be necessary to obtain this peak).

❖ Measurements:

What are the oscilloscope settings?

- for vertical amplification :
- for horizontal sweep :

Select the distance a between the 2 transmitters and the radius D of the receiver displacement arc.

- $D=$
- $a=$

Move the receiver and observe the amplitude peaks and minima of the signal sensed by the receiver and observe on the screen of the oscilloscope.

What explanation can you give about this observation?

How to proceed to perform correct measurements?

Note: The minima are more easily identified than the peaks. However, the minima are not zero and they are less and less marked. What explanation can you give about this finding?

❖ **Record the measurement results in the table below:**

The angles are given in algebraic notation relative to the graduation zero.

θ minimum									
minimum difference									
θ maximum									
minimum difference									

What is the measured variation for the angle θ from one maximum to the adjacent maximum (and from one minimum to the adjacent minimum)? In what angular domain is this verified?

What is the number of the observed peaks and minima:

- peaks
- minima

What can be said about the difference between the peaks and the minima and about their alternation?

3.2 Step 3: Exploitation of the results

❖ **Theoretical study of the interference phenomenon**

① **The superposition model of the waves from the two synchronous sources is used.**

The wavelength is λ .

At an arbitrary point M, the waves from the sources S_1 and S_2 exhibit a course difference $\delta = d_2 - d_1$.

What do you observe if the two waves arrive *in phase* at point M?

The relationship that renders this condition is: $\delta =$

What do you observe if the two waves arrive *in opposition of phase* at point M?

The relationship that renders this condition is: $\delta =$

② Calculation of the course difference

In the above setup schematic, the course difference is expressed as $\delta = d_2 - d_1$.

d_2 and d_1 are the hypotenuses of two rectangular triangles (MhS_1 and MhS_2).

Express the squares of these hypotenuses:

$$d_2^2 =$$

$$d_1^2 =$$

Write the difference:

$$d_2^2 - d_1^2 = (d_2 - d_1)(d_2 + d_1) =$$

verify that:


$$\delta = \frac{2aD \sin(\theta)}{\sqrt{D^2 + \frac{a^2}{4} + aD \sin(\theta)} + \sqrt{D^2 + \frac{a^2}{4} - aD \sin(\theta)}}$$

③ Approximations

The computerised study of the denominator $\Phi(\theta)$ of δ by means of ACTILAB yields its variations in the experimental conditions.

Plot the representative curve of $\Phi(\theta)$ versus θ :

(θ varies at most from -60° to $+60^\circ$ ($-\pi/3$ to $+\pi/3$))

Calculate the relative difference $(\Phi(\theta) - 2D)/2D$ for the various positions:

$$\theta = 0 \quad :$$

$$\theta = +\pi/6 \text{ or } -\pi/6 \quad :$$

$$\theta = +\pi/3 \text{ or } -\pi/3 \quad :$$

What conclusion can you make?

What is the approximate expression of the course difference: $\delta \approx$

④ Calculation of the number of peaks

What is the expression of δ for the peaks? (with $\lambda = 0.84$ cm).

Since $\frac{\lambda}{a}$ is constant, what can be said about $\sin(\theta)$?

Calculate the numerical values of k for $-\pi/3 \leq \theta \leq +\pi/3$.

The possible values of k for the amplitude peaks are:

$$k =$$

⑤ Number of minima

What is the expression of δ for the minima? (with $\lambda = 0.84 \text{ cm}$).

Since $\frac{\lambda}{2a}$ is constant, what can be said about $\sin(\theta)$?

Calculate the numerical values of k for $-\pi/3 \leq \theta \leq +\pi/3$.

The possible values of k for the amplitude minima are:

$$k =$$

⑥ Case of small angles

When $\sin(\theta)$ is small (for $\theta =$ about 25 deg. about on both sides of position M_0), the angle and its sine may be used indifferently: $\theta \approx \sin(\theta)$

- calculate the angles corresponding to the observable peaks (note the index of θ):

$$\theta = \quad \theta = \quad \theta = \quad \dots$$

- calculate the angles corresponding to the observable minima (note the index of θ):

$$\theta = \quad \theta = \quad \theta = \quad \dots$$

❖ Comparison of the experimental and theoretical results

Do the experimental measurements obtained agree with the theoretical results?

Do the number of peaks and the number of minima correspond to those observed?

Record and compare in the table below the experimental and theoretical values corresponding to of the small angles.

Conclusion:

Experimental peak (in degrees)					
Theoretical peak (in degrees)					
Experimental minimum (in degrees)					
Theoretical minimum (in degrees)					

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux