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## *Electron beams*

Réf :  
232 005

Français – p 1

English – p 6

Version : 9101

**Déflectron**

***Deflectron***

## 1. Principes et description

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**0 825 563 563\***

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## 1. Principles and description

The DEFLECTRON allows for the study of the behaviour of electrons in an electromagnetic field, and for the realisation of measurements on the deflection of these electrons. It contains:

- 1 stand with:
  - The pedestal that is made of alloy cast under pressure.
  - 1 loading arm used as a handle for transportation.
  - 1 fork support of the tube made of high mechanical and electrical resistance plastic. 2 clamps efficiently maintain the two axes of the tube.
- 1 set of Helmholtz coils: fixed on a rod made of stainless steel finished with a plastic protective tube, they allow for quick fixation to the stand. To remove the coils, lift the plastic rings and remove the rods and coils from their compartments.
- 1 electron tube:

It is a spherical tube that is under vacuum. It contains:

- The electron gun with its filament connected to the terminals Ø 4 mm, and an anode punctured with a slot that accelerates and focuses the electrons into a rectangular probe. The anode comes out laterally on a male banana pin Ø 4 mm.
- A calibrated screen covered with a fluorescent substance is placed sideways with respect to the electron beam in order to materialise its trajectory.
- Two electrical deflecting electrodes are connected to the outer part of the tube through 2 male plugs Ø 4 mm.

To remove the tube, slide the clamps of the arm backwards and push the tube forward while pressing with the palm of your hand on the base of the tube.

## 2. Technical characteristics

### 2.1 The tube

- Filament power supply: alternating 6 V
- Acceleration voltage: direct from 2 to 5000 V
- Voltage between the deflecting electrodes: direct from 0 to 5000 V
- Anode I max.: 1 mA
- Fluorescent screen, 90 mm length inclined at 15°
- Deflecting electrodes, 90 mm length, 52 mm distance
- Dimensions: Ø 136 mm – L = 260 mm

#### 2.1.1 The pedestal

- Collar made of alloy that was cast under pressure. Dim: 240 x 180 x 50 mm
- Weight: 3,2 kg
- A loading arm used as a handle for the transportation of the deflectron
- A fork support for the tube, made of insulating plastic

#### 2.1.2 The Helmholtz coils

Every coil has 320 wire turns made of enamelled copper

- average Ø: 136 mm.
- Rod made of stainless steel, finished with a protective plastic tube for fixation on the pedestal.

- Sidewise connection, female plugs Ø 4 mm.
- I max: 1 ampere in continuous rating (1,5 A for 10 minutes).

## 2.2 Accessories

To carry out experiments with the DEFLECTRON, the following sources and measuring devices are required:

### 2.2.1 For the Helmholtz coils

D.C supply from 6 to 12 V – 1 A (can go up to 18 Volts / 1,5 A for a few minutes).

### 2.2.2 For the tube

- 6 V / 2 A power supply for the filament,
- 5000 V power supply for the acceleration of the electrons and the deflection between the electrodes.

For certain experiments, it is worth having 2 high-voltage power supplies: one for the electron gun and one for the study of the deflection of the electrodes.

### 2.2.3 Spare the for DEFLECTRON Réf: 233 008

## 3. Experiments

### 3.1 Electrostatic deflection

#### 3.1.1 Theoretical calculation

For an electron of mass  $m$ , charge  $e$ , and speed  $v_0$ , the relation  $F = ma$  is applied. Hence:

$$x = v_0 t$$

$$y = -\frac{1}{2} \frac{eEt^2}{m}$$

The trajectory of the beam is a parabola of equation:

$$y = -\frac{1}{2} \frac{e}{m} \frac{E}{v_0^2} x^2$$

Otherwise, it can be written that:

$$\frac{1}{2} m v_0^2 = e v a$$

$$\text{with } E = k \frac{va}{d}$$


$K$  is a factor due to the geometry of the deflecting electrodes that present strong edge effects.

$$y = \frac{k}{4d} x^2$$

#### 3.1.2 Experiment

Connect one plate of the capacitor to the cathode and the other plate to the anode. Invert the + and the – (see figure 2). Observe the trajectories

Figure 2 :


**Comments:**

- To preserve the tube, do not permanently leave in operation, or alternatively decrease the filament voltage to  $\approx 5$  V.
- Obscure the tube if bothered by ambient light.
- To enhance the precision on the measurements, select the parameters in order to obtain a beam going through a node on the squared screen.
- If the same voltage is applied between the electrodes and the gun, the deflection obtained will not change because while the electrons are accelerated, the deflecting electrical field increases.
- If you have 2 power supplies, the parameters can be varied independently (see figure 3).

## 3.2 Electromagnetic deflection

### 3.2.1 Theoretical calculation

Field generated by two flat coils in Helmholtz position.

It can be written that:


$$B = \frac{\mu_0 R^2 N I}{(R^2 + x^2)^{3/2}} \text{ with } \begin{matrix} R = 68\text{mm} \\ N = 320 \end{matrix}$$

$$\Rightarrow B = K I$$

Power up the two in series coils by means of a 0-12 Volts power supply. The maximum continuous current in the coils is 1 A which corresponds to a 12 V voltage. If the power supply is  $> 12$  V a 1,5 A current can be provided for  $t \leq 10$  minutes.

Example of assembly:

Figure 6


Trajectory of an electron in a field B normal to its speed:

$$\frac{Bev}{m} = \frac{v^2}{R}$$

R = radius of curvature of the trajectory  
V = speed of the electron

Figure 7

The trajectory of the electron in the system (figure 7) can be written:

$$R = \frac{x^2 + y^2}{2y} \quad \text{equation of a circle}$$

### 3.2.2 Experiment

Power up the tube. Observe the beam.

- . Use V(anode) = 4000 V.
  - . Set several values for I(coils). Change the direction of the current.
  - . With I(coils) being fixed, increase V(anode) from 3000 to 5000 Volts.
- Measure R for each experiment.

### 3.3 Determination of the e/m ratio

*J.J THOMSON's method*

Simultaneously apply an electrical field E between the electrodes and a magnetic induction B on the screen in such a way that the effects are opposite or equal. In this case it can be written that:

$$eE = ev \wedge B$$

$$\text{hence } v = \frac{E}{B}$$

Suppress E (put the two electrodes at the voltage of the anode) and measure R corresponding to the circular trajectory due to the magnetic deflection.

$$R = \frac{mv}{eb}$$

$$\frac{e}{m} = \frac{E}{B^2 R} = \frac{va}{d} * \frac{1}{K^2 I_b^2 R}$$

$d = 52 \times 10^{-3} \text{ m}$  being the distance between the electrodes.

#### Experiments and results:

Perform the assembly of the "Electrostatic deflection" chapter (see figure 2). The values thus found for  $e/m$  are too high. This is due to the edge effects of the deflecting capacitor: the field E is not uniform.

#### Alternative method:

Measure the radius R, and plot of the electron beam deflected by a field B ( $E = 0$ ).  
Speed of the electrons:

$$e va = \frac{1}{2} mv^2$$

$$R = \frac{mv}{eB}$$

$$\boxed{\frac{e}{m} = \frac{2 va}{B^2 R^2}}$$

**Experiment:** measure the radius R corresponding to B(anode) and V(anode).

**Comment:** With the Millikan's experiment, the charge e of the electron can be determined.  
With the DEFLECTRON, the ratio  $e/m$  will be determined. Therefore, m can be deduced.

### 3.4 Experiment of the electronic « mirror »

When there is no more equilibrium between the actions of E on the electrons, from the following equation

$$F = e(E + v \wedge B) = m$$

the trajectory is deduced:

$$y = -\frac{eEt}{mw} + \left( \frac{eE}{m} + w v_0 \right) \frac{\sin wt}{w^2}$$

$$z = -\frac{1}{w^2} \left( \frac{eE}{m} + w v_0 \right) (1 - \cos wt)$$

$$\text{with } w = \frac{eB}{m}$$

$$\text{if } v_{0x} \text{ is negligible compared to } \frac{eE}{mw} = \frac{E}{B}$$

The parameters of a cycloid equation are found:

$$x = -\frac{eE}{mw^2} (wt - \sin wt)$$

$$y = -\frac{eE}{mw^2} (1 - \cos wt)$$

$$y \text{ maximum: } y_M = \frac{2 m E}{e B^2}$$

### Interpretation:

(1)  $v \rightarrow 0 \Rightarrow e\mathbf{v} \wedge \mathbf{B} = 0$ . Action of B is null

$F = -eE$

The electron has an accelerated uniform motion.

(2)  $v \rightarrow 0 \Rightarrow F = e\mathbf{v} \wedge \mathbf{B}$  tends to curve the trajectory.

(3) If B is sufficient the electron will not reach the + electrode. v is maximum.

(4) E tends to slow the electron down and  $e\mathbf{v} \wedge \mathbf{B}$  decreases.

(5)  $v \rightarrow 0$  and the electron is rejected by the electrode.

The trajectory is a cycloid whose maximum deflection is:

$$|y_M| = 2 \frac{me}{eB^2}$$

$$\Rightarrow \frac{e}{m} = \frac{2E}{b^2 y_M}$$

### Experiment:


- The speed of the electrons is  $v_{ox}$ . Only the trajectory after (1) must be considered.  
 $V(\text{anode}) = V(p1p2) < 3500 \text{ V}$  will be used.

- Increase I(coils) until suppression of the deflection ( $I_b = 2 \text{ A}$  max during 1 minute can be reached).

It is not necessary to obtain the second turn in order to carry out the measurement.

Measure  $y_M$ . Deduce  $e/m$ .

Figure 8.

## 4. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.  
For any repairs, adjustments or spare parts please contact:

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## Direct connection for technical support

A team of experts  
at your disposal  
from Monday to Friday  
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request  
immediatly to provide you  
with the right answers regarding  
your activity field : Biology, Physics,  
Chemistry, Technology.

### Free service\*

**+33 2 32 29 40 50\*\***

\* Call cost not included.

\*\* Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux