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Material properties

Réf :
242 013

Français – p 1

English – p 4

Version : 9101

Sonde de pression

Pressure probe

1. Principe et description

1.1 Principe

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1. Principle and description

1.1 Principle

The pressure probe, complete with a water or electronic pressure gauge, is used to study the distribution of pressures existing inside a liquid.

1.2 Description

The pressure probe includes:

- A cylindrical capsule of Ø 50 mm for exploration of the liquid media. The pressure is exerted on an extremely flexible rubber membrane covering the capsule and held stretched by a removable flexible bushing
- A replacement membrane
- A system with two coaxial tubes made of plastic:
 - The outer tube in the form of a J. It is screen-printed with marks every 50 mm
 - The inner tube connecting the capsule and the tube opening onto the pressure gauge can be rotated with a knurled knob: the movement is transmitted to the capsule via a watertight flexible connection.

2. Commissioning - fixing

Install the pressure probe vertically by locking the horizontal metal rod with a tightening nut on a base with rod (not provided). Connect the pressure probe to a water or electronic pressure gauge.

3. Use

3.1 Equipment required

- A container with Ø 130 mm (Part number 713262)
- A base, a retaining nut and a rod
- E.g.: Modumontage bracket (701292) + Polynux (703452) + Initio pressure meter (251181)

To make best use of the pressure probe, it is recommended that a container with a minimum depth of 200 mm be used, so that the entire working depth of the J tube can be used.

3.2 Experimentation

Sliding easily in the nut that holds it, the tube can be raised or lowered; it can pivot on itself by exploring the entire surface of water by the capsule.

The J tube being stopped at a level, by rotating the top knurled bushing; the capsule can be pivoted around a diametral axis, which shows that in these conditions, the average pressure that it supports is independent of its position.

4. Maintenance

The pressure probe requires almost no maintenance, as none of its constituent parts are sensitive to water. It is however recommended that you wipe it after every use.

When not in use, gently powder the rubber membrane and close the open port of the pressure gauge to avoid evaporation of the liquid.

The membrane can be changed after a certain period (puncture, ageing). In view of this, a spare membrane is provided. To install it:

Remove the plastic bushing from around the capsule

Replace the bushing: the membrane is automatically stretched.

5. Spare parts

- Spare membranes with attachment device

Part no. 242012(for 10 installations)

6. After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux