# **Mécanique**

Dynamique

# ***Mechanics***

*Dynamics*

**Ref :  
243 010**


**Tube de Newton**  
***Newton's tube***

Français – p 1

English – p 3

Version : 9001

## 1 Description

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
## 1 Description

Plastic tube equipped with 2 nozzles and a vacuum valve, containing 2 falling elements: a steel ball and a paper washer.

2 aluminium plates placed inside the nozzles help limit the electrostatic phenomena.

### ***Additional equipment***

Membrane vacuum pump  
Vacuum pipe inner  $\Phi$  8 mm

**Part no. 701 290**  
**Part no. 253 044**

## 2 Set-up

To make assembly easier, it is advisable to immerse the nozzle in water heated to about 50°C in order to soften the plastic.

1. Place the disk and the ball at the bottom of the tube. Hold the tube with one hand and the nozzle with the other, then forcefully fix the nozzle till the second groove.
2. Press hard on the nozzle with both thumbs to complete the assembly. A small hammer (not supplied) may also be used to force the nozzle in.

### 3 Technical characteristics

<b>Nozzle:</b>	for pipe inner $\Phi$ : 8 mm
<b>Dimensions:</b>	$\Phi$ 50 x 1060 mm
<b>Mass:</b>	540 g
<b>Guarantee:</b>	2 years

### 4 Experiments

#### 4.1 In air

Position the tube in the vertical position with the falling elements centred at the bottom of the tube, on the side opposite the vacuum valve.

Rotate the tube by 180°.

We observe that the 2 elements don't reach the bottom of the tube at the same time.

#### 4.2 In vacuum

Connect the tube to a vacuum tube with a vacuum pipe of 8 mm inner  $\Phi$  (not supplied).

Open the valve, vacuumize, then close the valve again and disconnect the vacuum pump pipe.

Position the tube in the vertical position with the falling elements centred at the bottom of the tube, on the side opposite the vacuum valve.

Rotate the tube by 180°.

We observe that the 2 elements reach the bottom of the tube at the same time.

**Note:** In case of a friction problem, bring the falling elements in contact with the vacuum valve then tighten the metal part of the valve by hand to remove electrostatic charges.

### 5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
Rue Jacques Monod  
BP 1900  
27 019 EVREUX CEDEX FRANCE  
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your  
disposal from Monday  
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries