Mécanique

Forces Pressantes

Mechanics

Pressure Forces

Réf :
243 035

Maquette Pression Surface

Surface Pressure Model

Français – p 1

English – p 6

Version : 8108

1. Description

1.1 Généralités

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1. Description

1.1 Generalities


This model highlights experimentally the relation between the three following physical quantities:

- Pressure
- Force
- Surface

The student is able to evaluate and measure the influence of each of the parameters, and validate the relation $P=F/S$.

The measurement can be performed

- Qualitatively: we then compare the pistons movement through syringes graduations thus representing comparison visual cues
- Quantitatively: we then measure the different pressures using a Computer Assisted Experimentation sensor or an Initio2 pressure meter for example.

1.2 Practically

The model is made of 3 syringes of different diameters allowing to act on the surface value: in fact, the piston diameter is the reference value for the surface.

This reference diameter is indicated on the model.

The syringes are fitted with plastic pipes for a direct connection to a pressure measurement sensor (Computer Assisted Experimentation, Initio and Initio2 pressure meter sensor)

Using test weights allows acting on the force value: in fact, each of the syringes is fitted with a mounting platform containing the weights.

An additional accessory allows connecting the two biggest syringes together in order to address the notion of hydraulic press.

The model has a storing notch for setting this accessory and thus avoid its loss.

2. Qualitative Measurements (without Measurement Sensor)

2.1 Required Equipment

In the case of a so-called « qualitative » measurement, only the weights –will be needed.

Series of brass weights 500 g
Series of brass weights 1000g

2.2 Protocol

2.2.1 At Constant Force

Make sure that the base of the syringes pistons is aligned with the black mark of the model and that the valves are in locked position.

Place a weight in the weights mount on one of the syringes.

Use the same weight on a second syringe.

We note a difference in the distance lowered by the piston for each of the 2 syringes.

Warning:

Each syringe has a suitable weight range; below, the piston does not react, and above, it is tightened.

Large model syringe	10 mL:	from 200g to 1000g
Medium model syringe	3 mL:	from 100g to 800g
Small model syringe	1 mL:	from 50 g to 100g

The distance travelled by the piston (piston « fall ») is followed by a proportional increase in pressure within the syringe.

The piston fall is thus proportional to the pressure.

2.2.2 At Constant Surface

Make sure that the base of the syringes pistons is aligned with the black mark of the model and that the valves are in locked position.

Place a weight in a weight mount of a syringe, starting with the lightest weight.

Change the weight using a higher weight and so on.

We note an increase in the fall distance of the piston in accordance with the increase of the positioned weight.

The distance travelled by the piston is followed by an increase in pressure.

The piston fall is proportional to the pressure.

2.2.3 Results - Observations

We thus note from $P=F/S$ that

- For F constant: the pressure (fall distance) increases when the surface decreases
- For S constant: the pressure (fall distance) increases with the force

3. Quantitative Measures (with Measurement Sensor)

3.1 Necessary Equipment

Computer Assisted Experimentation Sensor	Primo / VTT Range
	Visio / Visio+ Range
	Tooxy / Foxy Range
Others	Initio 2
	Manual or digital pressure meter
Weights	Series of brass weights 500 g
	Series of brass weights 1000g

3.2 At Constant Force

Make sure that the base of the syringes pistons is aligned with the black mark of the model and that the valves are in locked position.

Note the pressure displayed on the measuring device.

Place a weight in the weights mount on one of the syringes.

Note the pressure value.

Use the same weight on a second syringe.

Note the pressure value.


We note a difference in the distance lowered by the piston for each of the 2 syringes.

Warning:

Each syringe has a suitable weight range; below, the piston does not react, and above, it is tightened.

Large model syringe	10 mL:	from 200g to 1000g
Medium model syringe	3 mL:	from 100g to 800g
Small model syringe	1 mL:	from 50 g to 100g

The distance travelled by the piston (piston « fall ») is followed by a proportional increase in pressure within the syringe.

3.3 At Constant Surface

Make sure that the base of the syringes pistons is aligned with the black mark of the model and that the valves are in locked position.

Note the pressure value.

Place a weight in the weight mount of a syringe, starting with the lightest weight.

Note the pressure value.

Change the weight using a higher weight and so on by noting the pressure value for each weight.

We note an increase in the fall distance of the piston in accordance with the increase of the positioned weight.

The piston fall is proportional to the pressure.

The following measurement values are advised:

Large model syringe	10 mL:	from 200g to 1000g	In graduations of 200g
Medium model syringe	3 mL:	from 100g to 800g	In graduations of 100g
Small model syringe	1 mL:	from 50 g to 100g	In graduations of 10g


4. Results - Observations

We thus note from $P=F/S$ that

- For F constant: the pressure (fall distance) increases when the surface decreases
- For S constant: the pressure (fall distance) increases with the force

4.1 Results Example

Measurements taken with an Initio2 pressure meter on a constant Surface:

We thus note an equivalence between P_{measured} and $F/S_{\text{calculated}}$ then allowing to verify

Pressure=Force/Surface

5. After-Sale Services

The warranty is of 2 years.

For all settings, contact out **Technical Support** on **0 825 563 563**.

For all repairs, spare parts or return to our workshops, kindly contact:

JEULIN – S.A.V.
468 Jacques Monod Street
CS 21900
27019 EVREUX CEDEX France

0 825 563 563*

* 0.15 € TTC/min. from a landline

Assistance technique
en direct

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Direct connection
for technical support