Mécanique

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Mechanics

Fluid Dynamics

Réf : 243 038

Maquette Bernoulli

Bernoulli's Model

Français – p 1

English – p 5

Version : 1206

1. Description

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0 825 563 563*

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1. Description

This device allows highlighting Bernoulli's principle, more precisely its specific case : the Venturi effect.

Composed of a blower, 3 different diameter tubes and a support, it will consist in measuring the pressure difference between 2 different diameter tubes to deduce the flow of the blower, for example; the flow being considered constant.

Each of the 3 tubes has an output allowing to measure using a pressure meter, and have a different color in order to facilitate their differentiation (white, orange and grey). Thus, 6 different configurations are possible.

The blower is reversible: it can be used as vacuum or blower.

It is advised to use the differential pressure meter.

2. Experimental Protocol

Each of the tubes has a side with a joint which allows to ensure maximum sealing between each tube.

Fix both tubes on the blower while securing the support on the second tube.

Connect the pipes of a pressure sensor to each measuring tip.

Supply the blower and then activate its switch.

Wait until the measurement has stabilized and then note the value of the measured pressure difference.

It will be necessary to repeat the process several times in order to use an average of values in the calculations that will follow.

Make sure to verify that no object nor cloth is blocking the free end of the blower in order not to result in wrong values.

Change the tubes and repeat the operation.

The flow value of the blower can be deduced with Bernoulli's principle.

Our experimental measurements have determined a flow value of 13 L/s.

Therefore, based on this value, the differences in diameter of each tube can be determined and compared with the manufacturer's values.

3. Results Example

For a frictionless flowing fluid, Bernoulli's theorem states that:

$$P_A + \frac{1}{2} \rho \times v_A^2 = P_B + \frac{1}{2} \rho \times v_B^2 = P_C + \frac{1}{2} \rho \times v_C^2$$

Knowing that:

$$Q_v = v_A S_A = v_B S_B \quad \text{and} \quad v_B > v_A \quad \text{hence} \quad P_A > P_B$$

Thus, the following formula can be used to deduce the flow of the blower:

$$P_A - P_B = \Delta P = \frac{1}{2} \rho \times \left(\frac{1}{S_B^2} - \frac{1}{S_A^2} \right) \times Q_v^2$$

$$\Delta P = \frac{1}{2} \rho \times \left(\frac{1}{S_B^2} - \frac{1}{S_A^2} \right) \times Q_v^2$$

$$Q_v = \sqrt{\frac{2 * \Delta P}{\frac{1}{S_B^2} - \frac{1}{S_A^2} \times \rho}}$$

With:

Q_v = volume flow rate of the fluid in m³/s

ΔP = differential pressure in Pa

S = pipe section in m²

ρ = air density in kg/m³

	Average values of ΔP in kPa	Calculation of Q_v in m ³ /s	% error with respect to the theoretical value
Between White and Orange	0.68	0.014	16 %
Between White and Grey	2.08	0.011	-8 %
Between Orange and Grey	1.87	0.012	0 %

We note that the values measured allow finding the flow values consistent with the previous experimental values.

Discussing the sources of error may lead to uncertainty of measurement calculations, particularly in connection with the measuring device.

4. Technical Specifications

White tube diameter	3.106 cm
Orange tube diameter	2.149 cm
Grey tube diameter	1.516 cm
Flow of the blower (experimental measurements)	12 L/s

5. After-Sales Service

The guarantee period is 2 years.

For any repairs, contact the **Technical Support** on **0 825 563 563**.

The equipment should be returned in our workshops and for any repairs or spare parts, kindly contact:

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux