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Material properties

Réf :
251 014

Français – p 1

English – p 3

Version : 9101

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Simplified calorimeter

1. Description

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1. Description

The apparatus consists of :

- An inner vessel: internal \varnothing 100 mm - height 130 mm - thickness about 7/10 mm, in bright, oxidised aluminium, with lip on the upper edge.
- An outer vessel: internal \varnothing 120 mm - height 160 mm - thickness 7/10 mm, in bright, oxidised aluminium, with seam on the upper edge.
- A removable coil in plastic material designed to cap the outer vessel and to locate the inner vessel.

2. Characteristics of the inner vessel

- Area of cross section: 78.5 cm² (hence there are 78.5 cm³ if liquid per cm depth).
- Total volume: 1,020 cm³
- Weight of vessel empty: about 80 g

Theoretical "water value": will be obtained in gram calories by weighing the inner vessel and multiplying the weight in grams by the specific heat of aluminium equal to 0.217. It will thus be found that the value is about 18 calories.

You can also measure it experimentally by the mixing method, which makes an excellent experiment to familiarise students with the apparatus.

3. Heat loss

We have established the heat loss curves experimentally for the calorimeter with 200 and 350 g of water, ambient temperature being 18°C.

The following table gives the mean temperature drop per minute for temperatures varying from 36 to 19 degrees.

Temperature	Average drop in temperature in 1 minute	
	200 g of water	350 g of water
36	0.23	0.19
34	0.21	0.18
32	0.19	0.13
30	0.16	0.11
28	0.14	0.10
26	0.12	0.07
24	0.07	0.06
22	0.05	0.04
21	0.04	0.03
20	0.03	0.03
19	0.01	0.01

You will notice that at 10 degrees above ambient, the heat loss is such that the drop in temperature is about 1/10 degree C, per minute, but that, on the other hand, it becomes very low when it is maintained within 3 or 4 degrees from ambient.

4. Accessories

- 251 015 : Cover in transparent plastic material with aluminium agitator, hole for thermometer and shielded opening for inserting the mass. The cover again limits heat loss during the experiment and avoids splashes of liquid.
- 253 019 : Thermometer in 1/10ths of a degree -10 + 35° C ;
Ø 5 mm ; length 35 cm.

5. After-Sales Service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux