Méthode ESAO®

ESAO® method

**Réf :
251 016**

Français – p1
English – p3

Version : 7008

Sismomètre Initio SA

Initio SA seismometer

1. Introduction

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
1. Introduction

The seismometer allows displaying and recording a wave generated by a mechanical shock.

- You can **display the signals in real time** in the form of an LED which lights up increasingly brightly, depending on the intensity of the signal.
- You can **visualise and save the signal** in the form of an LED which lights up more brightly, depending on intensity of the signal.
- You can **transfer the raw signal of the wave to a computer assisted experiment interface** via an analogue ± 5 V output.

2. Description

The Initio SA seismometer is composed of a piezometric sensor linked to an Initio unit powered by a 9 V battery (supplied), or with a 12V 0.6A power supply unit ref.281612 (not supplied).

- (2) The seismometer provides 2 modes:
- Upper position: real time display of signal.
 - Lower position: save signal.

3. Additional equipment

For a simple demonstration of the impact of distance on a wave, it is possible to use a seismic scale model (Ref. 507 036).

4. Start-up and first use

- Insert a battery in the housing of the casing located in the upper section of the seismometer,
- Start by using the on/off button located under the seismometer casing.
- If you want to use 2 seismometers for a quantitative approach to the phenomenon:
 - place the 2 seismometers side by side,
 - cause a shock which is equidistant from the 2 sensors,
 - adjust the response using the adjustment button on the casing.
- Select the use mode you want (recording of the wave signal or real time display of the wave signal).
- Place the red sensor flat on the surface where the shock is to be made.
- Cause a shock.

Since the sensor is extremely sensitive, a small shock (pencil falling from a height of 10cm or breaking a piece of plastic on the seismic scale model for example) is sufficient to obtain correct display of the signal.

Example of experimental set-up with the seismic scale model:

5. Precautions for use

- Remember to switch off the seismometer after use checking the ON diode has gone out.
- Do not tap too hard, this could cause rebounding of the sensor.
- If possible, attach the sensor to the support which must itself be stable to avoid rebounds.

6. After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux