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Mesure de point de fusion

Chemistry measurement

Determining the melting point

Réf :
251 027

Français – p 1

English – p 8

Version : 2107

Banc chauffant Kofler
Kofler heating bench

Sommaire

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0 825 563 563*

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Summary

1	Introduction and set up	8
1.1	Kofler heating bench features	8
1.2	Diagram of the apparatus	9
1.3	Electrical connection	9
1.4	Temperature increase	9
1.5	Reading the temperature	9
1.6	Ambient temperature compensation	10
1.7	Adjusting the reading device	10
1.8	Accuracy	11
1.9	Observation of the melting process	11
1.10	Control substances and test substances	11
1.11	Substances susceptible to decomposition	12
1.12	Volatile substances	12
1.13	Eutectic temperatures	12
1.14	Hydrates	12
1.15	Melting point of mixtures	13
1.16	Softening temperature of synthetic powders and resins	13
2	Care of the heating bench	13
3	Bibliography	13
4	Use precautions	14
4.1	Breakdown or malfunction	14
4.2	Heat and electrical hazards	14
5	After-Sales Service	14

1 Introduction and set up

1.1 Kofler heating bench features

The Kofler heating bench is a tool designed to determine the melting point for an object by a quick and easy method. The principle on which it operates and the ease in handling make this apparatus suitable for general purpose tasks.

With the heating bench, you can not only determine the melting point but also observe temperature-dependent properties. In effect, the heating bench can be used for determining the degree of purity of various substances, the presence and the behaviour of crystallization water, as well as for observing decomposition and volatilisation, and determining eutectic temperatures and softening points of synthetic resins and powders.

The greatest advantage in using the heating bench rather than other methods is that it is exceptionally time saving.

It only takes 1 to 2 minutes to determine a melting point, including the prior adjustment of the reading device. When you have a substance with an unknown melting point, spread a few small crystals in the midzone of the heating surface (i.e., in the 100 to 140 °C temperature range).

You will then see right away if you need to push the substance to the left or to the right of the bench (i.e., towards higher or lower temperatures).

The melting temperature must be read at the precise location on the heating surface where the separation line between the liquid and solid particles establishes itself. This separation line is generally distinct if you do not use an excessive amount of the substance.

1.2 Diagram of the apparatus


Fig. Kofler heating bench with electronic control.

1 – On-off switch	4 – Slide
2 – Heating surface	5 – Pointer
3 – Graduated scale from +50 to 260 °C	6 – Guide

1.3 Electrical connection

The heating bench can be plugged directly into the network and started with the on-off switch (1) located on the right side.

The apparatus runs on a 230 V- 50Hz AC current.

The built-in electronics compensate for household current fluctuations.

1.4 Temperature increase

The heating bench requires ca. 40 minutes to attain its operating temperature. If frequent use of the bench is anticipated, it is recommended to leave it plugged in all day. The warm up time is thus avoided and the bench is always ready for immediate use.

The maximum power absorbed by the apparatus is 150 W/h.

1.5 Reading the temperature

A temperature scale graduated in 2 °C increments, ranging from +50 to 260 °C (3), and a reading device are provided for reading the temperature.

The latter consists of a slide (4), which is moved with the pointer (5) along the scale, and a guide (6).

The tip of the pointer is located above the surface of the heating bench (2). The guide mounted on the slide can be moved perpendicular to the scale.

1.6 Ambient temperature compensation

The temperature reading device is designed to take the influences of variations in ambient temperature into account.

In effect, the temperature fluctuations in the environment in which the bench is used affect the equilibrium temperature of the heating surface.

Consequently, it is necessary to adjust the reading device relative to the ambient temperature prior to each measurement.

1.7 Adjusting the reading device

The reading device is best adjusted in the following manner:

With the spatula (included with the device), place a small quantity of one of the control and test substances (see the following paragraph) on the heating surface (2);

for example, benzoine (melting temperature of 132-135 °C), on the heating surface near the 135 °C mark.

After at least 10 seconds, adjust the slide (4) along the scale so that the pointer (5) is located exactly on the border between the melted and non-melted fractions of the substance deposited on the bench.


In order to avoid parallax errors, be sure to take the reading in the exact position of coincidence between the pointer and its reflection on the heating surface. You can observe the liquid-solid separation with a magnifying glass in order to obtain a better view of the separation line.

Once the pointer is properly positioned, move the sliding guide (6) vertically until its indicator needle (7) is on the 135 °C value (between the 134 °C and 136 °C marks) of the temperature scale (3). The heating bench is now calibrated.

For any subsequent determination of a melting point, if you move the slide along the scale so that the pointer (5) is located on the separation line between the solid and liquid fractions of the substance, the indicator needle (7) of the guide will then indicate the corresponding temperature on the scale (3).

As long as the ambient conditions (especially the temperature) remain the same as those in which the reading device was set, the calibration will be valid.

Any change in these conditions requires a new calibration.

1.8 Accuracy

The accuracy of the measurements taken with the heating bench depends upon the care taken in the adjustment of the reading device (paragraph 1.6). If you need to obtain very accurate values, you must calibrate the bench immediately prior to examining the substance to be studied, using that substance among the control substances provided of which melting point is closest to that of the substance to be examined. Under these conditions, you can count on a precision of around 1 °C for the temperature measurement.

If however, the heating bench is only calibrated in the medium temperature range (with benzoine, for which the melting temperature is 132-135 °C), you cannot count on any greater precision than +2.5 °C in the remotest zones (ends of the bench). On the other hand, the accuracy of the measurements in the midzone of the bench is always +1 °C.

If you want to determine the melting point of an unknown substance on the heating bench and with maximum accuracy, proceed in the following manner:

- 1 Perform a quick determination of the melting point of the substance to be examined.
- 2 Recalibrate the bench with a control substance for which the melting point is closest to the melting point determined in step 1.
- 3 Now remeasure the melting point of the substance to be examined on the calibrated bench. Record the melting temperature, which will be correct and precise this time.

1.9 Observation of the melting process

You can observe the processes on the heating surface with your naked eye. We suggesting using a 2X magnifying glass for a more precise observation.

1.10 Control substances and test substances

A set of control substance and test substance samples are supplied with the heating bench for:

- calibrating the reading device of the heating bench,
- determining eutectic temperatures,

The characteristics of the substances are listed below:

Name of substance	Melting Temp.
Acide adipique	150-153 °C
Benzile	94,0 °C
Acétanilide	113-115 °C
Benzoine	132-135 °C
Benzanilide	161-165 °C
Dicyandiamide	208-210 °C
Saccharine	225-227 °C
Acide stéarique	68-70 °C

Note: The manufacturer of the device may substitute other substances with similar melting points for the substances described in the above table.

All of these substances, except dicyandiamid, may be used as control substances.

Dicyandiamid may not be used as a control substance because as it melts it decomposes in such a way that its melting point lowers rather rapidly. This substance is supplied with the heating bench only as a test substance for determining eutectic temperatures.

1.11 Substances susceptible to decomposition

In the case of substances that decompose as they are heated, a melting point higher than that indicated in the specialized bibliography is often measured on the heating bench. One must take into account that the decomposition temperature depends to a large extent on the speed of the heating. In the specialized bibliography, attention is thus drawn to the fact that the melting points of substances that readily decompose are only comparable when the temperature elevation occurs at a uniform rate during the course of the measurement.

These values were determined with rapid temperature elevation and in general they are the most reproducible. However, it is necessary to heat the substance to its melting or decomposition temperature on the heating bench much more rapidly than when simply determining a melting point.

For substances that readily decompose, reading the temperature after around 10 seconds will give the best results.

1.12 Volatile substances

While determining the melting point, certain volatile substances disappear from the heating bench by sublimation.

It is nevertheless possible to determine the melting point for this type of substance because most of them reach melting temperature on the bench before sublimating. The reading must therefore be taken immediately after spreading the substance. Sublimation without melting beforehand only occurs on the bench in the case of extremely volatile substances, such as hexachlorethane and borneol.

1.13 Eutectic temperatures

It is also possible to determine eutectic temperatures with the heating bench. To do so, mix the sample with a nearly equal quantity of a suitable test substance on a slide. Cover the mixture with a second slide and rub it between the two slides. Then spread the mixture on the heating surface with one of the slides so that it forms a narrow band oriented in the lengthwise direction of the heating element. On this band, note the temperature at which the band becomes "humid" and take the reading after ca. 1 minute. Taking a reading sooner is only recommended for highly volatile substances.

The bulletin "Thermomicroscopy in the Analysis of Pharmaceuticals" by Kuhnert-Brandstatter discusses in detail the considerable advantages of measuring eutectic temperatures to identify organic substances. The tables in this bulletin list not only the melting point and other physical constants for each substance but also the eutectic temperatures determined on the heating bench that in general correspond to the microscopically-determined temperature.

However, there will be differences for substance that melt slowly. The temperatures indicated by the heating bench are higher in this case.

1.14 Hydrates

In the case of hydrates, the presence of crystallization water manifests itself in different ways on the heating bench. When they are heated above a specific temperature, numerous hydrates melt without changing their state; however, they then solidify in the form of anhydrous substances after the water evaporates.

Below the melting temperature of the hydrate, one can often observe a change in the hydrate crystals after spreading them out and after separation of the crystallization water. In numerous hydrates, this phenomenon occurs in a relatively wide temperature interval.

After separation of the crystallization water, it is possible to determine the melting point of the melting point of the anhydrous substance immediately by moving the crystals on the heating surface from the zone of the lower temperatures to the zone of the higher temperatures.

1.15 Melting point of mixtures

To determine the melting point of a mixture, place the mixture of the two substances on the heating surface and observe whether the separation line between the melted and non-melted fractions of the mixture is produced at the same temperature as for the two individual substances and if the melting point of the mixture is clearly defined.

For a direct comparison, one can place the individual substances and the mixtures side by side on the heating surface and proceed with the observation.

1.16 Softening temperature of synthetic powders and resins

With the spatula, place a small amount of the synthetic resin powder to be examined on the heating surface in the temperature zone of the anticipated softening point. After about 1 minute, remove the powder with a fine-bristled brush. This operation should be performed either transversally relative to the heating surface or in the direction of the higher temperatures (to the right).

Removing the powder in the direction of the lower temperatures (to the left) would result in soiling the heating surface. The pointer must be set on the point where the powder just shows traces of the beginning of sintering.

2 Care of the heating bench

In order to prevent damage to the heating surface, you must use the spatula provided with the apparatus to deposit, spread, and move the substance samples on it. Using any other tool such as a scalpel or other steel or hard metal tool may scratch the heating surface.

After the analysis, the residues of the substances must be removed from the heating surface. This can be done by simply wiping it with a dry cloth, with water or even with an organic solvent. We recommend using a metal polish for residues that stick strongly to the heating surface.

Avoid, however, any cleanser that might scratch the stainless steel surface.

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4 Use precautions

4.1 Breakdown or malfunction

Do not take the apparatus apart if it is broken or not working properly. Return it to our after-sales service department for repairs. Any intervention performed on the apparatus while it is under warranty will void the latter.

Do not turn the apparatus on after alterations have been made to the heating device or if it appears not to be working properly.

4.2 Heat and electrical hazards

During use, this apparatus has hot surfaces of which the temperature can reach 260 °C in some places. It is important not to touch the heating surface (2) while the apparatus is operating or while it is cooling down. Any contact between this surface and the skin can cause burns.

Wait until the heating surface of the apparatus has cooled before putting it in a box or cabinet.

Never leave the apparatus running unattended.

Do not cover the Kofler heating bench while it is in use.

Do not use this apparatus near flammable substances or those that may explode upon contact with heat.

Make sure that the apparatus has cooled completely before cleaning it, especially with flammable and/or volatile organic solvents.

We strongly recommend working under an exhaust hood when taking melting point measurements on substances that are capable of liberating toxic gases or vapours.

Always use an electric current corresponding to the characteristics of the apparatus.

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX France

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries