Mesure

Initio® 2

Measure

Initio® 2

Réf :
251 181

Français – p 1

English – p 8

Version : 7105

Pressiomètre Initio 2
Initio 2 pressure meter

SOMMAIRE

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
1. Description

The Initio® 2 pressure meter is an instrument that measures pressures between 200 and 2000 hPa. It can also save and store data, and restore them later onto a computer.

A WIFI module (ref 451 470) is also available to display and save the measurement on a tablet or a PC equipped with WIFI.

The Initio® 2 pressure meter is fitted with an analogue output and can be used with an ExAO interface with direct inputs through banana connectors or a Groomy® interface using the jack.

2. Use

It is recommended that the Initio® 2 battery should be fully recharged when it is used for the first time. This should be done using the μ USB charger (7) delivered in a separate package. In general, it is recommended that complete charge/discharge cycles should be carried out (wait until the instrument is fully discharged before completely recharging it).

2.1 Start up

- Press the power on button (1) on the left side of the instrument.

The measurement is displayed directly on the instrument screen.

- To switch the Initio 2 off, do a short press on the power on button (1) and validate pressing the knob (3).

2.2 Browsing

Initio® 2 has three buttons for browsing :

Select-OK (3):

The central knob (3) is used to browse through the Initio 2 menus. Do a brief press on the OK button on the middle of the knob to validate.

Return (2):

The "Return" button (2) changes from one type of display to the other from the main screen in the following order:

- Measurement
- Table
- Graph

Within a particular menu, it acts as a Back button.

Menu (4):

This button displays the Initio 2 main menu.

2.3 Data acquisition

The Initio® 2 has an internal memory to save measurements. There are two available modes:

- Manual acquisition

Acquisition is made point by point. The user saves each measurement manually (see below)

- Time acquisition

Acquisition is made as a function of time, for a duration input by the user.

2.3.1 Manual acquisition

- Press the Menu « button » (4).
- Use the knob (3), to choose **Acquisition** and validate by pressing briefly on the knob (3).
- Select **Manual**.
- To save a measurement, choose the **Save** icon (white circle) and press knob (3) briefly to validate.
- To stop acquisition, choose the **Stop** icon (black square) and press knob (3) to validate. Validate the choice in the dialogue box.

2.3.2 Time acquisition

- Press the Menu « button » (4).
- Use the knob (3), to choose **Acquisition** and validate by pressing briefly on the knob (3).
- Select **Time**.
- Use the knob (3) to select the number defining the duration and press the knob (3).
- Modify this duration by turning the knob (3) and then validate by pressing the centre of the knob.
- For a loop acquisition (restart the acquisition at zero after the defined time) select **In loop**.
- Validate the final choice by selecting **OK**.

Acquisitions as a function of time last between 60 seconds and 50 hours. The interval between two measurements is adjusted automatically.

2.3.3 Save data

Once the acquisition is complete (time or manual), you need to save it in order to retrieve points on a computer.

- Press the Menu « button » (4)
- Select **File** and validate
- Select **Save**

You can only save one acquisition. Select **Open** to open the most recently saved acquisition.

2.3.4 Transfer data to a computer

Once the acquisition is complete, you can transfer the data onto a computer for processing:

- Install the Initio® 2 software available at address www.jeulin.com/Initio_2
- Connect the Initio® 2 to the computer using the USB-μUSB cable
- Open the Initio® 2 software
- Load the saved acquisition by clicking on the icon:

Data are then displayed on the screen. Do a right click on the graph to access the different graphic tools.

3. Data display

You can display data in the form of a table or a graph at any time. Simply press the Return key (2). You thus change in sequence from the measurement to the table and then to the graph. You can continue the acquisition regardless of what is displayed.

The table or graph will be displayed in full screen after 5 seconds without taking any action. Turn the knob if you want to obtain the associated functions once again.

The arrows at the bottom right in the table are used to move through the rows of the table.

4. Analogue outputs


The Initio® 2 has analogue outputs on banana connectors and on the 3.5 mm jack.

4.1 Banana connectors

- Connect the Initio® 2 to an ExAO interface provided with direct inputs using banana cables
- Configure the direct input as follows:
 - o 0 hPa corresponds to -5V
 - o 2000 hPa corresponds to +5V

Example with the Foxy® console

- Start the Generalist
- Put the icon for a direct input  on the required axis
- Select the range $\pm 5V$

- Fill in the customisation tab as follows:

- Click on the diskette to save the parameter settings

The parameter settings are complete, the displayed magnitude is now the pressure with symbol P and its unit is hPa.

4.2 Jack

The 3.5mm jack is used for example to connect the Initio 2 to a Groomy® interface (ref 181 000). The output voltage varies between 0 and 5V. We have:

- 0 hPa corresponds to 0V
- 2000 hPa corresponds to 5V

5. WIFI communication

The WIFI module (ref 451 470) is used to display Initio® 2 data, for example on a tablet.

5.1 Installing the WIFI module

The WIFI module is installed when the Initio® 2 is switched off.

- Raise the Initio® 2 stand
- Remove the back by sliding it downwards
- Put the WIFI module into place, replacing the back
- Push the module all the way into place until it click fits. It must be possible to retract the stand without forcing.

5.2 Display on tablet

Once the WIFI module has been installed, the tablet has to be configured :

- Press the « Menu » button (4)
- Select Parameters
- Select Network


Information about the module is displayed on the screen.

- On the tablet, choose the name of the WIFI network displayed on the Initio® 2.
- Open an internet browser
- Enter the address displayed by the Initio® 2 in the browser address bar.

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
Enter the reference value of the pressure.

6.2 Initio® 2 parameters

The various information on the Initio® is accessed by pressing Menu and selecting Setup → Informations.

This page displays the serial number of the instrument and the software version, the level of the battery charge and the working language. The working language can be modified using the knob.

6.3 Update

On our website at the address <http://en.jeulin.fr/en-resources/english-resources/software/free-software.html>, the last update for your Initio® 2 is available. Its version is indicated in the designation.

In order to know the version of your Initio® 2, press the “Menu” button and then “Setup” “Informations”.

7. Technical characteristics

Pressure	200 hPa +2000hPa
Resolution	1 hPa
Precision	$\pm 2\% \pm 4\text{hPa}$
Maximum pressure	4000 hPa
Usage limitation	Non-corrosive gases. For liquids, use pressure sensor ref 242 013
Input	Coupling for 3 mm inside diameter tube
Analogue output	-5 V to +5 V on 4 mm safety banana sockets 0 to +5 V on 3.5 mm Jack
Communication	USB, WIFI optional
Power supply	Li-ion battery
Dimensions	95 x 123 x 44 mm

8. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

+33 (0)2 32 29 40 50