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Laboratory apparatus

**Ref :
253 012**

Français – p 1

English – p 7

Version : 7004

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+33 (0)2 32 29 40 50

1 Description

The complete kit for recovering gasses by water displacement is a device for performing the techniques of recovering the gas dissolved in carbonated water or soda water and then identifying it with limewater. All necessary components for the activity are included in this kit and it is therefore not necessary to resort to expensive, cumbersome support and holding devices. Furthermore, the kit is constructed of durable, resistant materials.


Photo 1

- 1 – Sampling syringe.
- 2 – Abduction tube with quick coupler.
- 3 – Double suction cup.
- 4 – Single suction cups with slots.
- 5 – Rubber stopper with quick coupler.
- 6 – Bubble tube with quick coupler.
- 7 – Cylindrical gas recovery bottle.
- 8 – Unbreakable polycarbonate tub.

2 Contents

The kit contains:

- 1 plastic 60 ml syringe with a female quick coupler,
- 1 PVC abduction tube with an outer/inner diameter of 6x4 mm, length 850 mm, with a male quick coupler,
- 1 PVC bubbling tube with an outer/inner diameter of 6x4 mm, length 120 mm, with a male quick coupler,
- 1 cylindrical, transparent polystyrene bottle, diameter. 35 mm, h 70 mm,
- 1 rubber stopper with a channel and a female quick coupler,
- 2 suction cups with split rings,
- 1 double suction cup,
- 1 unbreakable polycarbonate tub.

The various parts of the complete kit for recovering dissolved gasses by water displacement must not come in contact with organic solvents (chlorinated or aromatic solvents, ketones or esters) or strong acids and bases. These substances may cause the parts to deteriorate.

3 Instructions for use

3.1 Introduction

The performance of the activity requires the complete kit for recovering gasses by water displacement as well as a bottle of soda or carbonated water.

In order to achieve good results, it is preferable not to open the bottle until the beginning of the exercise in order to ensure that a sufficient quantity of gas is liberated during the gas liberation step.

The activity involves three steps that are described in detail in the following pages.

This kit was designed to enable the student to figure out for him or herself how to assemble the apparatus in order to successfully complete each step of the activity. Therefore, none of the elements composing this kit will guide the student towards the solution. The purpose of the following sections is only to provide examples of various assemblies that have been tested.


This kit is designed to perform the basic techniques for recovering and analyzing a gas dissolved in carbonated or effervescent water. The principle of recovering a gas by water displacement, for example, can itself be the subject of a study, e.g., of the air being introduced into the cylindrical gas recovery bottle with the syringe.

3.2 Liberating and recovering the gas

Fill the tub two thirds full with tap water. Assemble the apparatus as shown in *Schema 3*. Shake the bottle containing the carbonated water, observe the release of the gas and its accumulation in the gas recovery bottle.

Photo 2

Schema 3

- 1 – Rubber stopper with quick coupler.
- 2 – Bottle of carbonated or soda water.
- 3 – Cylindrical gas recovery bottle.
- 4 – Double suction cup.
- 5 – Abduction tube with quick coupler.
- 6 – Single suction cups with slits.
- 7 – Unbreakable polycarbonate tub.

3.3 Sampling the recovered gas

After recovering a sufficient quantity of gas, disconnect the abduction tube from the stopper in the bottle of carbonated or soda water. Connect the abduction tube to the sampling syringe after having first connected the syringe to the female quick coupler (*Schema 6*). By pulling the plunger of the sampling syringe, suction off the gas collected in the gas recovery bottle as shown below...



Photo 4

Schema 5

Schema 6

- 1 – Male quick coupler for the abduction tube.
- 2 – Female quick coupler.
- 3 – Sampling syringe...

3.4 Identifying the gas sample

Disconnect the gas-filled syringe from the abduction tube. Then connect the syringe to the male quick coupler of the bubbling tube. Insert the latter into a container of limewater and depress the plunger of the syringe in order to make the gas bubble through the limewater (*Photo 7, Schéma 8*). Using the syringe enables the quantity of gas injected into the limewater to be precisely controlled and the formation of the white precipitate to be readily observed.

Photo 7


Schéma 8

- 1 – Bubbling tube.
- 2 – Male quick coupler for bubbling tube.
- 3 – Limewater.
- 4 – Container.

4 Maintenance and use precautions

Clean the parts of the kit in soapy water. Never use organic solvents, as they may cause the parts to deteriorate.

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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Service gratuit * :
+ 33 (0)2 32 29 40 50

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Rubrique FAQ

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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