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Material properties

Réf :
253 081

Français – p 1

English – p 5

Version : 4101

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Dewar flask calorimeter

1 Description

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0 825 563 563*

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1 Description

The calorimeter is made up of a double-walled glass flask bright dipped in vacuum, covered with a plastic outer cover, an internal plastic vessel and a closing cover.

The calorimeter's cover consists of the following elements:

- 1 guided non-removable plastic stirrer, with a rectangular cross-section
- 1 $\varnothing$ 11 mm port for a thermometer or a probe
- 1 $\varnothing$ 6 mm guide for an expansion thermometer or a temperature probe
- 1 central circular opening that can be plugged by a cover to insert resistances
- Two 2 and 4 Ω resistances connected in series on sheathed aluminium rods, connected to an insulating cover with $\varnothing$ 4 mm safety plugs.

**Noted: if the introduction of the thermometer is difficult, to use a little Vaseline to facilitate the passage and avoid the risks of accidental breakage.*

To make temperature measurements, the calorimeter can be used with:

- A fairly long glass expansion thermometer*
- A temperature probe with a diameter less than 7 mm, such as the **Initio® Thermometer** part no. 251 040 or the **Initio® SA Thermometer with analog output** part no. 251 045 (a retaining clip for the temperature probe is provided for these 2 devices).

2 Technical characteristics

- | | |
|--|---------------------------------|
| • Effective volume: | 950 ml |
| • Internal dimensions of the vessel: | $\varnothing$ = 85 mm x 180 mm |
| • External dimensions: | $\varnothing$ = 130 mm x 270 mm |
| • Heat capacity:
(with protective vessel) | negligible, < 30 J/C° |
| • Power supply: | 12 V max DC |

Caution!!

The resistances must always be submerged when they are powered on.

The value of the resistances is only given for information. A measurement using an ohmmeter is necessary.

You must not put the double-walled glass flask bright dipped directly under cold water following experiment to avoid thermal shock.

3 Adiabatic properties

3.1 Thermal inertia


The melting of an ice cube in 250 ml of water lowers the temperature in the enclosure.

Ambient temperature: 20°C

Temperature in the calorimeter after fusion: 13.3°C

The initial temperature difference between the surroundings and the calorimetric enclosure is 6.7°C.


time	0	1 h	2 h	3 h	4 h	5 h	7 h	10 h
T (°C)	13.3	14.2	14.9	15.5	16	16.3	17	17.5
ΔT (°C)	0	0.9	1.6	2.2	2.7	3.0	3.7	4.2

As calorimetric experiments last only a few minutes, we can consider that the calorimeter is perfectly adiabatic (insulated) for a temperature difference of $\pm 5^\circ\text{C}$ between the ambient and internal temperatures.

3.2 Heat loss

We have established the heat loss curve of the calorimeter with 200 g of water, the ambient temperature being 25°C.

Between 85 and 72 °C, the decrease in temperature is 1.8 °C per minute, on an average.

Below 70°C, the decrease remains almost constant at 0.27 °C per minute over the entire measuring range.

The loss of heat therefore remains low considering the experimental parameters:

- We are located more than 20 °C above the ambient temperature
- The volume of water represents only 1/5th of the calorimeter's effective volume.

4 Experiments

4.1 Heat capacity of the calorimeter


We pour 250 mL of water (cold) at 20.7°C in the calorimeter

We then pour 250 ml of hot water at 55.2°C in the calorimeter. The final temperature is 37.7°C.

According to the following formula, it is possible to determine the heat capacity of the calorimeter knowing that $C_{\text{water}} = 4187 \text{ J/Kg/}^\circ\text{C}$.

$$C = \frac{m_2 c_{\text{water}} (\theta_2 - \theta_{\text{eq}})}{(\theta_{\text{eq}} - \theta_1)} - m_1 c_{\text{water}}$$

We show that the heat capacity of the calorimeter with the plastic protective jacket is 32.05 J/°C.

4.2 Variation of the quantity of heat with the intensity

Place a mass m (about 300 g) of water in the vessel of the calorimeter.

Immerse the resistance and the accessories of the calorimeter and record the temperature of the assembly on the thermometer.

Prepare a circuit that includes a current source, an ammeter, a rheostat and the resistance in series.

Switch on the power and start marking time. Note down the current on the ammeter.

Let the current flow for a time t . Repeat the experiment several times by changing the current.

The quantity of heat is expressed by: $Q = m.C.(\theta_2 - \theta_1)$

with:

- m the mass of water
- C the specific heat of water
- θ_1 the initial temperature
- θ_2 the final temperature

If I_1 , I_2 , and I_3 are the different intensities measured, we can verify that the ratios Q_1 , Q_2 , etc. of the quantities of heat released are proportional to the squares of the measured currents.

4.3 Variation of the quantity of heat with time

Take several temperature readings θ_1 at regular intervals.

Calculate the quantities of heat released: $Q = m \cdot C \cdot (\theta_2 - \theta_1)$.

We observe that the quantities of heat released are proportional to the time intervals separating the start of the experiment from each of the readings.

4.4 Variation of the quantity of heat with resistance

Repeat the experiment under the same current and time conditions by successively using the 3 resistance values R_1 , R_2 , R_3 .

4.5 Determining the ratio W

Let the current flow for a time of t seconds. The quantity of energy released is $W = RI^2 \times T$, in Joules. By calculating the ratio W/Q , you will find a number that is close to 4 ($J = 4.18$).

4.6 Measuring the specific heat of liquids

An interesting application of Joule's law is to measure the specific heat of various liquids.

We obtain an increase in temperature θ_1 .

Replace water with another liquid and repeat the experiment with the same resistance, the same current and for the same duration. The rise in temperature is θ_2 .

The specific heat of the liquid M may be calculated using the formula:

$$C / m = \theta_2 / \theta_1.$$

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX France

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux