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0 825 563 563*

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1. Description

This test bench is designed to determine the materials thermal resistance. This device is based on the unit used for the standard measurement of the material thermal conductivity.

The sample with unknown thermal conductivity is placed between two metal blocks at a known and chosen temperature change. The forced steady state is stabilized by the energy coming from the upper block.

The interface between the device and the student allows the latter to manipulate all variables involved in the equation implicated by the unit. When the manipulation starts, the steady state which is necessary for obtaining correct values, is directly verified on the device.

The results are also read on the test bench through a third menu. It tracks the thermal resistance, the power supplied by the device equal to the heat flow, the thermal conductivity as well as the voltage and current.

The removable part of the device is designed to insert samples with a maximum thickness of 1 cm. They are less than the maximum limit thick and therefore allow material combinations in order to verify the status of addition of thermal resistances.

2. Principle

2.1 Technical Specifications

Batch of samples:

- depron plate of about 3 mm thick
- plywood plate of about 3 mm thick
- glass plate of about 3 mm thick
- two PVC plates of which one is porous (white).

12V, 5A power supply.

2.2 Theoretical Aspects


Figure: theoretical diagram

2.2.1 General Case on Flat Wall

Application of Fourier's law to the steady state thermal conduction. A flat and homogeneous wall is considered.

In a Cartesian system, Fourier's rapport thus applied is expressed by:

$$\Phi = \lambda \overrightarrow{\text{grad } \theta} \cdot \vec{n}$$

With \vec{n} unit vector oriented according to the flow
 θ : Temperature

Note: λ is always positive. The sign « - » translates the heat flows orientation towards low temperatures.

In our case, the wall is homogeneous, then: $\lambda = cste$

Also, the wall is flat, giving: $S=S'$ (see diagram).

Consequently, we can write:
 Equation 1

$$\Phi = \lambda S \frac{d\theta}{dx}$$

In a steady state case, the temperatures θ_1 and θ_2 are considered equal to constants; and θ is a constant physical quantity of the system.

By applying this so-called Fourier's law, we can determine the relation involved.

2.2.2 Determination of the Relation Involved

Application of the Fourier's law given by the above equation 1:

$$\Phi = \lambda S \frac{d\theta}{dx} \Rightarrow \Phi dx = -\lambda S d\theta$$

Applying the Fourier's law:

$$\Phi \cdot d(x) = -\lambda \cdot S \cdot d(\theta)$$

$$\Phi \cdot \int_{x_1}^{x_2} d(x) = -\lambda \cdot S \cdot \int_{\theta_1}^{\theta_2} d(\theta)$$

$$\Phi \cdot |x|_{x_1}^{x_2} = -\lambda \cdot S \cdot |\theta|_{\theta_1}^{\theta_2}$$

$$\Phi \cdot (x_2 - x_1) = -\lambda \cdot S \cdot (\theta_2 - \theta_1)$$

$$\Phi \cdot (x_2 - x_1) = \lambda \cdot S \cdot (\theta_1 - \theta_2) \text{ where } (x_2 - x_1) = e \text{ (wall thickness)}$$

Being $\Phi \cdot e = \lambda \cdot S \cdot (\theta_1 - \theta_2)$

The flow is then equal to:

$$\Phi = \frac{\lambda \cdot S \cdot (\theta_1 - \theta_2)}{e}$$

with $\frac{e}{\lambda \cdot S}$ (wall thermal resistance).

Then:

$$\Phi = \frac{(\theta_1 - \theta_2)}{R}$$

With: θ_2 and θ_1 in K (Kelvin)
 S in m^2
 e in m
 R in $K.W^{-1}$
 λ in $W.m^{-1}.K^{-1}$
 ϕ in W

2.3 Initiation and Navigation

2.3.1 Initiation


The device is directly initiated when the power cord is connected to the socket on its left side.

Once connected, the screen illuminates.

2.3.2 Caution

The upper and lower jaws must be set on each other before each measurement. The aim is to verify if the temperatures of these jaws are equivalent and close to the ambient temperature using the **grph** tab.

2.3.3 How to place the sample?

Simply lift the upper part of the thermoblock (see photo) to correctly place the sample in the jaws.

Place it in the dedicated space and then replace the jaw on the sample.


2.3.4 Navigation

The human-machine communication is carried out using the above described functionalities.

A rotary encoder allows switching from one upper tab to another.

Prm tab

The **Prm** tab for "parameters" allows entering the different parameters and starting the experiment.

Press the SET button to access the menu.

The output voltage varies from -5V to 5V for the following ranges :

T1: 0 to 40 °C

T2: 0 to 40 °C

Ta: 0 to 40 °C


U: 0 to 15 volts

I: 0 to 1.5 ampère

P: 0 to 15 watts

Rth: 0 to 100 K.W-1

λ: 0 to 5 W.m-1.K-1

The rotary encoder allows navigating between the different parameters provided. Press SET to change the selected parameter value.

The desired value is obtained with the rotary encoder.

Press SET to validate it.


In the **Prm** tab, select the Exp. line and use the rotary encoder to select "On" and then validate by pressing SET.

In a different tab, press and hold SET and select "On", then validate by pressing SET.


Manipulation starts.

Use either the Back button or the rotary encoder for a different tab.

Grph tab


The **Grph** tab for "Graph" allows following the temperatures progress up to the expected steady state. T_a represents the ambient temperature, T_1 that of the upper block and T_2 , the lower block.

To monitor the temperature evolution in full screen: Set to T_1 or T_2 and validate using the SET button. Switch from one temperature to another by pressing SET. Go back to the menu by pressing back.

The third navigation tab **Res.** for "Results" allows reading the values when the system is in steady state.


All results are given in the international system units.


U represents the voltage, I the current, P the power which is also the flow, R_{th} the thermal resistance and λ the thermal conductivity.

Time is expressed in minutes / seconds.

The fourth and last tab **Opt** for "Option" is used to set the analog outputs.

If you wish to associate a result to an analog output, press SET to enter the tab menu and use the rotary encoder to set S1, S2 or S3. Then press SET.

Use the rotary encoder to select the desired parameter, then press SET to validate

3. Objectives

The thermal resistance test bench allows verifying one of Fourier's relations in thermodynamics as seen above.

The maintained steady state allows working on a single variable, the thermal resistance, and verifying the tested material thermal conductivity.


The system also provides the opportunity to verify the thermal resistances addition during the layering of several materials, or the relation between thickness and thermal resistance of a material.

4. Manipulations


4.1 Materials Classification

We wish to verify the thermal conductivity of different supplied materials.

The most repetitive results are obtained by applying the following program:

Each sample thickness will be measured at the ruler for greater accuracy.

Values can be read after the steady state has been reached for 2 minutes. It is recommended to wait 5 minutes of stabilization.

In general, suppliers declare thermal conductivities included within a measuring range. The more commonly found are the following:

Material	Depron	PVC	Plywood	Glass
Thermal conductivity (W.m ⁻¹ .K ⁻¹)	0.03 to 0.04	0.12 - 0.25	0.12 - 0.15	0.8 – 1.3

Both PVC samples supplied show the need for manufacturers to declare a value range rather than a single value.

We can also note this conclusion:

$$\lambda_{\text{Depron}} < \lambda_{\text{Plywood}} < \lambda_{\text{PVC}} < \lambda_{\text{Glass}}$$

Then:

$$R_{\text{Depron}} > R_{\text{Plywood}} > R_{\text{PVC}} > R_{\text{Glass}}$$


4.2 Thermal Resistance Improvement

In this experiment, we wish to verify the resistance increase by adding material layers.


We know that: $R_{1+2} = R_1 + R_2$ where $e_{1+2} = e_1 + e_2$

Example of application:


A 3 mm porous PVC plastic thermal resistance is evaluated.

A 3 mm expanded PVC plastic thermal resistance is evaluated.

The two material layers are added to verify the above relation. The resistance of a wall composed of two material layers is then evaluated. The total thickness is 6 mm.

We note an increase in the thermal resistance when adding both material layers.

5. After-Sale Services

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux