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Étude des combustions

Chemistry teaching

Study of combustion

Réf :
253 140

Français – p 1

English – p 18

Version : 6106

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Combustion chamber

Contenu

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NOTES

Contents

1	Description	16
1.1	General	16
1.1.1	Composition and description	16
2	Precautions for use	17
2.1	Warning	17
2.2	Dangerous reactions	17
2.3	Use of oxygen	17
2.4	Use of the chamber	18
2.5	Protection of persons	18
3	Rapid start-up	18
3.1	Electric ignition	18
3.2	Rapid connections	18
3.3	Inlet and outlet valves	19
3.4	Connecting the injection syringe	19
3.5	Filling the chamber with oxygen	19
3.5.1	Water displacement method	19
3.5.2	Air displacement method	20
4	Experiments	20
4.1	Combustion in oxygen	20
4.1.1	Combustion of carbon (charcoal)	20
4.1.2	Combustion of iron (steel or iron wool)	23
4.2	Analysis of smoke	25
5	Upkeep and maintenance	27
5.1	Cleaning	27
5.2	Maintenance	27
5.2.1	Wearing parts	27
5.2.2	Heat screen	28
6	After-Sales Service	28

1 Description

1.1 General

The combustion chamber is a device for carrying out combustion in oxygen. It also allows investigating the combustion reactions of carbon and iron.

The chamber makes it possible to filter and analyse the smoke generated by the combustion of Armenia paper in air. The body of the chamber is in glass for easy observation of the phenomena which occur during combustion reactions.

The cover and its fittings incorporate many accessories to facilitate the experiments.

The accessories supplied with the combustion chamber make it easy to take samples for analysis of the reaction products.

The combustion chamber is also delivered with a solid cover used when filling the chamber with oxygen, in advance preparation for experiments.

1.1.1 Composition and description


Photo 1

- 1 – Cover with fittings
- 2 – Body of chamber in glass
- 3 – Solid cover
- 4 – Injection syringe
- 5 – Rapid male connector
- 6 – Rapid female connector

- 7 – Sampling outlet valve
- 8 – Sockets for connecting electric ignition
- 9 – Injection inlet valve
- 10 – Sample holding rods
- 11 – Removable sample holding clips
- 12 – Sampling syringe

2 Precautions for use

2.1 Warning

This appliance is designed for combustion reactions in the presence of oxygen, with a risk of fire, violent reactions or burning in the event of improper use.

The company JEULIN cannot be held liable for accidents if the combustion chamber or one of the accessories described in these instructions has been modified or converted by the user.

Also, the company JEULIN cannot be held liable in the event of an accident caused by non-compliance with the safety instructions, as set out below.

To guarantee the safety of users, it is vital to comply with the following instructions.

2.2 Dangerous reactions

The combustion chamber is designed **exclusively** for the operations described in these instructions.

It is formally prohibited to carry out the experiments and reactions described below. Non-compliance with these instructions may cause serious accidents.

No substances in any form whatsoever other than those specified in these instructions must ever be inserted in the combustion chamber.

- Never combust reducing metals such as aluminium, zinc, magnesium, calcium, sodium etc.
- Never carry out experiments using detonating mixtures of gas, liquids or solids (hydrogen/oxygen mixtures in stoichiometric proportions, explosive powders and solids etc.).
- Never carry out combustion of flammable solids, liquids or gases or those which can release toxic or corrosive gases.
- Never carry out any reaction involving gases other than oxygen (chlorine, hydrogen etc.).

2.3 Use of oxygen

Oxygen is a dangerous gas, particularly if highly concentrated. Combustible materials (wood, paper, plastic and metal) in an oxygen-rich atmosphere (> 20% by volume) may catch fire and burn fiercely.

All the usual precautions must be adopted to avoid fires.

- Never use the combustion chamber close to flammable or combustible materials.
- Never work on or near a combustible surface (table, wooden furniture etc.).

The combustion chamber is designed to carry out reactions in pure oxygen or the oxygen contained initially in the chamber. Combustion reactions stop when all the oxygen in the chamber at the start of the reaction is consumed. The quantity of heat released by combustion reactions is thus limited.

For this reason, the combustion chamber should never be used with an oxygen supplying device. Non-compliance with these instructions may cause fires and accidents.

Filling the chamber requires the use of pressurised oxygen bottles or cartridges which may be dangerous to handle. To prevent the risks intrinsic to this equipment, you are advised to read the Safety Data Sheets carefully before using these products.

The Safety Data Sheets can be downloaded free from the site:

www.jeulin.fr

2.4 Use of the chamber

Never use the combustion chamber after modifying it or linking it to any device other than the accessories supplied with the product.

Never use a combustion chamber with signs of wear or damage.

All parts which appear worn or damaged must be replaced before using the combustion chamber. The glass tank must be replaced if it has impact marks or cracks. Refer to the upkeep and maintenance section.

2.5 Protection of persons

When using the combustion chamber, it is necessary to comply with the safety instructions on personal protection.

- Wear a cotton overall during use
- Tie back your hair
- Operate on a clear, uncluttered work surface
- Access to and the area around the work bench must be kept clear.

3 Rapid start-up

3.1 Electric ignition

The cover and its fittings incorporate an electric ignition device by short circuit which triggers the combustion reaction of iron wool with the oxygen contained in the chamber, without any need to open it up.

The sample holding rods are linked to safety sockets, diameter 4 mm for establishing the short circuit. The sockets may be linked either to a very low voltage power supply ($U_{\text{max}} = 12$ V) or to a 4.5 V electric battery.

Use of the electric ignition is explained in paragraph 4.1.2.

3.2 Rapid connections

The combustion chamber is equipped with a rapid connection system which ensures hermetic sealing by simple engagement. The rapid connections are assembled as indicated below (*Diagram 2*).


Diagram 2

It is not necessary to apply force to the rapid connections to ensure the assembly is leak-tight.

Disassembly is by pulling the rapid connections apart. Do not pull on the flexible hoses, since this could damage the equipment.

3.3 Inlet and outlet valves

The injection inlet and the sampling outlet of the chamber are fitted with valves to isolate the contents of the chamber or to allow inserting the reagent or to sample the gases it contain.

Opening and closing of these valves is controlled simply by operating them as indicated below (*Diagram 3*).


Diagram 3

3.4 Connecting the injection syringe

The injection syringe is connected to the injection inlet by simple engagement as indicated below (*Diagram 4*).

Filling the chamber with oxygen

There are two possible methods for filling the combustion chamber with oxygen.


Diagram 4

3.4.1 Water displacement method

To fill the combustion chamber with oxygen, fill the glass tank with water and turn it over in a crystallizer or a tank, itself filled with water (*Diagram 5*).

Insert the oxygen supply pipe in the tank and open the pressure relief valve.



Diagram 5

When the oxygen has driven out the water it contains, close the oxygen supply, close the tank with the solid cover, keeping it upside down in the crystallizer or tank. Then remove the chamber from the crystallizer or tank.

3.4.2 Air displacement method

Open the chamber and insert the oxygen supply hose (Diagram 6).

Open the pressure release valve and allow the oxygen to escape for a few seconds so that it replaces the air initially contained in the chamber.

Cut off the oxygen supply and then close the chamber with the solid cover after removing the oxygen supply pipe.

This method, simpler to use than the water displacement one, has the disadvantage of not allowing precise monitoring of the oxygen concentration of the chamber after filling. However, it avoids the presence of water - inevitable when using the first method.

4 Experiments

4.1 Combustion in oxygen

The combustion chamber allows combusting carbon (stick of charcoal) and iron (iron wool or steel wool) in oxygen.

As stipulated in chapter 2, it is formally forbidden to perform combustion of materials other than those referred to above.

4.1.1 Combustion of carbon (charcoal)

This experiment allows studying the combustion of carbon (in the form of charcoal sticks) in oxygen, observing the reaction and describing the products.

a) Equipment required



Designation	P/N
Combustion chamber	253 140
Small beaker or glass with base	-
Stick of charcoal	-
Lime water	107 341
Oxygen bottle or cartridge	120 015

b) Procedure

The chamber must have previously been filled with oxygen (*paragraph 3.5*).

The valves on the cover with fittings must be in the closed position.

Attach a strip of charcoal, two to three centimetres in length, in one of the cover clips (*Diagram 7*). It is not necessary to place the two clips on the sample holders for this experiment.

Observe combustion of the strip of charcoal: it catches light as soon as it is placed in the oxygen in the chamber (*Diagram 8*).

Combustion is complete and there is almost no residue.

Diagram 8

c) Describing the combustion products

The combustion chamber and its accessories allow describing the combustion products in several different ways.

It is possible for the students themselves to suggest an analysis protocol and use it. The students can first of all observe that combustion of charcoal does not release any water vapour, because there is no condensation on the walls of the tank. This is true if the charcoal and the chamber are dry before the experiment (no humidity in the charcoal and no water on the walls of the chamber).

Characterisation of carbon dioxide is performed during the lime water test.

Analysis by sampling:

Characterisation of carbon dioxide by sampling gas present in the chamber.

- ① Open the two valves. Fit the rapid connection of the sampling syringe to that of the chamber sampling outlet (*Diagram 9*).
- ② Suck up a sample of the gases contained in the chamber with the sampling syringe (it is necessary to fill the syringe completely to perform the lime water test).
- ③ Disconnect the end piece of the sampling syringe from the combustion chamber and place the contents of the syringe in a beaker or glass with a base containing a small amount of lime water.


Diagram 9

The presence of calcium carbonate precipitate identifies the combustion product of carbon in oxygen as carbon dioxide.

Analysis in situ:

Characterisation of carbon dioxide by injecting lime water in the chamber.

- ① Open the two valves. Connect the injection syringe previously filled with lime water to the combustion chamber injection inlet (Diagram 10).
- ② Inject the contents of the injection syringe into the combustion chamber.
- ③ Gently shake the chamber until the precipitate can be observed.


Diagram 10

Diagram 11

As with the previous method, the precipitate of calcium carbonate identifies the combustion product.

The in situ analysis can also be performed by pouring lime water into the chamber before inserting the incandescent charcoal. Precipitate of calcium carbonate appears during combustion of the charcoal in oxygen.

d) End of the experiment

At the end of the experiment open the chamber by unscrewing the cover with fittings (Diagram 11).

Parts in contact with the charcoal (crocodile clip and sample holding rod) are heated to a high temperature during combustion, take care to cool them under cold running water before touching them.

4.1.2 Combustion of iron (steel or iron wool)

This experiment allows studying the combustion of iron (in the form of steel or iron wool) in oxygen, observing the reaction and the products.

a) Equipment required

Designation	Part No.
Combustion chamber	253 140
Power supply or 4.5 V battery	-
Two banana jack cables	-
Stainless steel or iron wool	107 464
Oxygen cartridge or bottle	120 015
Sand or water	-

b) Preparing for the reaction

First fill the chamber with oxygen (*paragraph 3.5*).

The valves on the cover with fittings must be in the closed position.

It is essential to place a layer of sand or water on the base of the chamber before carrying out the experiment for combustion of iron in oxygen.

In fact, drops of molten iron oxide are sprayed during the reaction and may break the chamber if they come into contact with the walls.

① Attach a few strands of stainless steel wool to the clips in the cover with fittings (*Diagram 12*). Ensure you do not attach too many stainless steel wool strands to each clip to ensure the electric ignition works properly.

② ③ Remove the solid cover and rapidly fit the cover with fittings and the strands of stainless steel wool.


Diagram 12


c) Igniting and observing the reaction

① Have an electricity supply to hand. A 4.5 V battery or other source of supply can be used.


The voltage supplied must not exceed 12 V.

- ② Connect the ignition sockets to the power supply with the cables fitted with banana jacks (*Diagram 13*).

- ③ Close the power supply circuit, or start up the power supply.


- ④ Observe the combustion reaction of iron in oxygen.

d) Preparation for the reaction

- ① At the end of the reaction the combustion products (iron oxides) will fall to the bottom of the chamber into the water or layer of sand (*Diagram 15*).
- ② To observe the combustion products, open the chamber by undoing the cover with fittings (*Diagram 15*).

! Parts in contact with the stainless steel wool (crocodile clip and sample holding rod) and the combustion products are heated to a high temperature during combustion, take care to cool these under cool running water before touching them.

4.2 Analysis of smoke

Another possible use of the combustion chamber is to analyse smoke. It is possible to produce smoke inside the chamber using a smoke generating product and then sample the smoke which can be studied by filtration.

The main advantage of the combustion chamber is confining the smoke and avoiding it being dispersed throughout the classroom.

a) Equipment required

Designation

Combustion chamber
 Armenia paper
 Filter paper or cotton

Part No.

253 140
 107 163
 -

b) Operation


The valves of the cover with fittings must be in the closed position.

Attach a piece of Armenia paper to the clips in the cover with fittings (*Diagram 16*).

- ① Light the sheet of Armenia paper using an electric heater or cigarette lighter and put out the flame by blowing on it.
- ② Rapidly fit the cover with fittings with the piece of burning Armenia paper in the chamber.

This experiment must be performed in air, do not fill the chamber with oxygen.

- ③ Observe combustion of the Armenia paper in the chamber. The smoke accumulates gradually in the chamber until its contents are totally opaque. *Diagram 17*).

Combustion of the Armenia paper ends naturally on account of lack of oxygen.


Diagram 17

c) Sampling and analysis of smoke

- ① Open the two valves. Assemble the rapid connection of the sample syringe to that of the chamber sampling outlet (*Diagram 18*).
- ② Suck up a sample of the smoke contained in the chamber using the sampling syringe.


Diagram 18

- ③ Disconnect the end of the sampling syringe from the combustion chamber and empty the contents of the syringe through a piece of cotton or filter paper (*Diagram 19*).


Diagram 19

- ④ Then observe the deposit on the piece of cotton or filter paper.

d) End of the experiment

At the end of the experiment open the chamber by loosening the cover with fittings (*Diagram 20*).

The parts in contact with the piece of Armenia paper (crocodile clip and sample holding rod) are heated to a high temperature during combustion. Take care to cool them under cold running water before touching them.

To avoid any discomfort caused by releasing the combustion smoke when opening the combustion chamber, it is advisable to perform this operation under a suction hood or close to an open window.


Diagram 20

5 Upkeep and maintenance

5.1 Cleaning

The combustion chamber and its accessories can be cleaned with water and a dishwashing-type detergent.

Dry the metal parts carefully after cleaning to avoid any rust forming.

During the smoke analysis experiment, tar deposits from condensation of smoke may occur. The deposits of tar can be eliminated (brown traces) by cleaning the syringe, the hoses and the sampling valve using denatured ethanol. Take a pipette full of denatured ethanol and run denatured ethanol into the soiled parts. The tar enters into solution with the ethanol and is eliminated.

Never perform a combustion experiment immediately after cleaning with ethanol or before you are sure there are no traces of liquid or ethanol or ethanol vapour in the chamber and its accessories.

Failure to comply with this instruction may result in fires, or even explosions.

Under no circumstances can other organic solvents be used, such as aliphatic hydrocarbons, aromatics or halogens or oxygen solvents (esters, ketones etc.). These products can irreversibly damage the plastic components and the rubber seals.

5.2 Maintenance

5.2.1 Wearing parts

The combustion chamber incorporates various wearing parts which must be replaced after several usages.

Reactions in oxygen are highly exothermic and parts in contact with the combustion reagents are subject to harsh conditions.

So these are designed to be replaced easily and economically.

a) Sample holding clips

The removable sample holding clips must be replaced at regular intervals.

b) Glass body of chamber

During the reaction with combustion of iron wool in oxygen, the glass body of the chamber is liable to be sprayed with molten iron oxide and then crack. It is essential to replace it if this occurs.

c) Leaktightness seal

The inlet and outlet valves are assembled on the cover with fittings using rubber seals. These must be replaced when they are damaged. To replace them proceed as shown below (*Diagram 21*).


Diagram 21

Assembly of the new seal is performed as indicated below (*Diagram 22*).


Diagram 22

Leaktightness seals (pack of ten)

PN: 253 045

5.2.2 Heat screen

Never remove the aluminium heat screen lining the cover with fittings. This protects the cover from radiation and projections which occur during the combustion reactions.

Failure to comply with this instruction may result in fire and accidents.

6 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX France

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux