Thermodynamique

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Thermodynamics

Infrared radiation

Réf :
253 174

Français – p 1

English – p 31

Version : 0000

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Infrared test bench

Sommaire

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Summary

Introduction	28
1.Theoretical part: Principle of infrared thermography	28
1.1 Definition	28
1.2 Principle	28
1.3 Infrared radiation	29
1.3.1 The electromagnetic spectrum	29
1.3.2 Elements of radiometry	30
1.4 Thermography	32
1.4.1 Measurement reliability	32
1.4.2 Detection	32
1.4.3 Thermal resolution.....	33
1.4.4 Thermal accuracy.....	33
1.4.5 Spacial resolution	34
1.4.6 Measurements conditions	34
1.5 Thermal image analysis techniques.....	35
1.5.1 Thermal framing	35
1.5.2 Isotherm	35
1.5.3 Palettes	36
2. Practical part.....	38
2.1 Comparative study of the emissivity of different materials	40
2.2 Calibrating emissivity on the thermal camera.....	42
2.3 Determining the emissivity of an unknown material	43
2.4 Influence of the material's surface condition on its emissivity.....	45
2.5 Study of the influence of angle on measurement	46
2.6 Study of the spatial resolution parameter using the slot.....	47
2.7 Special case : transparent material.....	49
2.8 Demonstrating the reflection of materials, reflected temperature	50
3. After sales service	51

Introduction

In this practical course, we'll be carrying out a number of manipulations and measurements using an infrared camera. The aim is to characterize the thermal behavior of materials, to adjust the camera to the case under study, and to take a critical look at the measurements made using infrared thermography..

1. Theoretical part: Principle of infrared thermography

Infrared thermography is a non-contact temperature measurement technique used to diagnose the state of an electrical, mechanical or thermal system in operation. It is widely used to optimize maintenance tasks without interrupting production flow, and to minimize maintenance costs..


1.1. Définition

Infrared thermography is a technique that uses appropriate equipment to obtain a thermal image of a thermal scene observed in the infrared spectral range. In other words, infrared thermography makes it possible to visualize temperature gradients that cannot be seen in the visible range.

1.2. Principle


Any object whose temperature exceeds absolute zero (i.e. $0^\circ\text{Kelvin} = -273.15^\circ\text{C}$) emits electromagnetic radiation. This radiation depends directly on the object's surface temperature.

The human eye can “see” heat at high temperatures: sun (6000°C), tungsten filament (2200°C), red iron (600°C), but for low-temperature objects ($<500^\circ\text{C}$), which emit little radiation, the human eye is inadequate (see

Planck's Law graph). The use of infrared radiation sensors makes up for this shortcoming: infrared cameras, which capture radiation between 3.6 and 13 microns in wavelength, convert this radiation into a thermal image, and allow temperatures to be integrated.


An infrared camera is a radiometer whose input signal is radiation power (W/m^2), and whose output signal is an electrical voltage proportional to the input signal. This voltage is translated into visual and digital indications. Thermographic cameras are designed to provide an image, called a thermogram, of static or dynamic phenomena in the infrared spectral range..

1.3. Infrared radiation

1.3.1. The electromagnetic spectrum

The electromagnetic spectrum is divided more or less arbitrarily into ranges called bands, depending mainly on the methods used to emit or detect radiation. There is no

fundamental difference between radiation in the different bands of the electromagnetic spectrum. They all respect the same laws, with nuances resulting from differences in wavelength.


Visible light, radio waves, TV and X-rays are all electromagnetic radiation. The visible range covers wavelengths from 0.4 to $0.8 \mu\text{m}$. The infrared band [figure 7.3] extends from 0.8 to $1000 \mu\text{m}$. In infrared thermography, we generally work in a spectral band that extends from 2 to $15 \mu\text{m}$ and more specifically in the $2\text{-}5 \mu\text{m}$ and $7\text{-}15 \mu\text{m}$ windows.

Wavelengths are most often expressed in μm (micrometers) or nm (nanometers), but an older unit is still often used: the Ångström (\AA). The relationship between these different units is: $10,000 \text{ \AA} = 1,000 \text{ nm} = 1 \mu\text{m}$, or $1 \text{ \AA} = 1.10^{-9} \text{ m} = 0.1 \text{ nm}$.

1.3.2. Elements of radiometry

a. Principle

Incident radiation refers to all radiation striking an object from outside. In the figure, the incident radiation, noted W_{INCID} , is all the radiation striking the target object from one or more sources. A certain proportion of the radiation, noted W_{α} , will always be absorbed, and the target object will retain its energy. A certain amount, W_{ρ} , will be reflected. The latter will not affect the target object in any way. Finally, a certain proportion of the radiation, W_{τ} , may pass through the target object. Like the quantity reflected, it will not affect the object. If we write this as a mathematical equation, we obtain:


$$W_{\alpha} + W_{\rho} + W_{\tau} = W_{\text{INCID}} = 100\%$$

An object has a certain capacity or ability to :

- Absorb: this is called absorptivity, α
- Reflect: so-called reflectivity, ρ
- Transmit: transmissivity, τ

The sum of the three is always equal to 1: $\alpha + \rho + \tau = 1$.

The resulting radiation captured by an infrared camera is the sum of all the radiation leaving the surface of an object, whatever its source. It comes from three types of source. In the figure, we have three sources of radiation: the target object itself, a source in front and a source behind. If we look at the target object from the right, the resulting radiation is a combination of the radiation from the target object itself, the radiation from the reflection on the object and the radiation from the source behind.

radiation from the reflection on the object of the front heat source (located on the right) and radiation from the rear heat source (located on the left) passing through the target object.

b. Emissivity

In infrared thermography, the most important part of the resulting radiation is generally the emitted part. An object has a certain capacity or ability to emit, characterized by its emissivity ε .

This is denoted ε and varies with temperature and wavelength.

$$\varepsilon = \frac{\text{Quantity of rays emitted by the real object}}{\text{Quantity of rays emitted by the black body}}$$

b. Notion de corps noir

A black body is an ideal radiating element that absorbs 100% of incident radiation, meaning it neither reflects nor transmits any radiation. This does not exist in reality, as there is always a little something reflected, for example. For a black body, the emissivity $\varepsilon = 1$ (et $\rho + \tau = 0$).

c. Notion of gray body

The term “gray body” is an abstract notion that enables relative measurements to be made. The basic assumption is that in the 3 to 5 μm band, the object being measured is a gray body.

d. Glossy bodies

hiny bodies, such as mirrors, are highly reflective. Referring to Drapper's equation, we have : $\varepsilon(\lambda) = 0$.

1.4. Thermography

1.4.1. Measurement reliability

One of the main problems with infrared thermography is the reliability of the measurements. This depends on :

- control of absolute temperature drift ;
- control of noise, i.e. the disturbance generated by rising sensor temperatures.

To avoid this phenomenon, the sensor can be cooled, but the latest-generation cameras no longer need to be cooled, as a software system for correcting measurement drift is included in the device. Regularly (every 30 s), temperature sensors integrated into the camera provide temperature information. The processing unit can then take this information into account to recalibrate the measurement.


1.4.2. Detection

One of the main problems with infrared thermography is the reliability of the measurements. This depends on :

- control of absolute temperature drift ;
- control of noise, i.e. the disturbance generated by rising sensor temperatures.

To avoid this phenomenon, the sensor can be cooled, but the latest-generation cameras no longer need to be cooled, as a software system for correcting measurement drift is included in the device. Regularly (every 30 s), temperature sensors integrated into the camera provide temperature information. The processing unit can then take this information into account to recalibrate the measurement.

Infrared cameras are equipped with a vanadium oxide matrix containing a multitude of sensors that transform the radiation received into electrical voltage.

1.4.2. Thermal resolution

This is the smallest temperature variation discernible by the device. It's also known as the M.R.T.D. (Minimum Resolvable Temperature Difference).

In other words, it's the apparent temperature contrast required to separate the details of an image. Example: 0.1°C at 30°C .

It depends on: - detector surface ;

- electronic bandwidth (noise increases with bandwidth) ;
- specific sensor detectivity (minimum required) ;
- the amount of flux emitted

1.4.3. Thermal accuracy (trueness)

This is the absolute value of the measurement, known with varying degrees of precision. Example: 0.5°C error for 15°C , i.e. 3%. This precision is also called N.E.T.D.: Noise Equivalent Temperature Difference. Noise appears like the snow effect when a camera works with too little light. Precision must not be confused with resolution.

1.4.4. Spatial resolution

It's the ability to accurately measure the temperature of small objects, “small” being defined in relation to the total dimension of the image. It is also the number of separable points. Spatial resolution is linked to factors such as the quality of the optics, the bandwidth of the electronics, and the scanning characteristics of the detector. It defines the size of the phenomenon that can be observed on a body as a function of the recording distance. We know that:

$$L_{\text{received}} = L_{\text{object}} + L_{\text{environnement}}$$

1.4.5. Measurement conditions

The following conditions must be met:

- measurement angle no greater than 45° to the object ;
- emissivity > 0.6 ;
- object size consistent with spatial resolution;
- consideration of the external environment.

1.5. Thermal image analysis techniques

The analysis of thermal images often involves the search for singular zones whose behavior differs from the rest of the image. For this reason, a number of functions are available to artificially enhance local contrasts. Handling these functions makes it easier to find what you want to highlight in an image. The three most important functions for enhancing thermal images are thermal framing, isotherm and palettes.

1.5.1. Thermal framing

Thermal scaling is the adjustment of image scale to optimize contrast for analysis purposes. Thermal scaling requires the use of the camera's level and gain controls. Once the area of interest in the image has been selected, these two parameters need to be adjusted so that the colors in the palette cover this part of the image as closely as possible. Less interesting areas of the image may be outside the scale. In this case, they will generally be represented in black or white.

In the example shown in figure 5.1, in the right-hand image, the gain is smaller and the level significantly higher. The result is better contrast. This image is easier to analyze, as the thermal zone is better highlighted..



Fig.5.1 Same image, automatically adjusted (left) and thermally registered on the component (right)

1.5.2. Isotherm

Some thermal imaging cameras feature an Isotherm option that allows certain colors in the scale to be replaced by opposite tones. It marks an interval of equal apparent temperature. The isotherm replaces only one color range with another that contrasts more with the colors used in the image. The isotherm can be moved up and down the

scale, and enlarged or reduced as required..

1.5.3. Palettes

A palette assigns different colors to mark defined apparent temperature levels. It can be more or less contrasty, depending on the colors used.

Fig.5.2 High-voltage transformer head with isotherm

A thermal image can generally be displayed with a maximum of 256 shades of color or gray simultaneously. On a grayscale, the black color is at one end of the scale, and becomes progressively lighter with each of the 256 steps until it becomes white. This means that there will be very little contrast between, say, the 93rd and 94th shades of gray. Color images, on the other hand, use a wide range of colors to achieve better contrast. The colors must be carefully matched to each other, to give the image a semblance of smoothness, and so that it is not painful to observe. An example of a thermogram of a storage tank taken with different palettes is shown in figure 5.3..


Fig.5.3 Examples of different palettes


There's a golden rule to apply when choosing palettes:

- Use high-contrast palettes on low-contrast target objects
- Use low-contrast palettes on high-contrast target objects

2. Experimental part

Equipment required - included in "IR test bench".

- Set of 8 samples of different materials: aluminum, steel (1 raw side, 1 polished side), copper, brass, steel (1 raw side, 1 black-painted side), stainless steel, transparent PP, transparent PET.
- 1 heating resistor,
- 1 adjustable-width slot,
- 2 material supports,
- 1 angular support.

Necessary complements

Optical stands or optical bench with jumpers

Alimentation 12V 5 A

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Before the test

Mount the camera on the stand or bench and identify the various disks.

*Steel : 1 side
painted/1 side
unfinished*


Copper

Brass

Aluminium

*Stainless steel : 1
raw face / 1 polished
face*

PET Lumex clear

PP translucent

Stainless steel

2.1. Comparative study of the emissivity of different materials

Remember: *The capacity or ability of an object to absorb incident radiation is always identical to its capacity to radiate its own energy. If an object is capable of absorbing incident radiation, it will also be able to emit its own energy in the form of radiation. The reverse is of course true: an object with a low absorption capacity will necessarily have a low emissivity. An opaque object with this characteristic is actually a good reflector. Consequently, a good reflector will have a low emissivity. In infrared thermography, the most important part of the resulting radiation is generally the emitted part. An object has a certain capacity or ability to emit, characterised by its emissivity ϵ . This is denoted ϵ and varies with temperature and wavelength.*

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
A two-sided steel sample is used, with one side treated with black paint (1) and the other untreated (2). Hold the material against the hot plate. Measure the temperature of the sample on each of the 2 sides.

When heated by the same plate, the two sides of the material have the same temperature. But because the black paint changes the emissivity of the material, the two sides have different emissivities. As a result, the temperatures observed on the camera are different. Side 1 has the same emissivity as the heating plate, so in the infrared image their colours are the same. You can therefore measure different temperatures with an infrared camera, but you need to be wary of the results you obtain, as over- or under-estimating emissivity can influence the measurements.

(NB: the marks on the plate are due to traces of paint on the sample of material used).


2.2. Calibrating emissivity on the thermal camera

Reminder: *Emissivity is one of the fundamental parameters of infrared thermography. A temperature measurement by radiation on an object with an emissivity of less than 0.6 will be inaccurate.*

We use the same material as for experiment 1. Side 1 is treated with black paint, which may represent the black body, and side 2 is untreated. Heat the sample and measure the temperature on each side. Note the two temperatures and compare. Side 1 has an emissivity of 0.95, so it has the same observed temperature as the hot plate, which is its actual temperature. However, the observed temperature of side 2 is much lower than that of side 1, so its emissivity will have to be adjusted. On side 2, change the emissivity until you obtain the same temperature as side 1. To do this, simply click and drag on the area to be set directly on the image.

A window will then open, allowing you to enter an emissivity value, and then confirm..


Based on tests 1 and 2, we can conclude that it is not possible to measure an exact temperature if we do not know the emissivity of the materials observed and if we have not properly calibrated the camera.

2.3. Determining the emissivity of an unknown material

Reminder: *The emissivity of a material is the ratio of the energy it radiates to the energy that a black body would radiate at the same temperature. It is therefore a measure of a body's ability to absorb and re-emit radiated energy. Emissivity depends on several factors: the temperature, the direction of the radiation and its wavelength.*

It is difficult to measure the temperature of highly reflective materials using infrared thermography. In order to obtain the correct values and avoid the influence of reflection, the surface must be modified, for example by using Chatterton. Black Chatterton is very often used by thermographers in industrial maintenance and in the building industry. When they carry out thermograms in electrical cabinets, they place these small strips of tape on the insulation of the conductors (whose actual emissivity is unknown) so as to be close to an emissivity of 1 and thus be certain of measuring the correct temperature.


Place a rubber tape large enough to cover the measurement spot on the material. Then measure its temperature, having set the emissivity to 0.95. Using the software, measure the temperature of a neighbouring area on the object and change the emissivity until you obtain the same temperature. The exact emissivity of the measured material.

The first step is to open the image in the software and select two points, one on the rubber tape, the other on the material. Although they have the same actual temperature, you can see that their observed temperatures are different.

Select a zone on the image and set the emissivity to a lower value (0.8, for example).

Repeat the operation until you obtain a zone temperature equal to the temperature of the heating resistor. You have now estimated the emissivity of the material under study.


2.4. Influence of the material's surface condition on its emissivity

Reminder: The energy incident on an object is divided into 3 parts: reflection, absorption and transmission. We can therefore write an energy balance, assuming that the incident energy is equal to 1. At a wavelength λ , we obtain $A(\lambda) + R(\lambda) + T(\lambda) = 1$ (A : absorption, R : reflexion, T : transmission). And emissivity is equal to absorption $A(\lambda) = \varepsilon(\lambda)$


The emissivity of metals and glass also varies with temperature. Surface oxidation of molten metal or liquid/solid phase change can cause emissivity to fluctuate greatly.

Using a stainless steel disc with different surface treatments (1 rough side, 1 polished side), we can observe the influence of the material's surface condition on its emissivity.

Hold the polished and unpolished stainless steel samples against the hot plate, measure and record the temperature for each of the 2 surfaces concerned, and compare the results..

Brut 0.5mm

Polished 2mm

We can see that their observed temperatures are different. So the emissivity of these two surface conditions is not equal, even though they are made of the same material. From this example, we can conclude that the state of a surface has a significant influence on the emissivity of a material. We also note that the polished surface has a lower emissivity than the other surfaces..


To take this further, you can do the same measurements as in experiment 3. Determine their emissivities using the duct tape and compare.

2.5. Study of the influence of angle on measurement


The aim of this experiment is to highlight the importance of the angle of observation in relation to the material in question.

In this experiment, the heating plate is placed on the angled support to produce different angles of incidence from the same fixed observation point. In order to demonstrate this influence, the temperature range of the camera must be changed manually. Here, as the highest temperature is around 46°C, we set the range from 34°C to 47°C.

Change the angle of the plate from -90° to 90° and observe the colour change in the images. Note the temperature for each 10° and plot the angle-temperature curve..


Infrared image at 0°

The 80° infrared image

In the image, we can see that the temperature remains constant at the beginning, so the emissivity is constant. However, when the temperature drops to 80°, the plate emits much less. We can see that the emissivity ϵ of an object is a function of the incident angle of observation.

The curve below gives a better understanding of this phenomenon:


2.6. Study of the spatial resolution parameter using the slot


Reminder : *Spatial resolution is a measure of the fineness of detail in an image, for a given dimension. It is used to determine the size of the smallest observable detail..*

To determine the influence of the camera's spatial resolution, we use an adjustable slit and the heating plate, which represents a near-perfect black body. Place the adjustable slit approximately 30 cm from the heating plate. Place the camera 20 cm from the adjustable slit. Adjust the slit to a width of 3 mm..

At the start of the test, the slit is thin and the camera is close to the slit, so the temperature is that of the hot plate. Note the distance between the camera and the slit " d ". Then move the camera back gradually until a temperature different from that of the hot plate is observed. At this distance, widen the slit and find the temperature.

$d=20cm$


$d=160cm$

It is also particularly interesting to analyse the temperature signal using the software. Select the **Temperature Profiles** tab, and draw a line over an area that includes part of the slot. The shape of the curve can then be seen in the small graph on the left.

To save this curve, enter its name and the abscissa unit. If a length marker is available, use it when drawing the profile.

At a fairly long distance, when the slit is large, the temperature signal shows a plateau with a constant temperature, as can be seen in Figure 1 below. However, when the slit is narrow, the signal no longer shows this horizontal plateau, but just a precise peak on the temperature profile..

2.7. Special case : transparent material

Transparent PP and transparent PET samples are used for this manipulation. The slit is also used to fix them. Look at the surface of the 2 materials. One allows better visibility through than the other. Adjust the width of the slot to the maximum. Using the two magnetic holders, fix the transparent PP sample against the slit. Observe the infrared image.

For rays in the visible range, the disc is transparent and the shapes are visible through it, but without detail. In the field of view of the infrared camera, infrared rays pass through more clearly than the eye can perceive, since the surface temperature of the heating resistor is clearly perceived.

The real image


The infrared image

Replace the PP sample with PET.

To the eye, this material appears much more transparent than its predecessor. However, infrared rays do not penetrate this material and the surface temperature of the heating resistor cannot be distinguished..

The real image

The infrared image


The composition/quality of transparent materials influences their ability to transmit infrared radiation.

2.8. Demonstrating the reflection of materials, reflected temperature

Reminder : *The resulting radiation, captured by an infrared camera, is the sum of all the radiation leaving the surface of an object, whatever its source. It comes from three types of source. The resultant radiation is a combination of the radiation from the target object itself, the radiation from the reflection of the front heat source on the object and the radiation from the rear heat source passing through the target object.*

To do this, we use a polished stainless steel sample, which has very good reflective properties..

Attach the sample to the angled support and measure the temperature of the surface of the sample. A student can then stand next to the disc. Change the angle and draw a light propagation diagram to find a suitable angle so that the student's image can be seen in the camera. Observe the change in temperature

In our example, we simply placed a hand close to the sample.

The real image

The infrared image

We can see that with reflection, the temperature is no longer the ambient temperature, but that of the hand. To take this further, less reflective materials can be used to compare the resulting thermal images.

1. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

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