Électricité

Magnétisme

Electricity

Magnetism

Réf :
263 006

Français – p 1

English – p 2

Version : 3109

Aimant en U à entrefer variable
U-shaped magnet with adjustable
air gap

1 Description

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
0 825 563 563 *

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1 Description

The magnet is made up of a 30 x 3 mm steel plate folded in a U-shape over each branch on which are placed 3 magnetized ferrites of which 2 are fixed and surround a third one that is removable, as shown in the diagram below. The N branch is marked by a red disc placed on its outer face; the other branch is not marked.

- A Ferrites magnetised in the direction of their thickness glued to the branches of the U.
- B Ferrites identical to A but removable and that can be placed in position B': maintaining one position or the other is done only by magnetic attraction.

Dimensions of magnetised ferrites: 42 x 20 x 9 mm

Important note:

The magnetised ferrite with the red disc must always remain on the marked branch of the magnet (branch N), the disc glued to this ferrite must always remain visible, in position B as well as in position B'.

The ferrite with the yellow disk must always remain on the unmarked branch (branch S) and its disc must always be visible.

2 Use

When the removable ferrites are in position B, we have a uniform field over a length of 126 mm that allows us to perform, among others, Laplace's rail and wire experiments, or the magnetic spectra experiment with iron filings.

With the removable ferrites in position B', the air-gap is reduced and a very strong field is created, which allows us to perform, for example, the Barlow's wheel experiment with a current of an intensity that is much lower than that with a standard U-shaped magnet, or the magnetic brake experiment (eddy current). On the other hand, the field at the end of the branches of the U is very weak compared to that of a standard U-shaped magnet.

3 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX France
+33 (0)2 32 29 40 50

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux