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Réf :
272 019

Français – p 1

English – p 5

Version : 8109

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0 825 563 563*

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1. Introduction

The Van De Graaff Generator can be used to conduct a wide variety of electrostatic

1.1 Principe of operation

Some insulating materials when separated from the surface of others leave those surfaces electrically charged each with the opposite sign of charge and with a high potential difference (p.d.).

1.2 Belts and rollers

A flexible belt made from an insulating material and running continuously over two rollers can, by the same process, produce a supply of charge where the surfaces separate.

The two rollers have to have different surfaces and together with the belt rubber, are chosen by experiment.

2. Combs

Charges are “sprayed” on to and removed from the belt by “combs” situated adjacent to the rollers. Actual contact between the combs and the belt is not essential because of the high p.d.’s. Combs can be simply a stretched wire, or a sharp or serrated edge: action depends on very high potential gradients due to their small radii (similar action to lightning conductors).

The lower comb is maintained at or close to earth potential and is a drain for negative charge, leaving the belt with positive charges that are carried up to the top comb.

Comb adjustments

The top and bottom aluminium mesh combs can be easily adjusted by 2 screws provided on the bracket to ensure good performance.

3. Collecting sphere

The top comb is connected to a sphere which, having inherent electrical capacity (proportional to its radius) will collect and store the charge on its outer surface until discharged either by breakdown of the surrounding air as a spark, or by conduction to an adjacent earthed object.

4. The whole apparatus

The main body of the apparatus consists of a motor unit which houses a transformer and a 12 V motor a charge generator an insulating support column enclosing a belt and spun aluminium charge collecting dome which stores the charge.

The apparatus is also supplied with a discharge dome connected to a Perspex insulating handle and a 4 mm hole for earth connection.

5. Operation



The control panel is located on the front back of the motor housing.

- On the front there are 2 switches one for on/off and the other to run the motor on 2 speeds.
- On the back there is IEC socket for connecting the power cord with built in fuse.
- Additional information for optimum use of the Van De Graaff Generator.

Only run the generator for short lengths of time and switch off immediately after each experiment. This will prolong the life of the belt as this component becomes worn with continuous use. Replacement Van De Graaff Generator belts are available from your usual dealer.

Store the generator in a dust free environment or under a polythene cover.

On wet and humid days dry the air around the generator and electrostatic items of apparatus being used by means of an electric heater or hair blow direr.

6. Accessories

Certainly the De Graaff Generator can produce striking demonstrations. The usual experiments are:

Faraday's cylinder to show electric charge resides on the outer surface of a charged hollow conductor.

Bouncing ball a conducting ball suspended on a non-conducting thread touching the charging sphere will, once it itself becomes charged, be repelled away from the sphere. If then it reaches an earthed surface, it will discharge and be attached once more to be recharged and for the process to continue

The head of hair is another demonstration of repulsion. Real hair or shredded paper strips bunched at one end used and provide a sensitive means of detecting charge.

The electric wind is produced by release of ions at the end of a pointed conductor enough to deflect a candle's flame.

To demonstrate the function of a lightning conductor'.

Hamilton's Mill utilizes the electric wind at the pointed ends of four arms to cause rotation about a pivot.

Kinetic Theory Model random motion is shown by metallic balls moving by repulsion and loss of charge within a transparent vessel.

Neon Indicator shows luminous discharge from the gaseous excitation by the high electric fields near the generator.

Caution

- Disconnect power before removing cover.
- Fuse located within IEC socket.
- Do not use in the vicinity of operating computer equipment or any other electronic equipment.
-

7. Replacing belts

Caution

Only qualified or competent personnel should carry out the following procedure:

- Ensure that mains supply is disconnected when removing the cover.
- Carefully remove the pulley with shaft on the lower roller by loosening the allen key on the lower roller (the allen key is fixed in the middle of the lower roller. The top roller can easily be lifted up from the side bearings and the belt can be replaced.

8. Cleaning

Caution

Always disconnect apparatus from the mains supply before cleaning.

General dusting should be sufficient for most cleaning requirements. However occasionally the unit may require a more through cleaning and a suitable solvent for this is methylated spirit applied in conservative amounts with a clean cloth. Parts of the generator which will require cleaning are the collector dome belt and brass top comb. When cleaning these components the following procedures must be carried out carefully.

Collector dome Apply a conservative amount of methylated spirit to all parts of the dome with a clean cloth. Ensure that solvent has evaporated before connecting to the mains supply.

Top Comb- Release comb by removing the self tapping screws and apply methylated spirit with a clean cloth. Allow any residue to evaporate before replacing the comb. Position of the comb should be such that it is just touching the belt or a finite air gap exists.

Belt- Ensure that the mains supply is disconnected before removing the cover.

Remove the belt following part 7 of these instructions and noting all the safety points. The belt can then be cleaned using methylated spirit applied with a clean cloth. Allow any solvent to evaporate before replacing the belt into the unit as per part 7 levelling of the top roller will be necessary.

Note

After removing cover it is good safety practice to test apparatus for earth bonding and insulation.

9. After-sales service

The device is under a 2 years guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

+33 (0)2 32 29 40 50