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Measurement and instrumentation

**Ref :
281 003**

Français – p 1

English – p 4

Version : 8009

Alimentation Evolio F6F12 / 1 A
Evolio F6F12 / 1 A power supply

1 Description

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1 Description

- 1 – Plastic box with double insulation
- 2 – On/off switch
- 3 – Operating indicator
- 4 – 6-12 V DC / 6-12 V AC switch
- 5 – **DC** output on safety sockets
- 6 – **AC** output on safety sockets

2 Setting up

- Connect the power supply to the 230 V mains supply.
- Select the desired voltage by rotating the switch to the desired position: 6 V =, 12 V =, 12 V ~ or 6 V ~.
- Connect your circuit (lamp, motor, voltage divider, etc.) to the selected output sockets.
- Toggle the master switch so that the power-on indicator lamp comes on.

Notes:

In **DC**, if the current flow is greater than the maximum permissible value (1 A), the current remains at about 1 A and the voltage gradually drops. Remove the source of the overload or the short circuit to return to normal operation.

In **AC**, if the current flow is greater than the maximum permissible value (1 A RMS), the power supply trips. This protection is ensured by a thermal protector which resets automatically after a few seconds to a few minutes once the source of the overload or the short circuit is removed.

In **AC**, the power supply delivers a maximum voltage that can fluctuate around 12 V in function of the load of the electrical circuit and the mains voltage (between 11.9 V and 13.8 V).

Under certain operating conditions, the thermal protection (thermal switch) located in the transformer may emit some crackling noise without any danger to the user or the power supply.

3 Precautions for use

- Avoid using the power supply in direct sunlight or close to a heat source.
- Never block the ventilation vents when the power supply is operating.
- If this device is used in a manner not mentioned in these instructions, the safety ensured by the device may be compromised.

4 Technical characteristics

- **Voltage:** 6 or 12 V fixed, DC or AC
- **Outputs:** Safety sockets \varnothing 4 mm
- **Maximum current:** 1A in DC and 1A RMS in AC
- **Stabilisation:** 1 %
- **Residual ripple:** less than 10 mV
- **Protection:** Primary: Auto-reset thermal protection in the transformer.
Secondary: Electronic in DC, thermal with auto-reset in AC
- **Power supply:** 230 V \pm 10 %, 50-60 Hz
- **Power consumed from the mains:** 30 VA
- **Dimensions:** 115 x 140 x 80 mm
- **Mass:** 950 g
- **Guarantee:** 2 years

This device complies with:

- *The low Voltage directive 73/23/CEE; applied standard: EN 61010-1*
- *The CEM 89/336/CEE directive; applied standards: EN 50081-1 and EN50082-1*

4.1 Environmental conditions for use of the device

- For indoor use
- Altitude less than 2000 meters

- Operating temperature from 15°C to 40°C
- Maximum relative humidity of 80% up to 31°C and 50% up to 40°C
- Mains supply variations of +/-10%
- Installation category II, pollution level 2.

4.2 Sectioning

It is the power plug that is used as a sectioning device.

4.3 Symbols used

: Caution, refer to accompanying documents.

: Caution, hot surface.

5 Care and maintenance

In case of a malfunction:

- Check the mains voltage (230 V)
- Check the open circuit voltages, without connecting any load
- Wait 2 to 3 minutes for the protection against over-currents to be turned on again.

Note: In case of use in high ambient temperatures, the thermal protection located in the primary of the transformer may get activated in exceptional cases. In this case, disconnect the power supply and wait for 1 hour so that the device cools down.

If the device still does not work, contact the After Sales Service.

In all cases, **do not open the device**.

Store the original packaging; it can be used to bring the device to our workshops for maintenance, if necessary.

6 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
0825 563 563

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries