Instrumentation

Alimentation

Instrumentation

Power supply

Réf :
281 034

Français – p 1

English – p 2

Version : 3109

Alimentation réglable 6000 V
Adjustable power supply 6000 V

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2 Use

This power supply is designed to power electronic tubes such as the deflectron ref. 232 005 or the E/M Thomson tube.

It is recommended that the high voltage safety cords provided with this power supply be used for the 0-6000 V output and the normal safety cords (1000 V) be used for the 6 V output.

When making the connections, do not connect or disconnect when live: connect the load (the vacuum tube) with the power off, then turn on the power supply.

Refer to the instructions of the device to be powered (deflectron, E/M tube, etc.) to make the necessary connections.

Do not operate the power supply in a confined area, leave some free space above and behind the device for ventilation.

In case a fuse needs to be changed, use only those models that are approved according to the CEI 127.

If this device is used in a manner not mentioned in these instructions, the safety ensured by the device may be compromised.

In case of any problem during operation, never try to open the device but contact our after sales service.

Environmental conditions for using the device:

- For indoor use
- Altitude less than 2000 meters
- Operating temperature from 15 to 40 °C
- Mains supply variations of ± 10 % for 230 V

This device complies with:

- *The low voltage directive 2006/95/CEE, applied standard: EN 61010-1*
- *The CEM 2004/108/CEE directive, applied standards: EN 61326-1*

3 After sales service

The guarantee is for 2 years, the equipment must be returned to our workshops. For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX France

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediately to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux