**Mécanique**

***Mechanics***

**Réf :  
282 005**

**Français – p 1**

**English – p 6**

**Version : 8006**

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# 1 General


## 1.1 Purpose of the instrument

The hydroelectric power station model explains the production of electrical energy through a bicycle alternator, by transformation of mechanical energy released by the action of water falling on a turbine.

Its design is such that it demonstrates the influence of the height of falling water, on the electrical signal generated by the rotating alternator. The position energy of water can then be introduced and it can be demonstrated that the loss of position energy of the water is reflected by a gain (creation) of electrical energy.

## 1.2 Description

### 1.2.1 Composition

### 1.2.2 Operation

A bicycle alternator is fixed to the bladed water wheel (turbine). The assembly is placed in a transparent bottomless tank which is provided with an orifice in its upper part with a connection ring onto which a reservoir and two extensions are adapted.

At the back of the tank, two safety bushings connected to the alternator are used to connect a dipole (LED, light emitting diode), the power level indicator reference 282 009 (see illustration above), an oscilloscope or an interface connected to the computer to display the electrical signal generated.

## 1.3 Technical characteristics


- Anti-splash tank with transparent PVC front face and top face. Dimensions  $L \times w \times h = 250 \times 147 \times 223.5$  mm.
- Connection ring on the top.
- 300 mL PVC reservoir with a rubber plug fixed to a plastic rod.

- Two 300 mm long transparent PVC tubes with a PVC connection sleeve.
- 1 waterproof bicycle alternator.
- 1 bladed water wheel with 8 Ø 185 mm blades, fixed onto the alternator spindle.
- 1 electrical output through 2 Ø4 mm safety bushings and at a pitch of 38 mm compatible with the Eco UME dipoles on the back of the tank.

## 2 Application

### 2.1 Electricity generation

The hydroelectric power station is supported on two supports placed at the bottom of a tank in which the flowing water will enter. The reservoir, closed off and filled with water, is connected to the tank, firstly directly in case ①, then through an extension in case ② and finally through two extensions in case ③:

The power level indicator reference 282 009 is connected to the alternator output, adjusted to its maximum sensitivity (using the potentiometer) corresponding to two extensions.

The number of LEDs lit up in each situation is noted, which can lead to a conclusion about the height of the falling water.

**Note :** In order to demonstrate generation of electrical energy, a light emitting diode (LED) can be connected to the alternator output. In this case, it is impossible to compare the powers involved. The only way to enable a qualitative comparison is to use the power level indicator.


### 2.1.1 Display of voltages produced

The sequence of the above three experiments is repeated, replacing the power level indicator by an interface connected to the computer (mountain bike on which a Primo or Visio voltmeter sensor is fitted).

A synchronization is fixed in the software used (Secondary school scientific workshop) to trigger acquisition at the same voltage level in the three experiments, which simplifies the comparison of the three voltages obtained.

The following records are obtained :


Observation of variable voltages in time demonstrates that:

1. The maximum value of the voltage increases with the drop height.
2. The period  $T$  (or the pseudo-period) reduces when the drop height increases.

The analysis of a record over a duration of 2 seconds shows that :

3. The maximum value of the voltage increases when the speed increases (and vice versa).
4. The period  $T$  (or pseudo period) reduces when the speed increases (and vice versa).

It can then be deduced that the voltage is periodic and the period is constant if the water flow through the turbine is controlled (constant flow).

The falling water rotates the turbine that accelerates. The maximum value of the voltage increases.

The water flow is terminated. The turbine decelerates. The maximum value of the voltage reduces.

### 3 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your  
disposal from Monday  
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries