Mécanique

Mechanics

**Réf :
282 006**

Français – p 1

English – p 6

Version : 8006

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1 General Presentation

1.1 Purpose of the apparatus

The air and water turbine kit is intended to demonstrate the principle of electricity generation and the working of power stations. The students can very quickly construct an electricity generation system (electricity power station).


The alternator is the component common to all electricity power stations. It is rotated using a turbine, which is itself mechanically driven.

Two examples are presented in the kit:

- Action of falling water (hydro-electric power station)
- Action of an air current (wind power plant or wind turbine)

The kit, in which all the parts are separate – and not in order – is designed so the students, split into work groups, can construct the experiments and compare their results depending on the method chosen.

1.2 Composition

1.3 Technical characteristics and description

1.3.1 Air turbine

- One electric motor 1.5 V; 70 mA; 1200 rpm
- One four-blade propeller diameter 56 mm

(The motor shaft engages in an orifice in the propeller centre)

1.3.2 Water turbine

- One bicycle alternator
- One water wheel with 8 paddles, diameter 185 mm
(The mobile part of the bicycle alternator engages in the paddle wheel centre)

1.3.3 Displaying the electrical signal

- Two cables L = 500 mm with crocodile clips on each end to connect the motor or the alternator to an LED. An oscilloscope, an interface linked to the computer, or possibly the power level indicator (ref. 282 009).
- One 3 mm bare red LED

2 Instructions for use

2.1 Preliminary remark

The 1.5 V motor is used as an alternator. This DC motor is equipped with a 'commutator' which allows it to rotate, since at each half turn the voltage at its terminals is inverted. Used as an alternator, the presence of the collector means the voltage generated, variable over time, is always positive or negative, depending on the direction of rotation of the paddles.

2.1.1 Use with the power level indicator (reference 282 009)

To demonstrate the energy produced by the alternator using the power level indicator, you must check the direction of the connection. The indicator power levels are displayed using red leds equipped - like all diodes – with a 'right direction'.

2.1.2 Use with an LED

In the case where the student wishes to demonstrate the voltage generated by the alternator using an LED, it is necessary to check the polarity. In fact, when the alternator is rotated by blowing on the propeller fixed to its shaft, it generates variable but non alternating voltage of maximum value below the LED voltage threshold. To light up the LED, a current caused for example by releasing compressed gas (air, oxygen etc) is necessary.

2.2 Hydro-electric power station

The turbine is linked to the bicycle alternator attached to a tightening yoke on a base.

A LED is connected to the alternator terminals. The power level indicator is also connected.

A jet of water rotates the paddle wheel.

Observe that the LED flashes (alternating voltage) and the power level can be marked on the level indicator (maximum on the illustration).

2.3 Wind power plant (wind turbine)

The propellor is linked to the 'motor', which becomes an alternator, connected to the power level indicator (Caution: watch the polarity: see 2.1).

An LED may also be connected to the alternator terminals (Caution polarity: see 2.1).

Then blow in front of the propeller blades to rotate the motor.

Observe the power level of the indicator is below that achieved with the hydro-electric power station (five LEDs lit up on the illustration).

If an LED is connected in the correct direction, it can also be observed it does not light up.

It can be deduced the alternator generates electricity of lower voltage than obtained with the hydro-electric power station.


2.4 Display of the voltages generated by the alternators using an interface

To display the voltage generated by a rotating alternator, connect it to an input on the selected interface (Visio, Primo or VTT equipped with a volt meter adaptor).

2.4.1 Alternator linked to the water turbine

Rotate the alternator by pushing a turbine blade with a finger.

Capture of voltage gives the following results:


Observation of the alternating voltage over time demonstrates that:

1. The maximum value of the voltage decreases as the speed of rotation decreases.
2. The period (or pseudo-period) increases when the speed of rotation decreases.

2.4.2 Alternator linked to the air propeller

Rotate the alternator by blowing briefly in front of the propeller.

Recording of the voltage gives the following results:

Observation of the alternating voltage over time demonstrates that :

1. The maximum value of the voltage increases as the speed of rotation increases then decreases as it reduces.
2. The voltage is always of the same sign (here, positive). This is because of the presence of a “commutator” (see 2.1).

2.4.3 Comparison of two sources of voltage

With the bicycle alternator linked to the water turbine, impulsion allows generating voltage of which the maximum value exceeds 2 volts, sufficient to light up an LED.

With the ‘small’ alternator connected to the propeller, ‘energetic’ blowing generates voltage of which the maximum value is in the order of 1.5 V, insufficient to light up an LED.

The results can be compared with the reality (hydro-electric power station and wind power plant) and the technologies used discussed, pointing out the dimensions of air and water power stations. In particular, for wind power plants, the need for very large propeller blades and a major gear reduction system in the propeller/alternator link can be demonstrated.

3 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your
disposal from Monday
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries