TP de chimie

Chemistry tutorial


**Ref :
282 023**

Français – p 1

English – p 3

Version : 7001

**Pile Daniell
*Daniell cell***

1 Principe - description

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
1 Principle – description

The phenomenon of the polarisation of electrodes creates a dissymmetry between these and thereby transforms an electrolyser into a generator.

1.1 Description

The DANIELL cell is made up of a zinc electrode dipped into a saturated solution of zinc sulphate ($\text{Zn}^{2+}, \text{SO}_4^{2-}$), a copper electrode dipped in a saturated solution of copper sulphate ($\text{Cu}^{2+}, \text{SO}_4^{2-}$) and porous vase separating the two electrolyte solutions (see the drawing).

1.2 Composition

The apparatus delivered dry and without the salts necessary for the experiments, includes:

- 1 cylindrical beaker of 600 mL,
- 1 porous vase,
- 1 electrode in zinc,
- 1 tubular electrode in copper.

Other accessories (not delivered):

- Copper sulphate crystals or powder,
- Zinc sulphate crystals.

1.3 Procedure

Place the 100 g of copper sulphate in the porous vase, add 200 mL of water, agitate with the copper electrode and place it in the centre. (The copper sulphate, undissolved, holds it vertically).

Prepare the zinc sulphate solution by dissolving 25 g of the crystallised product in 500 mL of water and place it into the glass container.

Place the zinc electrode concentrically to the porous vase.

Unipolar cords with banana clips shall be used for connecting the equipment to the usage circuit.

WARNING:

For preparing the copper sulphate solution, you can use:

- Either, copper sulphate crystals (blue): proceed as explained in the instructions by placing 100 g of CuSO_4 , 5 H_2O in the porous vase with 200 mL of water.
- Or, anhydrous copper sulphate (white powder): it is then necessary to use another recipient for preparing the solution, by placing the powder in the water. If the porous vase is used, there is a risk of destroying it during the crystallisation of copper sulphate in contact with the water.

2 Application

The DANIELL cell cannot be polarised. It has therefore a constant e.m.f. of 1.08 V and can be used in all the experiments in which we need a constant source with a very low output.

Note:

When the cell is working, the ions Cu^{2+} and SO_4^{2-} from the copper sulphate, each migrate to a pole, the concentration thus decreases little by little. It is therefore necessary to have in the porous vase undissolved surplus copper sulphate that enters into the solution gradually.

On the other hand, The SO_4^{2-} ions attack the zinc by producing zinc sulphate: the cathode thus dissolves slowly.

Because these cells are used most often in the circuits where there is very little output, their life is very long and in general, the question of replacing the cathode is not posed before several years of service.

3 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX France
+33 (0)2 32 29 40 50

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Service gratuit * :
+ 33 (0)2 32 29 40 50

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Rubrique FAQ

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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