Electricité

Electricity

**Ref :
283 040**

Français – p 1

English – p 3

Version : 6006

Pinces crocodiles sur support

Alligator clips on a stand

1 Principe et description

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1 Principle and description

The “alligator clips on a stand” product was designed to help students perform experiments to study conductors and insulators, as well as short circuits. The device consists of 2 alligator clips mounted on a plastic stand; each clip is connected to a 4mm diameter banana safety socket located on the side. A metal plate located on top of the box protects it from any incandescent residue from the steel wool that might fall on it during the short circuit study.

2 Technical specifications

$U_{\max} = 24 \text{ V}$

$I_{\max} = 3 \text{ A}$

Dimensions of the box: 105 x 65 x 33 mm

3 Experiments

3.1 Studying conductors and insulators

Other materials needed: various everyday objects, one 4.5 V battery, 3 cords and 2 other alligator clips, a lamp (or an LED or a motor or a buzzer).

Principle & description: The student shall test various materials in order to determine if they are conductors or insulators.

In order to do this, connect the "alligator clips on a stand" device and a lamp in series to the terminals of a battery; and then test various materials by attaching them between the alligator clips of the device.

3.2 Studying short circuits

Other materials needed: one 4.5 V battery, 2 cords and 2 other alligator clips.

Principle & description: The purpose of this experiment is to illustrate the hazards resulting from a short circuit in a generator. To perform this experiment, attach a small piece of steel wool between the alligator clips of the "alligator clips on a stand" device. With the aid of 2 other alligator clips and 2 cords, attach the assembly to the terminals of a battery. After a few seconds, the steel wool becomes incandescent.

4 After sales service

This device is under a 2 year warranty and should be returned to our factory in the event of any defects.

For all repairs, adjustments, or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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Service gratuit * :
+ 33 (0)2 32 29 40 50

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Rubrique FAQ

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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