# **Laboratoire**

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## ***Laboratory***

*Multi-purpose electrode support*

Réf :  
283 153

Français – p 1

English – p 15

Version : 8007

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***Multi-purpose electrode support***

## 1 Description

### 1.1 Généralités

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# 1 Description

## 1.1 General

The multi-purpose electrode support is a device that fixes and maintains a wide range of electrodes to perform electrochemical manipulations (electrolysis, batteries, etc.) and pH-metric and conductimetric analyses. It is compatible with a large number of equipment types and is designed to work under various configurations. It can be used to create a large number of setups that were formerly difficult in practice.

## 1.2 Composition


Diagram 1

- 1 – Multi-electrode housing
- 2 – Spiral for thermometer
- 3 – Housing for Ø 12mm probes
- 4 – Passage for burette tip
- 5 – Clamping screw

- 6 – Housing for stand
- 7 – Cable clamp
- 8 – Electrode terminals
- 9 – Sliding attachment stops

## 2 Complementary equipment

### 2.1 Electrodes provided

The multi-purpose electrode support is delivered with two  $\varnothing 6 \times 100$  mm graphite electrodes.

These electrodes may be ordered separately:

$\varnothing 6 \times 100$  mm graphite electrodes – set of 2

Ref. 283 088

### 2.2 Other electrodes

The multi-purpose electrode support may be used with metallic electrodes to make batteries and perform electrolysis.

It will hold two electrode sizes:

#### 120 × 30 mm large electrodes

<i>Description</i>	<i>Reference</i>
Copper electrode - thickness 2 mm	107 033
Iron electrode - thickness 2 mm	107 051
Zinc electrode - thickness 1.5 mm	107 419
Aluminium electrode - thickness 2 mm	107 006

#### 100 × 10 mm small electrodes

<i>Description</i>	<i>Reference</i>
Copper electrodes - thickness 1 mm - set of 10	107 204
Iron electrodes - thickness 1 mm - set of 10	107 203
Zinc electrodes - thickness 1 mm - set of 10	107 205
Aluminium electrodes - thickness 1 mm - set of 10	107 206

## 3 Usage precautions

### 3.1 Warning

The JEULIN Company cannot be held responsible for an accident that occurs due to failure to respect the safety instructions described in this manual.

#### 3.1.1 Usage conditions

This instrument is designed for experiments carried out cold. Use close to heated products or heating devices in operation can irretrievably damage the instrument.

When it is used as an electrode support for electrolysis, the multi-purpose electrode support must be connected to a power supply or to a battery outputting a maximum voltage of 30 V.

**It is strictly prohibited to connect the electrode support to the electrical power supply network. Failure to respect this instruction can cause serious accidents.**

### 3.1.2 Risk of fire

Electrolysis of aqueous solutions can produce dangerous gases. Hydrogen and oxygen can be emitted, that can cause fire and explosions.

- Apply all precautionary measures to prevent accumulation of these gases by ventilating the room in which these experiments are carried out.
- Do not leave the instrument in operation without supervision.
- Do not work close to inflammable or oxidizing substances.

### 3.1.3 Protection of persons

- Do not undertake any experiment without firstly becoming familiar with the risks that arise with the chemicals used.
- Wear personal protection equipment appropriate for the chemicals used.
- Toxic gases (chlorine, etc.) may be produced during certain electrolysis manipulations. Make sure that these gases cannot accumulate by ventilating the room in which the experiments are carried out and do not prolong experiments unnecessarily.
- Make sure that aisles and accesses to work stations are not cluttered to allow free circulation in case of emergency.

## 4 Setup and installation

### 4.1 Assembly

The multi-purpose electrode support is delivered disassembled. Screw in the terminals for the electrodes (*Diagram 1-8*) and the clamping screw (*Diagram 1-5*) into their housing as shown in the following diagram.


Diagram 2


### 4.2 Placement

#### 4.2.1 Placement on a beaker

The multi-purpose electrode support is designed to be installed on a 100, 250 or 400 mL beaker. The sliding attachment stops guarantee stable and safe assembly.


Proceed as follows to make the assembly (*diagram 3*):

- 1 Place the multi-purpose electrode support on the beaker.
- 2 Slide the sliding attachment stops until they are in contact on the neck of the beaker.

*Diagram 3*

The support is then firmly fixed to the beaker. The position of the stops on the beaker neck can be seen in the following diagram.


*Diagram 4*

#### **4.2.2 Placement on a stand**


The multi-purpose electrode support may be installed on a stand. Proceed as follows (*Diagram 5*):

- ① Insert the support rod into the stand housing (*Diagram 1-6*) and adjust it to the required height.
- ② Fix the multi-purpose electrode support onto the stand, by tightening the clamping screw (*Diagram 1-5*).

*Diagram 5*

- ③ Adjust the position of the universal support above the support base by turning it around the centre-line of the stand (*Diagram 6*). Then fix the multi-purpose electrode support in place by tightening the clamping screw.

*Diagram 6*

## 4.3 Installing electrodes and probes

### 4.3.1 Placement of electrodes for electrolysis and batteries

The different electrodes for electrolyses and batteries mentioned in chapter 2 are positioned in the same way, the multi-electrode housing (*diagram 1-1*) being compatible with all electrode formats.

Proceed as follows:

- ① Slide the electrode in the multi-electrode housing.
- ② Fix the electrode by tightening the terminal for electrodes.

#### Small electrodes


Diagram 7

#### Large electrodes


Diagram 8

Ø 6 mm graphite electrodes


Diagram 9

#### 4.3.2 Placement of pH and conductivity probes

The universal electrode support can also hold Ø 12 mm pH or conductivity probes. Furthermore, the housing for Ø 12 mm probes (*Diagram 1-3*) has two positions, consequently up to two probes can be put into place simultaneously.

Proceed as follows to place a probe on the multi-electrode support (*Diagram 10*):

- ① Slide the probe into the support, allowing it to pass through the central part of the housing for Ø 12 mm probes. Adjust the probe to the required height.
- ② Fix the probe by « clipping » it into one of the two side compartments. The second probe is put into position in the same way in the second compartment.


Diagram 10

#### 4.3.3 Placement of a temperature probe

The universal electrode support can also be used to hold thermometers and temperature probes. The thermometer spiral can hold all thermometers with a diameter of between 1 and 6 mm by blocking (Diagram 1-2)

The attachment is made as follows:

- ① Slide the thermometer or the temperature probe into the thermometer spiral, at the wide end of the spiral.
- ② Adjust the thermometer or the temperature probe to the required height.
- ③ Fix the thermometer or the temperature probe in place in the spiral, pushing it sideways towards the narrow part of the spiral.


Diagram 11


#### 4.3.4 Fixing cables

Probe cables and wires used to create a circuit during electrochemical experiments often cause hindrance and congestion during manipulations. The universal electrode support can solve this problem, because they can be fixed to the back of the device. Two cable clamps located on each side and at the back of the universal electrode support are used to fix the cables (*Diagram 1-3*)


Proceed as shown in diagram 12, to fix the cables on the multi-purpose electrode support clamp.

Push the cable into the cable-clamp. It will get jammed in the notch corresponding to its diameter.

To remove the cable, simply pull the cable to extract it from the cable-clamp.

*Diagram 12*

*Diagram 13*  
*View of cable in position in the cable-clamp.*

## 5 Example of experimental applications

### 5.1 Electrochemical applications

#### 5.1.1 Battery mode

The multi-purpose electrode support can be used to make battery elements using two strips of different metals (for example Zinc and Copper) immersed in a conducting solution contained in a beaker.

The current produced can be demonstrated by connecting this battery to a LED (*Diagram 14*). Several elements can also be assembled using different multi-purpose electrode supports and they can be installed in series.


Diagram 14

#### 5.1.2 Electrolysis mode

The multi-purpose electrode support can be used as an electrolyser to perform electrolysis manipulations. Simply connect it to an electrical power supply using cables with a  $\varnothing$  4 mm safety banana pin. A voltmeter and an ammeter can be connected to study the current in the circuit and the voltage at the terminals of the electrolyser. It can be used with graphite or metallic electrodes depending on the manipulation performed.


Diagram 15

The power supply voltage shall never exceed 30 V.

### 5.1.3 Conductivity tester mode

The multi-purpose electrode support may also be used as a conductivity tester for liquids, solutions and powders. To achieve this, it is connected to a power supply and is mounted in series with a dipole to characterize the presence of current in the circuit (Lamp, LED, buzzer, etc.). Graphite electrodes should be used in this case.


Diagram 16


The power supply voltage shall never exceed 30 V.

## 5.3 Applications in analytic chemistry

The other main function of the multi-purpose electrode support is to hold pH or conductivity electrodes during analyses.

The multi-purpose electrode support was designed to facilitate the ergonomics of setups using a beaker, a Mohr burette, a magnetic stirrer and a pH or conductivity electrode.

*Diagram 17*

The multi-purpose electrode support is useful for positioning and fixing all elements and it makes manipulations easy and safe (*Diagram 17*). One of the main innovations of this support is that it provides an opening to enable a burette tip to pass through (*Diagram 1-4*). All elements required to make an analysis can thus cohabit in the same volume overhanging the beaker.

## 6 Care and storage

### 6.1 Care

Although it is composed of strong materials (polypropylene, stainless steel, etc.), it is strongly recommended that the multi-purpose electrode support should be cleaned after each manipulation, and particularly after having been used for electrolysis manipulations releasing corrosive gaseous compounds (chlorine, etc.) in order to prevent any risk of degradation and oxidation. It may also be necessary to clean it if it is brought into contact with aggressive solutions.

Cleaning is done by rinsing the multi-purpose electrode support with running water. Then dry it with a cloth or absorbent paper.

Never use an organic solvent to clean the multi-purpose electrode support, this could irretrievably damage the parts made of a plastic material.

## **6.2 Storage**

There are no particular storage requirements.

## **7 After-Sales Service**

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your  
disposal from Monday  
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries