**Electricité**

Magnétisme

***Electricity***

*Magnetism*

Réf :  
292 002

Français – p 1

English – p 2

Version : 1109

**Aiguille d'Oerstedt**

***Oerstedt's needle***

## 1 Principe – description

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
**0 825 563 563 \***

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## 1 Principal – description

This experiment has the purpose of showing the magnetic effects of current.

- The magnetised needle and the conductor thread are placed on a transparent baseplate enabling the projection of the experiment. The back side of the circuit can be easily removed in order to be used with an overhead projector et to allow a display for all the students
- The 12/10copper thread, has sockets for banana plugs  $\varnothing$  4 mm I max : 10 A.
- Magnetised needle length 65 mm.

## 2 Handlings

Necessary accessories:

- Oersted's needle,
- Low voltage power supply with rheostat (or variable power supply),
- Ammeter,
- Current inverter.

In order to put the needle on its pivot, under the conductor thread, orientate it so the wire and needle are parallel.

Supply the thread  $I < 10$  A : the needle deviates and has the tendency to displace itself normally to the thread when the voltage is increased. Several amperes is enough to obtain a clear deviation.

**N.B.:** The transparent plastic material is electrifiable, which disrupts the magnetised needle.

**BEFORE EACH MANIPULATION, DE-ELECTRIFY THE SUPPORT  
WITH A CLOTH MOIST WITH ALCOHOL**

## 3 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27019 EVREUX CEDEX France**

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your  
disposal from Monday  
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries