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Electricity
Electromagnetism

Réf :
292 007

Français – p 1

English – p 4

Version : 9109

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0 825 563 563*

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With this apparatus, it is possible to show evidence of magnetic induction actions on a current.

The large diameter bobbin (~ 50 mm) is made up of 200 copper thread spires glazed 5/10th mm.

It is suspended by a two coaxial conductor thread with great flexibility so as not to disrupt the mobility of this coil.

The electrical link operates by a bipolar female socket fixed to the end of a steel tube which furthermore allows the apparatus to be suspended from a support using a tightening nut.

This device is equipped with 2 adapters plugged into 4 mm secure sockets already mounted for safe connection to the low voltage supply. If you are not disconnected from the housing, a six-pencil key of 1.5 mm is needed to tighten it (not provided).

1. Principal – description

When a magnet is approached by the flat coil suspended by two thread conductors which supply the coil with current, it inserts itself on the magnet and is placed approximately in the middle, so that the flux enters by the south face maximum.

1.1 Set up

Directly connect the coil to a continuous 6 V source (accumulator or inverter) ; the low heating enables use which is sufficiently long to show a demonstration.

1.2 Useful accessories to carry out a manipulation

Ref. 263 007 : Ticonal straight magnet 12 x 70 mm,

Ref. 281 483 : F6F12 Evo2 power supply,

Ref. 701 030 : Plate with rod,

Ref. 703 099 : Tightening nut.

2. Manipulations

2.1 Showing evidence of the max flux rule


Classic experiment. The excited coil is polarised; following the direction of the passing current, it is repulsed or attracted by the magnet.

If a straight bar magnet is placed following the coil axis, it can be seen to run into it.

When the current direction is inversed, firstly the coil is violently repulsed, returns face to face to run in once again.

In the previous experiment, the bar magnet can be replaced by a 2nd excited coil. Following the direction in which the current passes first you can see the 2 coils placed side by side facing each other or repelled side by side when the 2 opposite faces being of different polarities.

This experiment therefore needs 2 identical coils.

2.2 Mutual action of 2 coupled coils

This very interesting experiment also needs 2 coils.

On 2 distinct rod supports, horizontally fix 2 straight bar magnets so that they pass through the centre of the suspended coils above them. Electrically connect the 2 coils. Oscillate one of the two coils. After several oscillations you can see the other coil take up a synchronous pendular movement.

3. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.
For any repairs, adjustments or spare parts please contact:

JEULIN – TECHNICAL SUPPORT
468 rue Jacques Monod
CS 21900
27019 EVREUX CEDEX FRANCE

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.

** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux