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Electricity
Electromagnetism

Réf :
292 014

Français – p 1

English – p 5

Version : 4107

Bobine plate – bobine de Helmholtz

Flat coil – Helmholtz coil

1. Bobine plate

1.1 Description (figure 1)

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1 Flat coil

1.1 Description (figure 1)

The flat coil has a mean radius $R = 6.5 \text{ cm}$ and a width $2L = 2.5 \text{ cm}$.

The coil has $N = 95$ turns in 5 layers of 19 turns of 1 mm^2 enameled copper wire.

The maximum admissible current is 5 A (7 A for a few minutes).


Fig. 1

1.2 Justification of coil dimensions

The coil dimensions are such that it can be assimilated to a theoretical "infinitely flat" coil.

❖ Calculations for a theoretical coil (figure 2)

a) Induction at coil center O:

$$B_0 = \mu_0 \frac{NI}{2R} = \mu_0 I \frac{95}{0,13} = \mu_0 I \times 731$$

b) Induction at point M on the axis OM = R

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$$B_0 = \mu_0 I \times 3,9 \cdot 10^3 \times 0,19 = \mu_0 I \times 740$$

b) Induction at point M on the axis OM = R (figure 4)

$$\theta_1 = 40^\circ ; \cos \theta_1 = 0,766$$

$$\theta_2 = 129^\circ ; \cos \theta_2 = -0,629$$

$$B = \mu_0 I \times 3,9 \cdot 10^3 \times \frac{1}{2} (0,766 - 0,629) = \mu_0 I \times 269$$


Fig. 3


Fig. 4

❖ Conclusion:

The difference between the values obtained with the coil P/N 292.014 and the theoretical coil is less than 2%, i.e. the order of magnitude of the measurement accuracy of the Teslameter.

2. Helmholtz coils

2.1 Description (figure 5)

The apparatus is comprised of two coaxial coils B1 and B2 identical to the coil described in Section 1. Coil B1 is stationary. Coil B2 is movable.



The pointers I1 and I2 corresponding to the mid-point of the coils, and the graduation G engraved on the apparatus stand are used to measure the distance D between the centers of the two coils.

The guide E carries the Teslameter probe holder M. The probe S is then positioned in the coil axis.

The graduation on holder M allows the direct reading of the distance d between the probe S and the centre of the stationary coil B1.

Each coil has two associated terminals marked E and S for energizing the coils.

The connections shown on figure 6 below give additive induction values.

3. After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
 468 Rue Jacques Monod
 CS 21900
 27019 EVREUX CEDEX France

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Direct connection for technical support

A team of experts
at your disposal
from Monday to Friday
(opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request
immediatly to provide you
with the right answers regarding
your activity field : Biology, Physics,
Chemistry, Technology.

Free service*

+33 2 32 29 40 50**

* Call cost not included.
** Only for call from foreign countries.

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SAS au capital de 1 000 000 € - TVA intracommunautaire FR47 344 652 490 - Siren 344 652 490 RCS Evreux