# **Electricité**

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## ***Electricity***

*Electromagnetism*

**Ref :  
292 043**

**Français – p 1**

**English – p 5**

**Version : 9003**

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## ***Transformer and coils***

## 1 But

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## NOTES

## 1 Purpose


- To show there is no connection between the coils of a transformer
- Importance of the core and "magnetic" armature
- Influence of the number of windings
- Reversibility of the transformer

## 2 Equipment required

- |                                    |              |
|------------------------------------|--------------|
| • 6 V DC and AC supply             | Ref: 281 083 |
| • 6 V lamp on support              | Ref: 283 044 |
| • Set of 25 bulbs 3.5 V 0.2 A      | Ref: 283 093 |
| • Bell push type switch (set of 6) | Ref: 283 164 |
| • BT direct current motor          | Ref: 282 052 |
| • Multimeter CL 2035               |              |

## 3 First use

Assemble as follows:

Without a core, the coil has low resistance, and the current is then large: the coil heats and the protective varnish can be destroyed if the current is maintained (particularly at 12 V). To avoid this problem it is recommended that you put a push type switch in the primary circuit.

## 4 Experiments

### 4.1 The role of the magnetic circuit

- without the core: the bulb does not light:  $U_{CF} = 0.03 \text{ V}$
- place the core in the coil: the bulb lights slightly
- Close the magnetic circuit (fix U to I, core in coil): the bulb lights well.

### 4.2 Role of the number of windings

Read the voltages between input  $U_E$  and output without load (without bulb), then under load (with bulb): complete:

$U_E$	Voltage without load			Voltage with load		
	$U_{CF}$	$U_{DF}$	$U_{EF}$	$U'_{CF}$	$U'_{DF}$	$U'_{EF}$
6						
12						

Make the ratios  $U_E/U_{CF}$ ,  $U_E/U_{DF}$ ,  $U_E/U_{EF}$  and compare with the ratio to the number of **windings**  $m$  (1, 1/2, 1/4).


### 4.3 Observing signal form on the oscillograph

Check that:

- the voltage has the same form
- the input and output frequencies are the same
- the voltage value is reduced in a ratio equal to  $m$ .

### 4.4 Reversibility

- Now supply the second coil from the 6 V AC generator

- Measure the voltage  $U'_S$  at the terminals of the second coil: you will notice an increase in this voltage, with a ratio of 2.
- Also measure with the AC voltage  $U_E=12$  supplied by the generator

## 4.5 Principle of the alternator

Effect of a magnet moved in front of a coil.

Refer to the instructions supplied with the motor with magnet. Ref: 282 052

## 5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

**JEULIN - TECHNICAL SUPPORT**  
**Rue Jacques Monod**  
**BP 1900**  
**27 019 EVREUX CEDEX FRANCE**  
**+33 (0)2 32 29 40 50**

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

## Direct connection for technical support

A team of experts at your  
disposal from Monday  
to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

### Free service \*

**+ 33 (0)2 32 29 40 50\*\***

\* Call cost not included

\*\* Only for call from foreign countries