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Electromagnetism

**Réf :
292 077**

Français – p 1

English – p 10

Version : 8006

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3.5 Vérification expérimentale

The graph shows voltage U (V) on the y-axis (ranging from -2 to 1) versus time t (ms) on the x-axis (ranging from 0 to 300). The signal starts at 0 V, remains flat until approximately 140 ms (point 1), then drops to a minimum of about -1.2 V at 180 ms (point 3), rises to a peak of about 0.1 V at 200 ms (point 5), drops to a minimum of about -2.0 V at 240 ms (point 7), and rises to a peak of about 0.1 V at 260 ms (point 8). The signal then drops to a minimum of about -2.5 V at 280 ms (point 9) and rises to a peak of about 0.1 V at 300 ms (point 10).

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
1 General

1.1 Purpose of the apparatus

The alternator demonstration scale model is designed to illustrate during experiments, the generation of electricity through the conversion of mechanical energy. It can also be used to study the influence of a magnet's polarity on the direction of the current induced in qualitative fashion, by switching on a LED.

1.2 Description

①	Reversible straight magnet
②	Armatures in soft iron
③	Coil
④	Collector formed by two copper half-rings, forming part of the coil, each linked to one end of the wire constituting the coil.
⑤	Safety sockets linked to two copper blades which rub on the collector
⑥	Handle linked to a pulley to rotate the coil

1.2.1 Principle

In the demonstration alternator, the coil is rotated using the handle. Transmission is by a belt. A straight magnet is placed on two soft iron armatures to create a practically uniform magnetic field in which the coil can rotate. The magnet is mobile to demonstrate the influence of polarity on the direction of current.


The ends of the wire constituting the coil are linked to a collector formed by two copper half-rings which themselves form part of the coil. Two copper blades linked to two safety sockets, rub on the collector.

It is possible to connect the following to the alternator: a power level indicator (reference 282 009), a dipole (for example a LED), an oscilloscope or an interface (Visio, Primo or VTT equipped with a voltmeter adaptor) linked to a computer.

1.2.2 Voltage at the safety socket terminals

If the voltage produced at the coil terminals as it rotates in the magnetic field is alternating, the voltage obtained at the safety socket terminals is rectified "double alternating".

1.3 Technical Characteristics

The coil, equipped with a soft iron hub, is mobile around a stud which on one side forms part of the pulley, and on the other side, of a collector formed by two copper half-rings.

The stud is fixed to a plastic support screwed into a base which is also plastic.

Two semi-circular soft iron armatures are screwed to the support. Their upper form allows receiving a straight magnet constituted of two ferrites of which the polarities are marked by a colour code.

A handle is used to rotate the coil via a belt. Demultiplication is 10.

2 Instructions for use

2.1 Generation of electrical energy by conversion of mechanical energy

Connect an incandescent lamp (3.5 V; 0.2 A for example) to the alternator terminals. As soon as the coil is rotated (that is mechanical energy is supplied to the alternator), it can be observed the lamp lights up (irrespective of the direction of rotation of the coil) which means it is receiving electrical energy.


It can then be concluded that the mechanical energy received by the alternator is converted into electrical energy.

2.2 Study of electrical power

Connect the power level indicator (reference 282 009) to the alternator terminals.

Observe the power increases as the speed of rotation increases.

The electrical power available increases depending on the voltage generated at the alternator terminal. By increasing the speed of rotation, the voltage is

increased (see paragraph 2.4) and as a result, also the available electrical power.

Remark : The LEDs of the power level indicator can only light up in the passing direction, that is, only if the current which passes through them is circulating in the “right direction”. For this reason, the power level indicator can be used for only one direction of rotation of the coil, which depends on the polarity of the magnet, and the direction in which the power level indicator is connected.

2.3 Demonstration of the effect of the direction of rotation and polarity of the magnet

2.3.1 Effect of the direction of rotation

Rotate the alternator connected to a LED, first in one direction and then in the other. It can be observed only one direction of rotation lights up the LED.


DEL

2.3.2 Effect of polarity of the magnet

The LED lights up when the direction of rotation of the coil is contrary to the previous one.

DEL


2.3.3 Explanation

Mark one of the sides of the coil (noted F in the diagrams). With the magnet in a certain position, rotate the coil in one direction then another.


The polarity of the magnet is then inverted and the coil turned in one direction then another.

The situations in which the LED does and does not lights up are summarised below.

2.4 Voltage at the alternator terminals depending on rotation

Rotate the alternator off-load, connected to an interface linked to the computer. The speed of rotation increases and recording of the voltage gives:

Observation of the voltage which is variable over time allows highlighting two important phenomena:

When the speed of rotation increases :

3. The maximum voltage value increases
4. The period (or pseudo-period) reduces

3 Theoretical study of the alternator

3.1 Magnetic flux

Let S_T be total of the area of the turns forming the coil,

The expression of flux is : $\Phi = S_T \vec{B} \cdot \vec{n}$.

3.2 Arbitrary choice of the normal to the coil

Select as the positive direction, the direction of winding of the turns. The rule of the right hand gives the direction of the normal (see diagrams).

3.3 Electromotive force (EMF) : e

Any variation in flux produces an electromotive force of which the expression



is: $e = -\frac{d\Phi}{dt}$.







3.4 Study of electromotive force voltage and voltage at the alternator terminals over time

Consider the coil is rotating in relation to the magnet:

In the following study, reference is made to the surface of the coil from which the normal emerges, then $e = U_{BA}$.

Mark several positions of the coil in rotation in relation to the magnet. The explanations refer to the magnetic flux and its variations.

①	②
	
<p>Magnetic flux. The flux, negative, is minimal $e = U_{BA} = 0 \Rightarrow U = U_{BB} = 0$</p>	<p>Magnetic flux. The flux is negative and increasing $e = U_{BA} < 0 \Rightarrow U = U_{BA} < 0$</p>


	
<p>Magnetic flux. The flux is nil and increasing. $e = U_{BA} < 0 \Rightarrow U = U_{BA} < 0$</p>	<p>Magnetic flux. The flux is positive and increasing. $e = U_{BA} < 0 \Rightarrow U = U_{BA} < 0$</p>
	
<p>Magnetic flux. The flux, positive, is maximum. $e = U_{BA} = 0 \Rightarrow U = U_{AA} = 0$</p>	<p>Magnetic flux. The flux is positive and decreasing. $e = U_{BA} > 0 \Rightarrow U = U_{AB} < 0$</p>
	
<p>Magnetic flux. The flow is nil and decreasing. $e = U_{BA} > 0 \Rightarrow U = U_{AB} < 0$</p>	<p>Magnetic flux. The flux is negative and decreasing. $e = U_{BA} > 0 \Rightarrow U = U_{AB} < 0$</p>

The variations in the electromotive force can be summarised in the following tables:

positions	①	②	③	④	⑤	⑥	⑦	⑧
Flux Φ	mini	\nearrow	0	\nearrow	maxi	\searrow	0	\searrow
$e = -\frac{d\Phi}{dt}$	= 0	< 0	< 0	< 0	= 0	> 0	> 0	> 0
		\searrow	mini	\nearrow		\nearrow	maxi	\searrow
u	= 0	< 0	< 0	< 0	= 0	< 0	< 0	< 0
		\searrow	mini	\nearrow		\searrow	mini	\nearrow

3.5 Experimental verification

The experimental verification gives :

Remark: the amplitude of the voltage and of the period varies on account of an increasing speed of rotation

4 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact :

JEULIN - TECHNICAL SUPPORT
 Rue Jacques Monod
 BP 1900
 27 019 EVREUX CEDEX FRANCE
 +(33) 0825 563 563

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Service gratuit * :
+ 33 (0)2 32 29 40 50

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Rubrique FAQ

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

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