# Mesure

Initio® 2

# Measurement

Initio® 2

Réf :  
292 087

Français – p 1

English – p 7

Version : 1205

**Joulemètre Initio® 2**

**Initio® 2 Joulemeter**

## 1. Description

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**0 825 563 563\***

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## 1. Description

The Initio® 2 joulemeter is an instrument which measures voltage, intensity (direct or effective), power consumed, energy and time. It also provides information on the phase shift between voltage and current.


It can be used to record data, store it and then restore it later on a computer.

A WIFI module (ref 451 470) is also available for visualising and saving measurements on a tablet or a PC equipped with WIFI.

Equipped with an analogue output, the Initio®2 joulemeter can be used with an ExAO interface with direct inputs via banana jacks or via a Groomy® interface and jack.

## 2. Use

Upon first use, we recommend that the battery of the Initio2 be fully charged. For charging use the  $\mu$ USB (7) charger supplied in separate packaging. In general, we recommend using the instrument in complete charge-discharge cycles (wait until it is completely discharged before fully recharging).

## 2.1 Putting into use

- Press the start button (1) situated on the left side of the instrument.

The measurement displays directly on the instrument screen.

- To switch off the Initio®2, briefly press the start button (1), turn the knurled knob to make your selection and confirm by pressing the knurled knob (3).

## 2.2 Navigation

Initio®2 has three navigation buttons:

- Selection-confirmation (3):

The central knurled knob (3) is for moving between the menus of the Initio®2. Briefly pressing the OK button in the center of the knurled knob confirms the selection.

- Back (2):

The "Back" button (2) makes it possible to change from one type of display to another from the main screen in the following order:

- Measurement
- Table
- Chart

Within a menu it is used to go back a step.

- Menu (4):

This button displays the main menu of the Initio®2.

## 2.3 Changing measurement type

The joulemeter measures various physical values

- Intensity
- Voltage
- Power consumed
- Energy consumed

### Changing value

To switch from one measurement to the next:

- Go to the **Measurement** screen
- Using the knurled knob, select the measurement symbol in the top left
- Briefly press the knurled knob to change the value

### TRMS - Direct Selection:

For measurement of current or voltage, it is possible to measure the direct or TRMS value.

To modify the mode:

- Go to the Measurement screen
- Select TRMS or Direct in the top left of the screen
- Briefly press the knurled knob to modify the mode

### Energy measurement:

In Table or Chart mode, select this symbol  and press the knurled knob to start energy measurement. Time measurement will also appear at the top of the screen.

Energy measurement lasts for 60 seconds by default. To modify this time:

- Press the Menu button
- Select Acquisition and time
- Select the duration box and modify the time using the knurled knob.
- Select OK to return to the main screen

## 2.4 Data acquisition

The Initio®2 has an internal memory which allows it to save measurements. Two modes are available:

- Manual acquisition

Acquisition takes place on a point by point basis. Each measurement is manually saved by the user (see below)

- Time acquisition

Data acquisition takes place as a function of time, according to a duration entered by the user

### 2.4.1 Manual acquisition

- Press the Menu button (4)
- Using the knurled knob (3), select **Acquisition** and confirm by pressing on the knurled knob (3)
- Select **Manual**
- To save a measurement, select the **Save** icon (white circle) and confirm by briefly pressing the knurled knob (3).
- To stop the acquisition select the **Stop** icon (black square) and confirm by pressing the knurled knob (3). Confirm the choice in the dialogue box.

As soon as you press the Save icon, the energy meter starts up

### 2.4.2 Time acquisition

- Press the Menu button (4)
- Using the knurled knob (3), select **Acquisition** and confirm by pressing on the knurled knob (3)
- Select **Time**
- Using the knurled knob (3), select the number defining the duration and press on the knurled knob (3).
- Modify this duration by turning the knurled knob (3) then confirm by pressing on the centre of the knurled knob.
- If the data acquisition is to loop (return of the acquisition to zero at the end of a defined period) select **Loop**
- Confirm the final choice by selecting **OK**

Time data acquisitions last for a time between 8 seconds and 48 hours. The interval between two measurements adjusts automatically.

### 2.4.3 Recording the data

Once the acquisition is completed (time or manual) it must be saved to be able to retrieve the points on a computer.

- Press the Menu button (4)
- Select **File** and confirm
- Select **Save**

Only one data acquisition can be saved. By selecting **Open**, you open the last saved data acquisition.

### 2.4.4 Transfer of the data to the computer

Once the acquisition is completed, it is possible to transfer the data to a computer to process it:

- Install the Initio®2 software available under the address [www.jeulin.fr/Initio2](http://www.jeulin.fr/Initio2)
- Connect the Initio®2 to the computer via the USB-μUSB cable
- Open the Initio® 2 software
- Load the saved acquisition by clicking the icon:

The data then appears on the screen. Right-clicking the chart provides access to various chart tools.

## 3. Data visualisation

It is possible to view the data at any moment in tabular or chart format. To do so, simply press the Back button (2). The display then toggles between measurement, table and chart. Whatever the current display, the data acquisition can be monitored.

After 5 seconds without activity, table or chart is output as a full screen display. To again obtain the associated functions, turn the knurled knob.

In the table, the arrows situated at the bottom right are used to move between the rows of the table.

## 4. Analogue outputs

The Initio2 has analogue outputs in the form of banana jacks and a 3.5 mm jack. The U and I analogue outputs deliver a voltage of -5V to +5V for full-scale grades, whether direct or TRMS measurement.

**Output 1 :** top and middle socket => voltage U based on a full scale of 0-30V i.e. S1 = -5V for U = 0V, S1 = 0V for U = 15V and S1 = 5V for U = 30V

We can assimilate the response of the analog output to an affine function  $S1 = (10/30) \times U - 5$

**Output 2 :** bottom and middle socket => current I based on a full scale from 0 to 5A i.e. S2 = -5V for I = 0A, S2 = 0V for I = 2.5A and S2 = 5V for I = 5A

We can assimilate the response of the analog output S2 to an affine function  $S2 = (10/5) \times I - 5$

## 5. Multi display mode

Select a quantity at the top left of the Joulemeter screen.  
Press and hold the OK button to activate this mode allowing to simultaneously display U, I, P and E (denoted W).

## 6. WIFI communication

Using the WIFI (451 470) module, it is possible to visualise the Initio2 data on a tablet or similar.

## 6.1 Installation of the WIFI module

The WIFI module is fitted with the Initio®2 switched off.

- Lift the Initio® 2 handle
- Remove the back by sliding it down
- Position the WIFI module instead of the back
- Push the module firmly into place until it clicks. It must be possible to retract the handle without forcing it.

## 6.2 Visualisation on a tablet

Once the WIFI module is installed, the tablet must be configured:

- Press the Menu button
- Select Parameters
- Select Network

Information concerning the module is output to the screen.

- On the tablet, chose the name of the WIFI network displayed on the Initio®2.
- Open a web browser
- Enter the address displayed by the Initio®2 in the address bar of the browser.

## 7. Parameters

Various information concerning the Initio®2 can be obtained by pressing on Menu and selecting Parameter → Information

This page displays the serial number of the instrument and the software version, the battery charge level and the operating language. The latter is changed by turning the knurled knob.

## 8. Update

On our website at the address <http://en.jeulin.fr/en-resources/english-resources/software/free-software.html>, the last update for your Initio® 2 is available. Its version is indicated in the designation.

In order to know the version of your Initio® 2, press the “Menu” button and then “Setup” “Informations”.

## 9. Technical characteristics

Function	Range	Resolution	Detail
Direct ammeter	± 5A	1mA	2%
TRMS ammeter	0/5A	1mA	3%
Direct voltmeter	±10V ±30V	1mV 10mV	2%
TRMS voltmeter	0/10V 0/30V	1mV 10mV	3%
Wattmeter	±100W ±150W	10 mW 100mW	Sum of U and I details
Joulemeter	0/13MJ	From 1J to 10kJ	

### Clock (display h:m:s)

Range	Resolution	Detail
0.0 - 2 s	1 ms	± 1 UR
2.00 - 20 s	10 ms	± 1 UR
20.0 - 200 s	100 ms	± 1 UR
200 - 2000 s	1 s	± 1 UR

### Protection:

Ammeter protected by an automatic resettable fuse (cut-off power: 40 A / 30 V)

Max. operating voltage: 30 V peak

Max. operating current: 5 A peak

### Supply:

Li-ion battery with sector adapter and micro-USB cable (supplied)

### Analogue output:

-5 V to +5 V via 4 mm safe banana jacks

0 to +5 V via 3.5 mm jack

### Sizes:

95 x 123 x 44 mm

## 10. After-sales service

The device is under a 2-year guarantee, it must be sent back to our workshops.

For any repairs, adjustments or spare parts please contact:

### JEULIN – TECHNICAL SUPPORT

468 rue Jacques Monod

CS 21900

27019 EVREUX CEDEX FRANCE

**+33 (0)2 32 29 40 50**