Electricité
Electricity

**Ref :
302 047**

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Teaching scale model for Ohm's law

Français – p 1

English – p 5

Version : 7007

1 Description

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NOTES

1 Description

The teaching scale model for Ohm's law is a complete unit for tracing the graph of the characteristics of an Ohmic dipole. Serigraphy allows simple connection of the ammeter and the voltmeter without risk of error.

A multi-purpose connector of the loud speaker type, is used to connect the ohmic dipoles.

A potentiometer is incorporated in the scale model to vary the voltage from 1.3 V to ~ 9 V.

Accessory for attaching a bulb with an E10 base.

2 Additional equipment

Recommended additions

- 12 V DC power supply
- 1 ammeter
- 1 voltmeter

Consumables

- Bulbs, base E10, 6 V
- Bare resistances : from 120 to 1000 Ω

3 Technical characteristics

Connection: Power supply for the ammeter and voltmeter on a safety socket diameter 4 mm.

Protection: Against reversal of polarity, against overheating: voltage regulator.

Accuracy: +/- 5%.

Pmax: 1W.

Dimensions: 105 x 65 x 50 mm.

4 Assembly and installation

Connection of the ohmic conductor

A multipurpose connector of the loud speaker type is used to connect all types of resistances.

A support is also delivered for screwing in a 6 V, 100 mA bulb.

Measuring

Follow the instructions on the unit:

The ammeter at **A**: Rating 200 mA for a resistance of ~ 120 Ω and 20 mA for resistances between 400 – 1000 Ω .

The voltmeter at **B**: Rating 2 V.

Power supply

Use a 12 V DC power supply.

The potentiometer allows varying the input voltage continuously between 1.25 and 9 V.

Example:

1. U constant, R varies → measure I versus the ohmic dipole:

Thanks to the dipole support, at fixed voltage (here adjusted to 5 V using the potentiometer), we measured the intensity with 3 different resistances. The calculation R (resistance value) x I (measured) gives us the original voltage.

| Resistance | U | I measured (mA) | R x I calculated (V) |
|--------------|-----|-----------------|----------------------|
| 120 Ω | 5 V | 41.2 | 4.94 |
| 470 Ω | 5 V | 10.8 | 5.07 |
| 1 k Ω | 5 V | 5.05 | 5.05 |


2. Variation of U → measurement of I:

Using the potentiometer we varied the voltage at the lamp terminals and measured the corresponding intensity.

| U (V) | I (mA) |
|-------|--------|
| 1,5 | 45.8 |
| 3 | 69.4 |
| 4.5 | 87 |
| 6 | 102.2 |
| 7.5 | 115.9 |
| 9 | 123.9 |

Result obtained: The graph of the voltage versus intensity gives a straight line:

Study of Ohm's law $U=R'I$

5 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

JEULIN - TECHNICAL SUPPORT
Rue Jacques Monod
BP 1900
27 019 EVREUX CEDEX FRANCE
+33 (0)2 32 29 40 50

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediatly to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service * :
+ 33 (0)2 32 29 40 50

** Call cost not included*